**Polyclonal Rabbit Anti-Human Kappa Light Chains/APC, Rabbit F(ab')<sub>2</sub>**  
**Polyclonal Rabbit Anti-Human Kappa Light Chains/FITC, Rabbit F(ab')<sub>2</sub>**  
**Polyclonal Rabbit Anti-Human Kappa Light Chains/RPE, Rabbit F(ab')<sub>2</sub>**

**Code C0222**  
**Code F0434**  
**Code R0436**

## ENGLISH

### Intended use

For in vitro diagnostic use.

C0222, F0434 and R0436 are intended for use in flow cytometry. In flow cytometry, antibodies to kappa light chains are useful for the demonstration of cell surface kappa light chains, and, thus, for the identification of monoclonality (clonal excess) in B-cell lymphoproliferative disorders together with a panel of other antibodies (1-4). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

### Introduction

Most B cells, with the exception of pre-B progenitors and pre-B cells, and mature plasma cells, express immunoglobulin on their surface. Each cell expresses only one light chain type. In normal peripheral blood and lymph nodes, there is a mixture of kappa-positive and lambda-positive cells, with two-thirds of the cells expressing kappa and one-third expressing lambda (5). Since lymphoid neoplasms are usually clonal expansions of a single cell, malignant cells uniformly express the same light chain isotype. Neoplastic B-cell lymphoproliferative disorders can frequently be suspected on the basis of the demonstration of a marked predominance of cells expressing a single light chain type (4). Rare cases of kappa-negative or lambda-negative B-cell non-Hodgkin's lymphoma have been described (4, 6).

### Reagent provided

The Anti-Kappa Light Chains conjugates, C0222, F0434 and R0436, have been produced from a F(ab')<sub>2</sub> fragment of affinity-isolated polyclonal rabbit antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Use 10 µL of the product for staining of up to 10<sup>6</sup> leucocytes from normal human peripheral blood.

Conjugate concentration g/L: See label on vial.

Antibody Code	Fluorochrome	Ig Reagent Code
C0222	APC (Allophycocyanin)	X0998
F0434	FITC (Fluorescein Isothiocyanate Isomer 1)	X0929
R0436	RPE (R-Phycoerythrin)	X0930

### Preparation

1. The antibody used for APC, FITC and RPE conjugation has been solid-phase absorbed with human plasma proteins to remove traces of contaminating antibodies.
2. The absorbed antibody has been further purified by affinity chromatography on a column with immobilized human kappa light chains.
3. The affinity-isolated antibody has then been degraded with pepsin and the F(ab')<sub>2</sub> fragment isolated by gel filtration.
4. Finally, the F(ab')<sub>2</sub> fragment has been conjugated with allophycocyanin, fluorescein isothiocyanate isomer 1 or R-phycoerythrin.

### Immunogen

Polyclonal immunoglobulin light chains of kappa type isolated from a pool of human sera.

### Specificity

Anti-Kappa Light Chains reacts with free kappa chains as well as kappa chains in intact immunoglobulin molecules. The specificity has been ascertained as described below. To obtain maximum sensitivity, the crossed immunoelectrophoresis specificity test was performed prior to affinity purification and pepsin degradation.

Crossed immunoelectrophoresis: Only kappa-related precipitates appear when the antibody is run against human plasma. Staining: Coomassie Brilliant Blue.

Flow cytometry: When Anti-Kappa Light Chains is applied as described in the staining procedure in combination with Anti-CD19/RPE, HD37, or Anti-CD19/FITC, HD37 on lysed human whole blood, a specific staining of a part of the CD19-positive B lymphocytes is seen corresponding to the expected range of kappa light chain expression.

Flow cytometric analysis of single cell suspensions from formalin-fixed, paraffin-embedded tissue samples demonstrated that Anti-Kappa Light Chains labels reactive hyperplastic lymph nodes (10/10 cases). In non-Hodgkin's lymphomas, the antibody labeled 4/10 cases (5). The remaining cases were positive for lambda light chains (5/10 cases) or showed no expression of light chains (1/10 case) (4).

Flow cytometric analysis of peripheral blood lymphocytes demonstrated that Anti-Kappa Light Chains labels a proportion of B-cell chronic lymphocytic leukemias. Thus, in one study of 121 cases (2), 37 were positive for kappa. In another study of 165 cases (3), 97 were positive for kappa.

### Precautions

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused solution should be disposed of according to local, State and Federal regulations.

### Storage

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. During storage a small precipitate may occasionally develop causing a fine granular non-specific staining. By a simple filtration (0.22 µm cellulose acetate filter), the original high quality of the conjugate will be restored. Conjugates should not be stored in diluted form. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Services.

### Procedural note

Before staining samples of peripheral blood, the mononuclear cells must be isolated by centrifugation on a separation medium or the blood sample washed to remove the soluble serum proteins. As human monocytes bind serum immunoglobulins via their surface Fc receptors, these cells should be removed by depletion or identified.

### Staining procedure

1. Collect venous blood into a test tube containing an anticoagulant.
2. Transfer 100 µL of the anticoagulated blood into the test tube.
3. Add 2 mL of 0.01 mol/L PBS, pH 7.4. Mix gently by using a vortex mixer.
4. Centrifuge at 300 x g for 5 minutes, then aspirate the supernatant leaving approximately 50 µL of fluid.
5. Repeat steps 3 and 4 two more times.
6. Add 10 µL rabbit immunoglobulin fraction, Code X0903, for blocking. Mix gently by using a vortex mixer and incubate in the dark at 37 °C for 30 minutes.

7. Add 10 µL of fluorochrome-conjugated Anti-Kappa Light Chains. Mix gently.
8. Use a non-reactive F(ab')<sub>2</sub> fragment of rabbit immunoglobulin, and conjugated with the same fluorochrome, as a negative control (see table).
9. Incubate in the dark at 4 °C for 30 minutes.
10. Add 1-2 mL of erythrocyte-lysing reagent to each tube and mix gently. Follow the reagent manufacturer's recommendations for time and temperature of incubation.
11. Centrifuge at 300 x g for 5 minutes.
12. Aspirate the supernatant, leaving approximately 50 µL of fluid.
13. Add 2 mL 0.01 mol/L PBS, pH 7.4. Mix gently by using a vortex mixer.
14. Repeat steps 11 and 12.
15. Resuspend pellet in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
16. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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

1. Johnson A, Olofsson T. Flow cytometric clonal excess analysis of peripheral blood, routine handling, and pitfalls in interpretation. *Cytometry* 1993;14:188-95.
2. Cartron G, Linassier C, Bremond JL, Desablens B, Georget MT, Fimbel B, et al. CD5 negative B-cell chronic lymphocytic leukemia: clinical and biological features of 42 cases. *Leuk Lymphoma* 1998;31:209-16.
3. Gandini D, Lanza F, Latorraca A, Levato F, Del Senno L, Castoldi G. Immunophenotypic and genotypic characterization of B-cell chronic lymphocytic leukemia patients from Northern Italy. *Haematologica* 1993;78:18-24.
4. Leers MPG, Theunissen PHMH, Ramaekers FCS, Schutte B, Nap M. Clonality assessment of lymphoproliferative disorders by multiparameter flow cytometry of paraffin-embedded tissue: an additional diagnostic tool in surgical pathology. *Hum Pathol* 2000;31:422-7.
5. Deegan MJ. B lymphocytes and plasma cells: their development and identification. In: Keren DF, editor. *Flow cytometry in clinical diagnosis*. Chicago: ASCP Press; 1989. p. 139-63.
6. Kaleem Z, Zehnauer BA, White G, Zutter MM. Lack of expression of surface immunoglobulin light chains in B-cell non-Hodgkin lymphomas. *Am J Clin Pathol* 2000;113:399-405.

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Agilent Technologies Singapore (International) Pte Ltd.