## Monoclonal Mouse Anti-Human CD20, Clone B-Ly1

Code No./ Code/ Code-Nr. **F 0799** FITC-Conjugated  
 Code No./ Code/ Code-Nr. **R 7013** RPE-Conjugated  
 Code No./ Code/ Code-Nr. **C 7132** RPE-Cy5-Conjugated

### ENGLISH

**Intended use** For in vitro diagnostic use.  
 F 0799, R 7013 and C 7132 are intended for use in flow cytometry. CD20 can be considered a pan-B cell antigen, as it is expressed on the surface of all mature B-lymphocytes but not in secreting plasma cells (1). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

**Synonyms for antigen** B1, Bq35 (2).

**Introduction** The CD20 antigen is a non-glycosylated 33-37 kDa membrane phosphoprotein (3). It contains three extensive hydrophobic regions which pass through the cell membrane four times. In addition, the long carboxyl and amino ends of the molecule are located intracellularly with only a minor portion of the molecule exposed outside the cell (4).  
 The CD20 molecule is involved in regulation of B cell activation, proliferation and differentiation (1). The CD20 antigen is expressed early during pre-B cell development, presumably just before the expression of cytoplasmic  $\mu$ -chains, and it persists until plasma cell differentiation (3). Weak expression of CD20 may be seen in a subpopulation of T cells but the antigen is not expressed in normal myeloid, erythroid, monocytic or mesenchymal cells (1).  
 The CD20 antigen is expressed in acute lymphoblastic leukaemia, B-cell chronic lymphocytic leukaemia, hairy cell leukaemia and Burkitt's lymphoma (5).

**Reagent provided** The Anti-CD20 conjugates, F 0799, R 7013 and C 7132, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 100 tests (10  $\mu$ L of conjugate for up to 10<sup>6</sup> leucocytes from normal human peripheral blood).  
 Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
F 0799	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927
R 7013	RPE (R-Phycoerythrin)	X 0928
C 7132	RPE-Cy5 (R-Phycoerythrin + Cyanine 5)	X 0955

**Specificity** Anti-CD20, B-LY1, was included in the Fourth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with the CD20 antigen (3). Cross-blocking studies indicate that a majority of antibodies in the CD20 cluster studied at the above workshop, defines a single epitope (4).

**Precautions** 1. For professional users.  
 2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.  
 3. As with any product derived from biological sources, proper handling procedures should be used.

**Storage** Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Staining procedure**

- Collect venous blood into a test tube containing an anticoagulant.
- Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
- Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
- Mix 100  $\mu$ L cell suspension with 10  $\mu$ L fluorochrome-conjugated Anti-CD20.
- Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
- Incubate in the dark at 4 °C for 30 minutes.
- Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
- Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

**Product-specific limitations** It has been observed that RPE-Cy5-conjugates may bind to monocytes, resulting in background staining (6).

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
- Leong A.S-Y, Cooper K, Leong F.J.W-H. Manual of diagnostic antibodies for immunohistology; London: Oxford University Press; 2003. p. 89-91.
- Tedder TF. CD Guide. CD20. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 766.
- Dörken B, Möller P, Pezzutto A, Schwartz-Albiez R, Moldenhauer G. B3. B-cell antigens: CD20. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p 46-8.
- Tedder TF, Penta A. B3.1. Structure of the CD20 antigen and gene of human and mouse B-cells: use of transfected cell lines to examine the Workshop panel of antibodies. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p 48-50.
- Pezzutto A, Behm F, Callard RE, Clark EA, Genetet N, Goodall AH, et al. B22.3. Flow cytometry analysis of the B-cell blind panel: joint report. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p 165-74.
- van Vugt MJ, van den Herik-Oudijk IE, van de Winkle JGJ. Binding of PE-CY5 conjugates to human high-affinity receptor for IgG (CD64). Blood 1996;88(6):2358-61.

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