Monoclonal Mouse Anti-Human CD19, Clone HD37

Code No./ Code/ Code-Nr. C 7224 APC-Conjugated
Code No./ Code/ Code-Nr. F 0768 FITC-Conjugated
Code No./ Code/ Code-Nr. R 0808 RPE-Conjugated
Code No./ Code/ Code Nr. C 7066 RPE-Cy5-Conjugated

ENGLISH

Intended use For in vitro diagnostic use.
 C 7224, F 0768, R 0808 and C 7066 are intended for use in flow cytometry. CD19 is the broadest lineage-specific surface marker for B cells and it is present on the surface of virtually all B lymphocytes, including early B progenitor cells (1). CD19 expression is maintained in B-lineage cells that have undergone neoplastic transformation (2). Antibodies to CD19 are considered essential for the initial evaluation of acute and chronic lymphoproliferative disorders (3). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Introduction CD19 is a 120 kDa transmembrane glycoprotein with a reduced Mr of 95 000 (2, 4). CD19 is a member of the immunoglobulin superfamily with two extracellular C2-type domains (2), and it is a critical signal transduction molecule that regulates B lymphocyte development, activation, and differentiation (4). CD19 expression is restricted to normal and neoplastic B cells, being absent from T cells, monocytes, and granulocytes (2). The CD19 antigen appears early during B cell maturation, probably at late pro-B cell stage around the time of Ig heavy chain rearrangement. It then persists during all stages of B cell maturation and is lost upon terminal differentiation to plasma cells (2). The advantage of immunodetection of CD19 expression is that B lineage leukaemias and lymphomas rarely lose the epitope (5). Expression of CD19 in some acute myelocytic leukaemia cells has been described (2).

Reagent provided The Anti-CD19 conjugates, C 7224, F 0768, R 0808, and C 7066, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).

Isotype: IgG1, kappa. **Conjugate concentration mg/L:** see label on vial.

Antibody Code No.	Fluorochrome	Negative Code No.	Control
F 0768	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927	
R 0808	RPE (R-Phycoerythrin)	X 0928	
C 7066	RPE-Cy5 (R-Phycoerythrin-Cyanine 5)	X 0955	
C 7224	APC (Allophycocyanin)	X 0968	

Immunogen Hairy cell leukaemia cells (6).

Specificity Anti-CD19, HD37, was included in the Second, Third, Fourth and Fifth International Workshops and Conferences on Human Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with CD19 (7, 8).

Anti-CD19, HD37, labels human B cells in peripheral blood, bone marrow and other tissues. Additionally, B-cell lymphoproliferative disorders gave positive reactions with the antibody, i.e. acute lymphoblastic leukaemia, chronic lymphocytic leukaemia, hairy cell leukaemia, lymphoblastic lymphoma (Burkitt type), centroblastic/centrocytic lymphoma, and diffuse non-Hodgkin's lymphoma. The antibody was shown to be unreactive with other cells in the human haematopoietic system and did not react with non-haematopoietic cells, e.g. in kidney, liver, breast or lung tissues (8). Anti-CD19, HD37, inhibits anti-Ig-induced B-cell activation and proliferation (1).

Precautions 1. For professional users.
 2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
 3. As with any product derived from biological sources, proper handling procedures should be used.

Storage Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Staining procedure 1. Collect venous blood into a test tube containing an anticoagulant.
 2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
 3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
 4. Mix 100 µL cell suspension with 10 µL of the fluorochrome-conjugated Anti-CD19.
 5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
 6. Incubate in the dark at 4 °C for 30 minutes.
 7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
 8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Product-specific limitations It has been observed that RPE-Cy5-conjugates may bind to monocytes, resulting in background staining (9).

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