**Monoclonal Mouse Anti-Human CD3/APC**, Clone UCHT1  
**Monoclonal Mouse Anti-Human CD3/FITC**, Clone UCHT1  
**Monoclonal Mouse Anti-Human CD3/PerCP**, Clone UCHT1  
**Monoclonal Mouse Anti-Human CD3/RPE**, Clone UCHT1  
**Monoclonal Mouse Anti-Human CD3/RPE-Cy5**, Clone UCHT1

**Code C7225**  
**Code F0818**  
**Code PR702**  
**Code R0810**  
**Code C7067**

## ENGLISH

### Intended use

For in vitro diagnostic use.

C7225, F0818, PR702, R0810 and C7067 are intended for use in flow cytometry. CD3 is a pan-T-cell-lineage-restricted antigen, which is a valuable marker for normal and neoplastic T cells. Antibodies to CD3 are considered essential for the initial evaluation of acute and chronic lymphoproliferative disorders (1). Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

### Summary and explanation

Human TCR/CD3 is a complex structure on the lymphocyte surface. It consists of the TCRαβ or TCRγδ heterodimer and the associated CD3 complex. The CD3 complex is composed of six polypeptides with usually four different transmembrane CD3 chains, γ (gamma), δ (delta), ε (epsilon), and ζ (zeta). Three different dimers, γε, δε, and ζζ, constitute the CD3 complex. The Mr of CD3ε is 20 000 (2).

The CD3 complex is crucial in transducing antigen-recognition signals into the cytoplasm of T cells and in regulating the cell surface expression of the TCR complex. Further it plays an important role in the differentiation of thymocytes (2).

CD3 is first detectable in early thymocytes and its appearance probably represents one of the earliest signs of commitment to the T cell lineage. In cortical thymocytes, during early stages of maturation, the CD3 antigen is predominantly present in the cell cytoplasm. In medullary thymocytes, the CD3 antigen is predominantly detected on the cell surface (3, 4).

The majority of T cell neoplasms also expresses CD3, but it is absent from non-T-cell lymphoid malignancies (5). Consistent with the pattern of synthesis of the antigen in normal thymocytes, the earliest site detectable within neoplastic cells is the cell cytoplasm (3).

### Reagent provided

The Anti-CD3 conjugates, C7225, F0818, PR702, R0810 and C7067, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10<sup>6</sup> leucocytes from normal peripheral blood).

**Isotype:** IgG1, kappa. **Conjugate concentration mg/L:** see label on vial.

Antibody Code	Fluorochrome	Negative Control Code
F0818	FITC (Fluorescein Isothiocyanate Isomer 1)	X0927
R0810	RPE (R-Phycoerythrin)	X0928
PR702	PerCP (Peridinin Chlorophyll Protein)	X7909
C7067	RPE-Cy5 (R-Phycoerythrin-Cyanine 5)	X0955
C7225	APC (Allophycocyanin)	X0968

### Immunogen

Human infant thymocytes and lymphocytes from a patient with Sézary disease (6).

### Specificity

Anti-CD3, UCHT1, was included in the First and Third International Workshops and Conferences on Human Leucocyte Differentiation Antigens (Paris 1982, Oxford 1986), and studies by a number of laboratories confirmed its reactivity with the CD3 antigen (7). Anti-CD3, UCHT1, reacts with the 20 kDa ε-chain of CD3 (8).

Anti-CD3 labels T cells in thymus, bone marrow, spleen, tonsil and blood (3-6). It also labels Purkinje cells in the cerebellum, which is the only other cell type known to bind antibodies to CD3 (9).

The antibody is able to induce in vitro proliferation of mature thymocytes and T cells in the presence of interleukin-2 (IL-2) (10).

### Precautions

- For professional users.
- This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
- As with any product derived from biological sources, proper handling procedures should be used.
- Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
- Unused solution should be disposed of according to local, State and Federal regulations.

### Storage

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

### Staining procedure

- Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.
- Add 10 µL of fluorochrome-conjugated anti-CD3 and mix gently by using a vortex mixer.
- Incubate the tube in the dark at 2-8 °C for 30 min utes or at room temperature (20-25 °C) for 15-30 mi nutes.
- Add 100 µL of Uti-Lyse™ Reagent A (Dako Code S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
- Add 1 mL of Uti-Lyse™ Reagent B (Dako Code S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
- Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
- Add 2 mL of PBS to the tube and resuspend the cells by using a vortex mixer.
- Repeat step 6.
- Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.
- Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples shoul d be analysed within 24 hours after staining.

Note that fluorochrome conjugates are light sensitive. Reagents should be protected from light during storage and samples should be protected from light during staining procedure and until analysis.

(102361-003)

C7225/F0818/PR702/R0810/C7067/EFG/LHP/2008.09.30 p. 1/4

Dako Denmark A/S | Produktionsvej 42 | DK-2600 Glostrup | Denmark | Tel. +45 44 85 95 00 | Fax +45 44 85 95 95 | CVR No. 33 21 13 17

### Procedural notes

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the isotype and fluorochrome of the conjugated antibody. Recommended control reagents are shown in the table above

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, such as Dako EasyLyse™, Code S2364, the PBS in step 9 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

### Product-specific limitations

It has been observed that RPE-Cy5-conjugates may bind to monocytes, resulting in background staining (11).

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
<p>References/ Références/ Literatur</p>												
<ol style="list-style-type: none"><li>Braylan RC, Orfao A, Borowitz MJ, Davis BH. Optimal number of reagents required to evaluate hematologyphoid neoplasias: results of an international consensus meeting. Cytometry 2001;46:23-7.</li> <li>Saito T, Yamazaki T. TC4. CD3 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 44-48.</li> <li>Campana D, Thompson JS, Amlot P, Brown S, Janossy G. The cytoplasmic expression of CD3 antigens in normal and malignant cells of the T lymphoid lineage. J Immunol 1987;138:648-55.</li> <li>Swerdlow SH, Angermeier PA, Hartman AL. Intrathymic ontogeny of the human T cell receptor associated CD3 (T3) antigen. Lab Invest 1988;58:421-7.</li> <li>Erber WN, Mynheer LC, Mason DY. APAAP labelling of blood and bone-marrow samples for phenotyping leukaemia. Lancet 1986;i:761-5</li> <li>Beverley PC, Callard RE. Distinctive functional characteristics of human "T" lymphocytes defined by E rosetting or a monoclonal anti-T cell antibody. Eur J Immunol 1981;11:329-34.</li> <li>McMichael AJ, Gotch FM. T-cell antigens: new and previously defined clusters. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sept 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p.31-62.</li> <li>Tunnacliffe A, Olsson C, Traunecker A, Krissansen GW, Karjalainen K, de la Hera A. T3.2. The majority of CD3 epitopes are conferred by the epsilon chain. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 295-6.</li> <li>Garson JA, Beverley PC, Coakham HB, Harper EI. Monoclonal antibodies against human T lymphocytes label Purkinje neurones of many species. Nature 1982;298:375-7.</li> <li>Denning SM, Tuck DT, Singer KH, Haynes BF. T2.11. Activation of human thymocytes via CD3 and CD2 molecules. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM et al, editor. Leukocyte Typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference. Oxford-New York-Tokyo: Oxford University Press; 1987. p. 144-7.</li> <li>van Vugt MJ, van den Herik-Oudijk IE, van de Winkel JGJ. Binding of PE-CY5 conjugates to the human high-affinity receptor for IgG (CD64). Blood 1996;88:2358-61.</li></ol>												
<p>Explanation of symbols/ Légende des symboles/ Erläuterung der Symbole</p>												
<table> <tbody><tr> <td><b>REF</b></td> <td>Catalogue number Référence du catalogue Bestellnummer</td> <td> <p>Temperature limitation Limites de température Zulässiger Temperaturbereich</p></td> <td> <p>Use by Utiliser jusque Verwendbar bis</p></td></tr> <tr> <td><b>IVD</b></td> <td>In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum</td> <td> <p>Keep away from sunlight (consult storage section) Conservser à l'écart du soleil (se reporter à la section conservation) Lichtgeschützt lagern (siehe Abschnitt zur Lagerung)</p></td> <td> <p>Manufacturer Fabricant Hersteller</p></td></tr> <tr> <td></td> <td>Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten</td> <td><b>LOT</b></td> <td>Batch code Code du lot Chargenbezeichnung</td></tr> </tbody></table>	<b>REF</b>	Catalogue number Référence du catalogue Bestellnummer	 <p>Temperature limitation Limites de température Zulässiger Temperaturbereich</p>	 <p>Use by Utiliser jusque Verwendbar bis</p>	<b>IVD</b>	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	 <p>Keep away from sunlight (consult storage section) Conservser à l'écart du soleil (se reporter à la section conservation) Lichtgeschützt lagern (siehe Abschnitt zur Lagerung)</p>	 <p>Manufacturer Fabricant Hersteller</p>		Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	<b>LOT</b>	Batch code Code du lot Chargenbezeichnung
<b>REF</b>	Catalogue number Référence du catalogue Bestellnummer	 <p>Temperature limitation Limites de température Zulässiger Temperaturbereich</p>	 <p>Use by Utiliser jusque Verwendbar bis</p>									
<b>IVD</b>	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	 <p>Keep away from sunlight (consult storage section) Conservser à l'écart du soleil (se reporter à la section conservation) Lichtgeschützt lagern (siehe Abschnitt zur Lagerung)</p>	 <p>Manufacturer Fabricant Hersteller</p>									
	Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	<b>LOT</b>	Batch code Code du lot Chargenbezeichnung									

(102361-003)	C7225/F0818/PR702/R0810/C7067/EFG/LHP/2008.09.30 p. 4/4
Dako Denmark A/S   Produktionsvej 42   DK-2600 Glostrup   Denmark   Tel. +45 44 85 95 00   Fax +45 44 85 95 95   CVR No. 33 21 13 17	