## Monoclonal Mouse Anti-Human CD4, Clone MT310

Code No./ Code/ Code-Nr. C 7226 APC-Conjugated  
 Code No./ Code/ Code-Nr. F 0766 FITC-Conjugated  
 Code No./ Code/ Code-Nr. R 0805 RPE-Conjugated  
 Code No./ Code/ Code-Nr. C 7069 RPE-Cy5-Conjugated

### ENGLISH

**Intended use** For in vitro diagnostic use.  
 C 7226, F 0766, R 0805 and C 7069 are intended for use in flow cytometry. Measurement of the relative as well as the absolute number of CD4-positive cells in peripheral blood is critically important in assessing the immunological status of a patient. Additionally, antibodies to CD4 have an important role in the immunophenotyping of reactive lymphocytes and in the investigation of lymphoproliferative disorders (1). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

**Introduction** CD4 is a 55 kDa glycoprotein with four immunoglobulin-like external domains, a transmembrane domain and a highly conserved intracellular domain (1). CD4 acts as a coreceptor in MHC class II-restricted antigen-induced T-cell activation. Further, CD4 has a role in the regulation of T-B lymphocyte adhesion in the absence of antigen recognition, and in thymic differentiation (2).  
 CD4 is expressed on thymocyte subsets, it marks helper/inducer T cells, and it is expressed in 55-65% of mature peripheral T cells. Subpopulations of suppressor or cytotoxic T cells can be identified among CD4-positive T cells. CD4 is also expressed on monocytes/macrophages, Langerhans cells and other dendritic cells. CD4 is not expressed on B cells.  
 The majority of peripheral T-cell lymphomas are derived from the helper T-cell subset, so most postthymic T-cell neoplasms are CD4-positive (1).

**Reagent provided** The Anti-CD4 conjugates, C 7226, F 0766, R 0805 and C 7069, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10<sup>6</sup> leucocytes from normal human peripheral blood).  
Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code No.	Fluorochrome	Negative Code No.	Control
F 0766	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927	
R 0805	RPE (R-Phycoerythrin)	X 0928	
C 7069	RPE-Cy5 (R-Phycoerythrin + Cyanine 5)	X 0955	
C 7226	APC (Allophycocyanin)	X 0968	

**Specificity** Anti-CD4, MT310, was included in the Second International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with the CD4 antigen (3).

**Precautions**

- For professional users.
- This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
- As with any product derived from biological sources, proper handling procedures should be used.

**Storage** Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Staining procedure**

- Collect venous blood into a test tube containing an anticoagulant.
- Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
- Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
- Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD4.
- Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
- Incubate in the dark at 4 °C for 30 minutes.
- Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 0.01 mol/L PBS, pH 7.4.

8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

**Product-specific limitations** It has been observed that RPE-Cy5-conjugates may bind to monocytes, resulting in background staining (4).

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