Monoclonal Mouse Anti-Human CD8, Clone DK25

Code No./ Code/ Code-Nr. **C 7227 APC-Conjugated**
 Code No./ Code/ Code-Nr. **F 0765 FITC-Conjugated**
 Code No./ Code/ Code-Nr. **R 0806 RPE-Conjugated**
 Code No./ Code/ Code-Nr. **C 7079 RPE-Cy5-Conjugated**

| ENGLISH | | | | | | | | | | | | | | | | |
|---------------------------|---|---------------------------|--------------|---------------------------|--------|--|--------|--------|-----------------------|--------|--------|-------------------------------------|--------|--------|-----------------------|--------|
| Intended use | For in vitro diagnostic use. C 7227, F 0765, R 0806 and C 7079 are intended for use in flow cytometry. Measurement of the relative as well as the absolute number of CD8-positive cells in peripheral blood is critically important in assessing the immunological status of a patient (1). Additionally, CD8 has an important role in the immunophenotypic analysis of reactive and neoplastic populations of T cells (2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. | | | | | | | | | | | | | | | |
| Introduction | CD8 is a 68 kDa disulphide-linked, transmembrane glycoprotein. The CD8 molecule serves as a receptor for MHC class I molecules and can mediate a function as a coreceptor in TCR-ligand binding and T-cell activation. Its cytoplasmic part is associated with the p56 ^{lck} tyrosine kinase (1). CD8 is expressed by the great majority of cortical thymocytes and approximately 30% of medullary thymocytes, and by class I major histocompatibility complex restricted, mature suppressor/cytotoxic T cells. In addition a proportion of $\gamma\delta$ T cells and NK cells express CD8 (1). Thus, CD8 is expressed by about 25-35% of peripheral T cells, and about 30% of NK cells express low levels of CD8 (2). Immunophenotyping of lymphocytes is widely applied as a component of the testing of peripheral blood for diagnosis of immunodeficiencies such as acquired immunodeficiency syndrome (AIDS) (3, 4). An increase in the percentage of CD8-positive T cells has been reported during cytomegalovirus infection (5). | | | | | | | | | | | | | | | |
| Reagent provided | The Anti-CD8 conjugates, C 7227, F 0765, R 0806 and C 7079, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 μ L of conjugate for up to 10 ⁶ leucocytes from normal peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> see label on vial. | | | | | | | | | | | | | | | |
| | <table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 0765</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> <tr> <td>R 0806</td> <td>RPE (R-Phycoerythrin)</td> <td>X 0928</td> </tr> <tr> <td>C 7079</td> <td>RPE-Cy5 (R-Phycoerythrin-Cyanine 5)</td> <td>X 0955</td> </tr> <tr> <td>C 7227</td> <td>APC (Allophycocyanin)</td> <td>X 0968</td> </tr> </tbody> </table> | Antibody Code No. | Fluorochrome | Negative Control Code No. | F 0765 | FITC (Fluorescein Isothiocyanate Isomer 1) | X 0927 | R 0806 | RPE (R-Phycoerythrin) | X 0928 | C 7079 | RPE-Cy5 (R-Phycoerythrin-Cyanine 5) | X 0955 | C 7227 | APC (Allophycocyanin) | X 0968 |
| Antibody Code No. | Fluorochrome | Negative Control Code No. | | | | | | | | | | | | | | |
| F 0765 | FITC (Fluorescein Isothiocyanate Isomer 1) | X 0927 | | | | | | | | | | | | | | |
| R 0806 | RPE (R-Phycoerythrin) | X 0928 | | | | | | | | | | | | | | |
| C 7079 | RPE-Cy5 (R-Phycoerythrin-Cyanine 5) | X 0955 | | | | | | | | | | | | | | |
| C 7227 | APC (Allophycocyanin) | X 0968 | | | | | | | | | | | | | | |
| Specificity | The specificity of Anti-CD8, DK25, is equivalent to that of the CD8-clustered antibodies OKT8 and Leu-2a (6). | | | | | | | | | | | | | | | |
| Precautions | <ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. | | | | | | | | | | | | | | | |
| Storage | Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services. | | | | | | | | | | | | | | | |
| Staining procedure | <ol style="list-style-type: none"> Collect venous blood into a test tube containing an anticoagulant. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. Mix 100 μL cell suspension with 10 μL of the fluorochrome-conjugated Anti-CD8. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). Incubate in the dark at 4 °C for 30 minutes. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. Analyse on a flow cytometer. | | | | | | | | | | | | | | | |

Product-specific limitations It has been observed that RPE-Cy5-conjugates may bind to monocytes, resulting in background staining (7).

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- Nakauchi H. TC9. CD8 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 65-7.
- Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 1999. p. 55-6.
- Watret KC, Whitelaw JA, Froebel KS, Bird AG. Phenotypic characterization of CD8+ T cell populations in HIV disease and in anti-HIV immunity. Clin Exp Immunol 1993;92:93-9.
- Hengel RL, Nicholson JK. An update on the use of flow cytometry in HIV infection and AIDS. Clin Lab Med 2001;21:841-56.
- Collier AC, Meyers JD, Corey L, Murphy VL, Roberts PL, Handsfield HH. Cytomegalovirus infection in homosexual men. Relationship to sexual practices, antibody to human immunodeficiency virus, and cell-mediated immunity. Am J Med 1987;82:593-601.
- Stein H, Lennert K, Feller AC, Mason DY. Immunohistological analysis of human lymphoma: correlation of histological and immunological categories. Adv Cancer Res 1984;42:67-147.
- van Vugt MJ, van den Herik-Oudijk IE, van de Winkel JGJ. Binding of PE-CY5 conjugates to the human high-affinity receptor for IgG (CD64). Blood 1996;88:2358-61.

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