Monoclonal Mouse Anti-Human CD45, Leucocyte Common Antigen, Clone T29/33

Code No./ Code/ Code-Nr. **C 7230 APC-Conjugated**
 Code No./ Code/ Code-Nr. **F 0861 FITC-Conjugated**
 Code No./ Code/ Code-Nr. **R 7087 RPE-Conjugated**
 Code No./ Code/ Code-Nr. **C 7099 RPE-Cy5-Conjugated**

ENGLISH

Intended use

For in vitro diagnostic use.

C 7230, F 0861, R 7087 and C 7099 are intended for use in flow cytometry. CD45 is one of the most abundant leucocyte cell surface glycoproteins and is expressed exclusively on cells of the haematopoietic system and their progenitors (1). In flow cytometry, anti-CD45, together with a panel of other antibodies, is considered essential for the initial evaluation of chronic lymphoproliferative disorders and acute leukaemias (2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Introduction

CD45 is a single-chain type I transmembrane protein typically expressed at high levels on all nucleated cells of haematopoietic origin (3). Thus, CD45 is present on T and B-lymphocytes, granulocytes, monocytes and macrophages, with the exception of maturing erythrocytes and megakaryocytes (1). There are five different isoforms of CD45, based on differential splicing of exons 4, 5 and 6, named ABC, AB, BC, B and 0. The Mr of the isoforms ranges from 220 000 for the ABC isoform to 180 000 for the 0 isoform. All the CD45 isoforms share the same intracellular segment which has been shown to have tyrosine phosphatase activity and which has a functional role in lymphocyte activation and differentiation (3).

Antibodies that recognize all five isoforms are known as anti-CD45 (3).

Reagent provided

The Anti-CD45 conjugates, C 7230 F 0861, R 7087 and C 7099, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).

Isotype: IgG1, kappa. Conjugate concentration mg/L: see label on vial.

| Antibody Code No. | Fluorochrome | Negative Control Code No. |
|-------------------|--|---------------------------|
| F 0861 | FITC (Fluorescein Isothiocyanate Isomer 1) | X 0927 |
| R 7087 | RPE (R-Phycoerythrin) | X 0928 |
| C 7099 | RPE-Cy5 (R-Phycoerythrin-Cyanine 5) | X 0955 |
| C 7230 | APC (Allophycocyanin) | X 0968 |

Specificity

Anti-CD45, T29/33, was included in the Third and Fourth International Workshops and Conferences on Human Leucocyte Differentiation Antigens. At the Third Workshop, the antibody was clustered as a CD45 antibody reacting with **all** the known isotypes of the CD45 family, also called the leucocyte common antigen family (4). At the Fourth Workshop, the expression of CD45 in a range of haematopoietic cell lines, and the lack of CD45 in non-haematopoietic cell lines, were demonstrated. Notably, CD45 was expressed at a very low level in myeloma cells and in the U-266 plasmacytoid cell line (5).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

Storage

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Staining procedure

1. Collect venous blood into a test tube containing an anticoagulant.
2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
4. Mix 100 µL cell suspension with 10 µL of the fluorochrome-conjugated Anti-CD45.
5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
6. Incubate in the dark at 4 °C for 30 minutes.
7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Product-specific limitations It has been observed that RPE-Cy5-conjugates may bind to monocytes, resulting in background staining (6).

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
1. Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 1999. p. 95-8.
2. Braylan RC, Orfao A, Borowitz MJ, Davis BH. Optimal number of reagents required to evaluate hematolymphoid neoplasias: results of an international consensus meeting. *Cytometry* 2001;46:23-7.
3. Sewell WA, Cooley MA, Hegen M. NL6. CD45 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 499-502.
4. Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1, and leucocyte common antigens: new and previously defined clusters. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leukocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 788-803.
5. Schwinzer R. N7. Cluster report: CD45/CD45R. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leukocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 628-34.
6. van Vugt MJ, van den Herik-Oudijk IE, van de Winkel JGJ. Binding of PE-CY5 conjugates to the human high-affinity receptor for IgG (CD64). *Blood* 1996;88:2358-61.

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Agilent Technologies Singapore (International) Pte Ltd.