Monoclonal Mouse Anti-Human CD34/APC Class III, Clone BIRMA-K3
Monoclonal Mouse Anti-Human CD34/FITC Class III, Clone BIRMA-K3
Monoclonal Mouse Anti-Human CD34/RPE Class III, Clone BIRMA-K3
Monoclonal Mouse Anti-Human CD34/PerCP-Cy5.5 Class III, Clone BIRMA-K3

Code C7238
Code F7081
Code R7125
Code PR706

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>C7238, F7081, R7125 and PR706 are intended for use in flow cytometry. CD34 is expressed on early lymphohematopoietic stem and progenitor cells and also on endothelial cells. Antibodies to CD34 can be utilized to quantitate and purify lymphohematopoietic stem/progenitor cells for transplantation and research (1-6). CD34 is a useful marker for stem and progenitor cells as well as for the subclassification of leukemias. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.</p>																				
Summary and explanation	<p>CD34 is a single-chain transmembrane protein of approximately 116 kDa, expressed on immature hematopoietic stem/progenitor cells, capillary endothelial cells, embryonic fibroblasts and rare glial cells in nervous tissue (1, 5).</p> <p>CD34 appears to be expressed at its highest level on the earliest progenitors, and to decrease progressively with maturation (2, 4). CD34 is a stage-specific, rather than a lineage-specific, leucocyte differentiation antigen. The most immature definable B-lymphoid precursors (CD19+/CD10+/TdT+) are CD34+ (2, 4, 5). Immature T-lymphoid precursors also express TdT and CD34 (5). Normal peripheral blood lymphocytes, monocytes, granulocytes, and platelets do not express CD34 (4). Many acute B-lymphoid leukemias and acute myeloid leukemias (AML), and few acute T-lymphoid leukemias express CD34 (4). Chronic lymphoid leukemias, lymphomas and multiple myelomas may be found to be uniformly CD34 negative (2, 4, 5).</p> <p>Monoclonal antibodies to CD34 can be confined to three main classes, class I, class II and class III, defined by the differential sensitivity of the corresponding CD34 epitopes to degradation by specific enzymes. As an example, class III monoclonal antibodies recognize a CD34 epitope resistant to neuraminidase, chymopapain and glycoprotease (1, 3).</p>																				
Reagent provided	<p>The Anti-CD34 conjugates, C7238, F7081, R7125 and PR706, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial of C7238 contains conjugate for 50 tests of up to 10⁶ KG-1a cells. Each vial of F7081, R7125 and PR706 contains conjugate for 100 tests of up to 10⁶ KG-1a cells. See table 1 for volume of conjugate per test.</p> <p><u>Isotype</u>: IgG1, kappa. <u>Conjugate concentration</u>: see label on vial.</p>																				
	<p>Table 1</p> <table border="1"> <thead> <tr> <th>Antibody Dako Code</th> <th>Fluorochrome</th> <th>Volume of conjugate per test</th> <th>Isotype Reagent Dako Code*</th> </tr> </thead> <tbody> <tr> <td>F7081</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>10 µL</td> <td>X0927</td> </tr> <tr> <td>R7125</td> <td>RPE (R-Phycoerythrin)</td> <td>10 µL</td> <td>X0928</td> </tr> <tr> <td>C7238</td> <td>APC (Allophycocyanin)</td> <td>10 µL</td> <td>X0968</td> </tr> <tr> <td>PR706</td> <td>PerCP-Cy5.5 (Peridinin Chlorophyll Protein Cyanine 5.5)</td> <td>5 µL</td> <td>No Dako Code*</td> </tr> </tbody> </table>	Antibody Dako Code	Fluorochrome	Volume of conjugate per test	Isotype Reagent Dako Code*	F7081	FITC (Fluorescein Isothiocyanate Isomer 1)	10 µL	X0927	R7125	RPE (R-Phycoerythrin)	10 µL	X0928	C7238	APC (Allophycocyanin)	10 µL	X0968	PR706	PerCP-Cy5.5 (Peridinin Chlorophyll Protein Cyanine 5.5)	5 µL	No Dako Code*
Antibody Dako Code	Fluorochrome	Volume of conjugate per test	Isotype Reagent Dako Code*																		
F7081	FITC (Fluorescein Isothiocyanate Isomer 1)	10 µL	X0927																		
R7125	RPE (R-Phycoerythrin)	10 µL	X0928																		
C7238	APC (Allophycocyanin)	10 µL	X0968																		
PR706	PerCP-Cy5.5 (Peridinin Chlorophyll Protein Cyanine 5.5)	5 µL	No Dako Code*																		
	<p>*It is recommended to use a non-reactive monoclonal antibody of the same isotype and conjugated with PerCP-Cy5.5</p>																				
Immunogen	KG-1a cells (3).																				
Specificity	<p>Anti-CD34, BIRMA-K3, was included in the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens (Kobe 1996), and studies by different laboratories confirmed its reactivity with CD34 and the class III epitope (1).</p> <p>In flow cytometry the antibody labels KG-1a cells (a primitive myeloid leukemia cell line, known to be CD34+) (3).</p>																				
Precautions	<ol style="list-style-type: none"> For in vitro diagnostic use. For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. Unused solution should be disposed of according to local, State and Federal regulations. 																				
Storage	<p>Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.</p>																				
Staining procedure	<ol style="list-style-type: none"> Collect venous blood into a test tube containing an anticoagulant. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. Mix 100 µL cell suspension with the volume of fluorochrome-conjugated Anti-CD34 stated in table 1. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). Incubate in the dark at 4 °C for 30 minutes. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 0.01 mol/L PBS, pH 7.4. Analyse on a flow cytometer. 																				

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Product-specific limitations

It has been reported that binding of CD34 monoclonal antibodies may be reduced by fixatives present in erythrocyte-lysing reagents, suggesting that the CD34 epitope is negatively affected. Hence, fixatives should be used with caution if the monoclonal antibody is to be used for enumeration purposes (7, 8). Additionally, it has been observed that RPE-Cy5-conjugates may bind to monocytes, resulting in background staining (9).

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
1. Nishio H, Tada J, Hashiyama M, Hirn J, Ingles-Esteve J, Suda T. MC7. CD34 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 974-84.
2. Civin CI, Trischmann TM, Fackler MJ, Bernstein ID, Bühring HJ, Campos L, et al. M7.1. Report on the CD34 cluster workshop. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 818-25.
3. Höffkes H-G, Lowe JA, Pedersen RO, Schmidte G, McDonald DF. BIRMA-K3, a new monoclonal antibody for CD34 immunophenotyping and stem and progenitor cell assay. J Hematother 1996;5:261-70.
4. Civin CI, Strauss LC, Fackler MJ, Trischmann TM, Wiley JM, Loken MR. Positive stem cell selection - basic science. Prog Clin Biol Res 1990;333:387-402.
5. Krause DS, Fackler MJ, Civin CI, May WS. CD34: structure, biology, and clinical utility (review). Blood 1996;87:1-13.
6. Campos L, Guyotat D, Archimbaud E, Devaux Y, Treille D, Laresse A, et al. Surface marker expression in adult acute myeloid leukaemia: correlations with initial characteristics, morphology and response to therapy. Br J Haematol 1989;72:161-6.
7. Macey MG, McCarthy DA, van Agthoven A, Newland AC. How should CD34+ cells be analysed? A study of three classes of antibody and five leucocyte preparation procedures. J Immunol Methods 1997;204:175-88.
8. Serke S, van Lessen A, Huhn D. Quantitative fluorescence flow cytometry: a comparison of the three techniques for direct and indirect immunofluorescence. Cytometry 1998;33:179-87.
9. van Vugt MJ, van den Herik-Oudijk IE, van de Winkel JGJ. Binding of PE-CY5 conjugates to the human high-affinity receptor for IgG (CD64). Blood 1996;88:2358-61.

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Agilent Technologies Singapore (International) Pte Ltd.