Monoclonal Mouse Anti-Human CD5, Clone DK23

Code No./ Code/ Code-Nr. **C 7242 APC-Conjugated**

Code No./ Code/ Code Nr. **F 0795 FITC-Conjugated**

Code No./ Code/ Code Nr. **R 0842 RPE-Conjugated**

ENGLISH													
Intended use	For in vitro diagnostic use. C 7242, F 0795 and R 0842 are intended for use in flow cytometry. In flow cytometry, anti-CD5 is considered essential for the initial evaluation of chronic lymphoproliferative disorders together with a panel of other antibodies. Anti-CD5 may also be useful in immunophenotyping of acute leukaemia (1). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.												
Synonyms for antigen	T1, Tp67 (2).												
Introduction	CD5 is a transmembrane glycoprotein monomer with an Mr of 67 000. It possesses a large cytoplasmic domain suitable for signal transduction (3). CD5 appears early in thymocyte development (3), and is expressed at low density on thymocytes and at high density on all mature T lymphocytes (2). The majority of T-cell malignancies (76%) expresses CD5, and almost 85% of T-cell acute lymphoblastic leukaemias are CD5 positive (4). Further, CD5 is expressed at low density on a small subset of mature B lymphocytes (B1a cells), which is expanded during fetal life and in several autoimmune disorders as well as in some B-cell derived lymphoproliferative disorders, notably chronic lymphocytic leukaemia (CLL) (2) and centrocytic leukaemia (5).												
Reagent provided	The Anti-CD5 conjugates, C 7242, F 0795 and R 0842, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ leucocytes from normal peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> see label on vial.												
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 0795</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> <tr> <td>R 0842</td> <td>RPE (R-Phycoerythrin)</td> <td>X 0928</td> </tr> <tr> <td>C 7242</td> <td>APC (Allophycocyanin)</td> <td>X 0968</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 0795	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927	R 0842	RPE (R-Phycoerythrin)	X 0928	C 7242	APC (Allophycocyanin)	X 0968
Antibody Code No.	Fluorochrome	Negative Control Code No.											
F 0795	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927											
R 0842	RPE (R-Phycoerythrin)	X 0928											
C 7242	APC (Allophycocyanin)	X 0968											
Specificity	The specificity of Anti-CD5, DK23, is equivalent to that of the CD5-clustered antibody Leu-1 (5). In chronic leukaemia and myeloma, Anti-CD5, DK23, has been shown to label 4/4 (100%) T cell chronic lymphocytic leukaemias, 2/2 (100%) T cell polyclonal lymphocytic leukaemias, 36/39 (92%) B cell chronic lymphocytic leukaemias, 5/6 (83%) B cell polyclonal lymphocytic leukaemias, 0/14 (0%) hairy cell leukaemias, and 0/25 (0%) myelomas (6).												
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.												
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.												
Staining procedure	1. Collect venous blood into a test tube containing an anticoagulant. 2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. 3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. 4. Mix 100 µL cell suspension with 10 µL of the fluorochrome-conjugated Anti-CD5. 5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). 6. Incubate in the dark at 4 °C for 30 minutes. 7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. 8. Analyse on a flow cytometer.												

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
1. Braylan RC, Orfao A, Borowitz MJ, Davis BH. Optimal number of reagents required to evaluate hematolymphoid neoplasias: results of an international consensus meeting. *Cytometry* 2001;46:23-7.
2. Lozano F, Simarro M, Calvo J. CD5. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference*; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 1112-3.
3. Lozano F, Calvo J, Roca A, Places L, Simarro M. TC6. CD5 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference*; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 56-8.
4. Leong AS-Y, Cooper K, Leong FJW-M. *Manual of diagnostic antibodies for immunohistology*. London: Oxford University Press; 1999. p. 51-2.
5. Stein H, Lennert K, Feller AC, Mason DY. Immunohistological analysis of human lymphoma: correlation of histological and immunological categories. *Adv Cancer Res* 1984;42:67-147.
6. Erber WN, Mynheer LC, Mason DY. APAAP labelling of blood and bone-marrow samples for phenotyping leukaemia. *Lancet* 1986;i:761-5.

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