Monoclonal Mouse Anti-Human CD61, Platelet Glycoprotein IIIa/APC, Clone Y2/51 Code/ Réf./ Code-Nr. C7280
Monoclonal Mouse Anti-Human CD61, Platelet Glycoprotein IIIa/FITC, Clone Y2/51 Code/ Réf./ Code-Nr. F0803

ENGLISH

Intended use For in vitro diagnostic use.
 C7280 and F0803 are intended for use in flow cytometry. CD61 is a selective marker of platelets and platelet precursors, and antibodies to CD61 may be of value for identification of megakaryoblastic leukaemias (1, 2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Synonym for antigen 9-o-acetyl-GD3 (1).

Introduction CD61 is a 110 kDa single-chain protein and is the common β -subunit of the CD41/CD61 and the CD51/CD61 complexes, both of which are noncovalently associated heterodimers. The complexes require divalent cations for complex integrity and receptor functions (1). The CD41/CD61 complex appears early in megakaryocyte maturation (3). The activated CD41/CD61 complex is a receptor for von Willebrand factor, soluble fibrinogen and fibronectin (1) and plays a central role in platelet activation and aggregation (4). Glanzmann's thrombasthenia is a well-defined but rare inherited disorder of platelet function. It is caused by deficiency or abnormality of the CD41/CD61 complex with bleeding due to defective platelet hemostatic plug formation. Glanzmann's thrombasthenia has been classified into two types: type 1, in which CD41/CD61 is nearly completely absent; and type 2 in which 5-20% of the complex is expressed (5). Clot retraction is absent in type 1 and normal or moderately decreased in type 2. CD51/CD61 complexes play a role in tumour metastasis and adenoviral infection (1).

Reagent provided The Anti-CD61 conjugates, C7280 and F0803, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (10 μ L of conjugate for up to 10⁶ platelets from normal human peripheral blood).
Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code	Fluorochrome	Control ReagentCode
C7280	APC (Allophycocyanin)	X0968
F0803	FITC (Fluorescein Isothiocyanate Isomer 1)	X0927

Immunogen Phytohemagglutinin (PHA)-stimulated peripheral blood cells (6).

Specificity Anti-CD61, Y2/51, was included in the Fourth International Workshop and Conference on Human Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with CD61 (7).
 Anti-CD61 detects platelets in peripheral blood and bone marrow and it also reacts with megakaryocytes (6).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

Storage Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Staining procedure

1. Collect venous blood into a test tube containing EDTA as an anticoagulant.
2. Within 5 minutes centrifuge at 200 x g for 5 minutes at room temperature.
3. Collect 100 μ L of the upper platelet-rich plasma (this is sufficient for 20 tests) and mix into 1 mL of 1% paraformaldehyde in 0.01 mol/L PBS, pH 7.4. Incubate at 4 °C for ½-1 hour.
4. Centrifuge at 1200 x g for 5 minutes at room temperature. Gently aspirate the supernatant and discard it leaving approximately 50 μ L of fluid.
5. Add 2 mL 0.01 mol/L PBS, pH 7.4, and resuspend the platelets by using a vortex mixer.
6. Centrifuge at 1200 x g for 5 minutes at room temperature. Gently aspirate the supernatant and discard it leaving approximately 50 μ L of fluid.
7. Add 1 mL 0.01 mol/L PBS containing 2% fetal calf serum and resuspend the platelets by using a vortex mixer.
8. Mix 50 μ L platelet suspension with 10 μ L fluorochrome-conjugated Anti-CD61.
9. Incubate in the dark at 4 °C for 20 minutes.
10. Wash with PBS. Resuspend the platelets in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
11. Analyse on a flow cytometer.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Procedural notes

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 8: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the isotype and fluorochrome of the conjugated antibody. Recommended control reagents are shown in the table above.

Product-specific limitation C7280 and F0803 will in conjunction with Dako lysing procedures form platelet-leukocyte aggregates (PLA). PLA can be reduced by employing a no-lyse-no-wash procedure (8).

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
1. Sun Q-H, Newman PJ. CD Guide. CD61. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 810.
2. Van Dongen JJM, Adriaansen HJ. Immunobiology of leukemia. In: Henderson ES, Lister TA, Greaves MF, editors. Leukemia. Philadelphia, London, Toronto, Montreal, Sydney, Tokyo: WB Saunders Company; 1996. p. 83-130.
3. Vinci G, Tabilio A, Deschamps JF, van Haeke D, Henri A, Guichard J, et al. Immunological study of *in vitro* maturation of human megakaryocytes. Br J Haematol 1984;56:589-605.
4. Nachman RL, Leung LLK. Complex formation of platelet membrane glycoproteins IIb and IIIa with fibrinogen. J Clin Invest 1982;69:263-9.
5. George JN, Caen JP, Nurden AT. Glanzmann's thrombasthenia: the spectrum of clinical disease. Blood 1990;75:1383-95.
6. Gatter KC, Cordell JL, Turley H, Heryet A, Kieffer N, Anstee DJ, et al. The immunohistological detection of platelets, megakaryocytes and thrombi in routinely processed specimens. Histopathology 1988;13:257-67.
7. von dem Borne AEG, Modderman PW, Admiraal LG, Nieuwenhuis HK. Joint report of the platelet section. P1: Platelet antibodies, the overall results. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p 951-66.
8. Hagberg IA, Lyberg T. Evaluation of circulating platelet-leukocyte conjugates: a sensitive flow cytometric assay well suited for clinical studies. Platelets 2000;11:151-60.

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Agilent Technologies Singapore (International) Pte Ltd.