**Polyclonal Rabbit  
Anti-Human  
Lysozyme EC 3.2.1.17/FITC  
Code No./ Code/ Code-Nr. F 0372**

**ENGLISH**

<b>Intended use</b>	<p>For in vitro diagnostic use.</p> <p>F 0372, Polyclonal Rabbit Anti-Human Lysozyme EC 3.2.1.17/FITC, is intended for use in flow cytometry. F 0372 has been optimized for the labelling of intracellular lysozyme when used in combination with the Dako IntraStain Fixation and Permeabilization Kit, code No. K 2311. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.</p>
<b>Introduction</b>	<p>Lysozyme, originally identified and named by Alexander Fleming in 1922 and defined as a bacteriolytic tissue element, is ubiquitous in its distribution among living organisms and has in humans been demonstrated in several organs and tissue fluids (1). Within the haematopoietic system, lysozyme is known as an intracellular pan-myeloid marker molecule that is selectively expressed by cells of the granulo-monocytic lineage. Normal monocytic precursors seem first to synthesize lysozyme and then later on myeloperoxidase, whereas myeloperoxidase synthesis seems to precede lysozyme in neutrophil development thus mirroring the expression pattern in acute myeloid leukaemia (AML). Mature monocytes are thought to continuously synthesize and secrete lysozyme while granulocytes do not synthesize but only secrete the preformed enzyme (2).</p>
<b>Reagent provided</b>	<p>F 0372 is the purified immunoglobulin fraction of rabbit antiserum conjugated with fluorescein isothiocyanate isomer 1 (FITC). The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10<sup>6</sup> leucocytes from normal human peripheral blood).</p> <p><u>Immunoglobulin concentration g/L:</u> See label on vial.</p> <p><u>F/P ratio:</u> E<sub>495 nm</sub>/E<sub>278 nm</sub> = 1.00 ± 0.15 corresponding to a molar FITC/protein ratio of 4.</p>
<b>Immunogen</b>	<p>Lysozyme isolated from urine of patients with monocytic leukaemia.</p>
<b>Specificity</b>	<p>The antibody reacts with human lysozyme. Traces of contaminating antibodies have been removed by solid-phase absorption with human plasma and urine proteins.</p> <p>The specificity has been ascertained as follows:</p> <p><u>Crossed immunoelectrophoresis:</u> Only the lysozyme precipitation arch appears when using 12.5 µL unconjugated antibody per square cm gel area against 2 µL of concentrated urine from patients with monocytic leukaemia. No precipitation arch is seen when the antibody is tested against 2 µL normal human plasma or 2 µL concentrated normal human urine. Staining: Coomassie Brilliant Blue.</p> <p>In immunocytochemistry the antibody labels the cytoplasm in polymorphs and monocytes but rarely in lymphocytes (3). The antibody has been recommended as part of a panel of antibodies for the immunophenotyping of acute leukaemia in bone marrow biopsy specimens as it aids in the differentiation of AML from acute lymphoblastic leukaemia (ALL) (4). The International Lymphoma Study Group (ILSG) recommends the antibody as part of a panel of antibodies for the immunocytochemical immunophenotyping of histiocytic and dendritic cell neoplasms (5).</p>
<b>Precautions</b>	<ol style="list-style-type: none"> <li>1. For professional users.</li> <li>2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>3. As with any product derived from biological sources, proper handling procedures should be used.</li> </ol>
<b>Storage</b>	<p>Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens.</p> <p>Conjugates should not be stored in diluted form. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.</p>
<b>Staining procedure</b>	<ol style="list-style-type: none"> <li>1. Transfer 50 µL (up to 10<sup>6</sup> cells) of the cell suspension to be analysed (whole blood or bone marrow) to a test tube.</li> <li>2. Add 100 µL Dako IntraStain Reagent A (Fixation), code No. K 2311. Mix gently by using a vortex mixer to ensure that the cells are in suspension.</li> <li>3. Incubate at room temperature for 15 minutes.</li> <li>4. Add 2 mL PBS and mix gently by using a vortex mixer.</li> <li>5. Centrifuge at 300 x g for 5 minutes, then aspirate the supernatant, leaving approximately 50 µL of fluid.</li> <li>6. Mix thoroughly by using a vortex mixer to ensure that the cells are in suspension and add 100 µL Dako IntraStain Reagent B (Permeabilization), code No. K 2311. Add 10 µL of F 0372. Mix gently by using a vortex mixer to ensure that the cells are in suspension.</li> <li>7. Use a non-reactive polyclonal rabbit antibody conjugated with FITC as a negative control.</li> </ol>

8. Incubate in the dark at room temperature for 15 minutes.
9. Repeat steps 4 and 5.
10. Resuspend the pellet in an appropriate fluid for flow cytometric analysis, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
11. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
1. Åström M, Bodin L, Hörnsten P, Wahlin A, Tidefelt U. Evidence for a bimodal relation between serum lysozyme and prognosis in 232 patients with acute myeloid leukaemia. *Eur J Haematol* 2003;70:26-33.
2. Scheinecker C, Strobl H, Fritsch G, Csmarits B, Krieger O, Majdic O, et al. Granulomonocyte-associated lysosomal protein expression during in vitro expansion and differentiation of CD34+ hematopoietic progenitor cells. *Blood* 1995;86:4115-23.
3. Krugliak L, Meyer PR, Taylor CR. The distribution of lysozyme, alpha-1-antitrypsin, and alpha-1-antichymotrypsin in normal hematopoietic cells and in myeloid leukemias: an immunoperoxidase study on cytocentrifuge preparations, smears, and paraffin sections. *Am J Hematol* 1986;21:99-109.
4. Toth B, Wehrmann M, Kaiserling E, Horny H-P. Immunophenotyping of acute lymphoblastic leukaemia in routinely processed bone marrow biopsy specimens. *J Clin Pathol* 1999;52:688-92.
5. Pileri SA, Grogan TM, Harris NL, Banks P, Campo E, Chan JKC, et al. Tumours of histiocytes and accessory dendritic cells: an immunohistochemical approach to classification from the International Lymphoma Study Group based on 61 cases (review). *Histopathol* 2002;41:1-29.

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Agilent Technologies Singapore (International) Pte Ltd.