**Polyclonal Rabbit
Anti-Human
Lactoferrin/FITC
Code No./ Code/ Code-Nr. F 0395
Edition/ Édition/ Ausgabe 01.08.03**

ENGLISH	
Intended use	For in vitro diagnostic use. F 0395, Polyclonal Rabbit Anti-Human Lactoferrin/FITC, is intended for use in flow cytometry. F 0395 has been optimized for the labelling of intracellular lactoferrin when used in combination with the DakoCytomation IntraStain Fixation and Permeabilization Kit, code No. K 2311. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.
Introduction	Lactoferrin is an 80 kDa cationic iron-binding glycoprotein with a proposed immunoregulatory, bactericidal and bacteriostatic activity (1, 2). Free lactoferrin is present in mammalian secretions, such as milk, tears and sperm (1). In the blood, lactoferrin is present in the secondary granules of neutrophil granulocytes and is rarely demonstrated in monocytes and lymphocytes (3). The expression is restricted among haematopoietic cells to the postmitotic maturation compartment of the granulocytic lineage starting from the myelocyte stage. Normal and malignant myeloblasts and promyelocytes are lactoferrin negative (4). In flow cytometry, antibodies to lactoferrin are useful for the subclassification of acute myeloid leukaemia (AML) (4, 5), as lactoferrin is preferentially present in the more mature myeloid subpopulation.
Reagent provided	F 0395 is the purified immunoglobulin fraction of rabbit antiserum conjugated with fluorescein isothiocyanate isomer 1 (FITC). The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Immunoglobulin concentration g/L:</u> See label on vial. <u>F/P ratio:</u> E _{495 nm} /E _{278 nm} = 1.00 ± 0.15 corresponding to a molar FITC/protein ratio of 4.
Immunogen	Lactoferrin isolated from human colostrum.
Specificity	The antibody reacts with human lactoferrin. Traces of contaminating antibodies have been removed by liquid-phase absorption with human plasma proteins. The specificity has been ascertained as follows: <u>Crossed immunoelectrophoresis:</u> Only the lactoferrin precipitation arch appears when using 12.5 µL unconjugated antibody per square cm gel area against 2 µL human milk, and no precipitation arch is seen when using 2 µL human plasma. Staining: Coomassie Brilliant Blue. In immunoelectronmicroscopy, the unconjugated antibody labels the secondary granules of neutrophil granulocytes, whereas primary granules, mitochondria, Golgi apparatus and nuclei are negative (6).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. Conjugates should not be stored in diluted form. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Staining procedure	1. Transfer 50 µL (up to 10 ⁶ cells) of the cell suspension to be analysed (whole blood, bone marrow or mononuclear cells) to a test tube. 2. Add 100 µL DakoCytomation IntraStain Reagent A (Fixation), code No. K 2311. Mix gently by using a vortex mixer to ensure that the cells are in suspension. 3. Incubate at room temperature for 15 minutes. 4. Add 2 mL PBS and mix gently by using a vortex mixer. 5. Centrifuge at 300 x g for 5 minutes, then aspirate the supernatant, leaving approximately 50 µL of fluid. 6. Mix thoroughly by using a vortex mixer to ensure that the cells are in suspension and add 100 µL DakoCytomation IntraStain Reagent B (Permeabilization), code No. K 2311. Add 10 µL of F 0395. Mix gently by using a vortex mixer to ensure that the cells are in suspension. 7. Use a non-reactive polyclonal rabbit antibody conjugated with FITC, as a negative control. 8. Incubate in the dark at room temperature for 15 minutes. 9. Repeat steps 4 and 5.

- Resuspend the pellet in an appropriate fluid for flow cytometric analysis, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
- Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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

- Német K, Simonovits I. The biological role of lactoferrin. Haematologia 1985;18:3-12.
- Brock J. Lactoferrin: a multifunctional immunoregulatory protein? Immunol Today 1995;16:417-9
- Butler TW, Heck LW, Huster WJ, Grossi CE, Barton JC. Assessment of total immunoreactive lactoferrin in hematopoietic cells using flow cytometry. J Immunol Methods 1988;108:159-70.
- Knapp W, Strobl H, Majdic O. Flow cytometric analysis of cell-surface and intracellular antigens in leukemia diagnosis. Cytometry 1994;18:187-98.
- Koniková E, Glasová M, Kusenda J, Babušíková O. Intracellular markers in acute myeloid leukemia diagnosis. Neoplasma 1998;45:282-91.
- Esaguy N, Águas AP, Silva MT. High-resolution localization of lactoferrin in human neutrophils: labeling of secondary granules and cell heterogeneity. J Leukoc Biol 1989;46:51-62.

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