**Polyclonal Goat
Anti-Mouse Immunoglobulins/FITC
Goat F(ab')₂**

Code F0479

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>F0479 is intended for use in flow cytometry.</p> <p>F0479 has been produced to give optimal performance in flow cytometry when used in combination with a monoclonal mouse antibody. Being the F(ab')₂ fragment of goat immunoglobulin, the reagent is particularly well-suited for work on unfixed cells exhibiting Fc receptors, and also for other applications where the Fc part of the antibody molecule could disturb. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.</p>
Summary and explanation	<p>Immunofluorescence (IF) is a common laboratory technique used in almost all aspects of biology. This technique, has been widely used both in research and clinical diagnostics. Applications include the evaluation of cells in suspension, cultured cells, frozen tissue, FFPE tissue, beads, and microarrays for the detection of specific proteins. In the IF techniques, antibodies are chemically conjugated to fluorescent dyes such as fluorescein isothiocyanate (FITC). These labeled antibodies can be use directly or indirectly to bind to the antigen of interest, which allows for antigen detection through fluorescence techniques. The degree of fluorescence can then be detected using a flow cytometer, array scanner, or automated imaging instrument, or visualized using fluorescence or confocal microscopy (1).</p> <p>F0479 is used in the indirect staining method as a secondary antibody for example in the testing of mouse hybridoma cell cultures for antibody production.</p>
Reagent provided	<p>F0479 is the purified and F(ab')₂ fragmented immunoglobulin fraction of goat antiserum conjugated with fluorescein isothiocyanate isomer 1 (FITC). The conjugate is provided in liquid form in 0.05 mol/L Tris-HCl, 15 mmol/L NaN₃, 1% bovine serum albumin, pH 7.2.</p> <p><u>Protein concentration g/L:</u> See label on vial.</p>
Preparation	<ol style="list-style-type: none"> 1. The immunoglobulin fraction has been solid-phase absorbed to remove antibodies cross-reacting with human immunoglobulins and fetal calf serum. 2. The absorbed immunoglobulin fraction has been further purified by affinity chromatography using agarose beads coupled with mouse immunoglobulins. 3. The affinity-isolated immunoglobulin fraction has then been degraded by pepsin and the F(ab')₂ fraction isolated by gel filtration. 4. Finally, the F(ab')₂ fraction is conjugated with fluorescein isothiocyanate isomer 1.
Immunogen	<p>Immunoglobulins, mainly IgG, isolated from mouse serum.</p>
Specificity	<p>The antibody reacts with all mouse IgG subclasses and mouse IgM. Reaction with other mouse immunoglobulins has not been investigated, but the antibody will probably react with these via their light chains.</p> <p>Cross-reaction with human immunoglobulins and fetal calf serum is very low as determined by ELISA, less than 0.5%. The cross-reaction with immunoglobulins of rabbit and swine is less than 0.5%.</p> <p>The cross-reaction with immunoglobulins of guinea pig is about 10%, and with rat immunoglobulins about 30%.</p>
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (NaN₃) a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	<p>Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user.</p> <p>Conjugates should not be stored in diluted form. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.</p>
Dilution guidelines	<p>In flow cytometry, F0479 may be used at a volume of 10 μL per test and in a dilution of 1:10-1:20. This is a guideline only. The optimal antibody concentration may vary depending on specimen and preparation method and should be validated by each individual laboratory. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.</p>

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
References/ Bibliographie/ Literaturnachweise

1. Taylor CR and Rudbeck L. Education Guide: Immunohistochemical Staining Methods. Sixth Edition. Dako, Carpinteria, California; 2013.

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Agilent Technologies Singapore (International) Pte Ltd.