**Monoclonal Mouse
Anti-Human
CD11c, Protein 150,95/FITC**
Clone KB90
Code No./ Code/ Code-Nr. F 0713

ENGLISH							
Intended use	For in vitro diagnostic use. F 0713 is intended for use in flow cytometry. Anti-CD11c is considered essential for the initial evaluation of B-cell lymphoproliferative disorders together with a panel of other antibodies (1-3). Interpretation of staining results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.						
Synonym for antigen	Leucocyte surface antigen p150,95, α X integrin chain, $\alpha^x\beta_2$, and CR4 (4).						
Introduction	The CD11c antigen is the 150 kDa α -chain of the leucocyte integrin complex p150,95. The protein belongs to a family of related heterodimeric proteins designated CD11 in the classification system for human leucocyte differentiation antigens. The CD11 family shares a common 95 kDa β -chain (CD18). In addition to the p150,95 protein, this protein family contains two other members, leucocyte function associated-1 antigen (LFA-1, CD11a) and iC3b receptor (CD11b) (4). On circulating leucocytes CD11c is expressed strongly on monocytes and weakly on neutrophil granulocytes (4). NK cells, activated B and T cells also express CD11c (5). CD11c is expressed on hairy cell leukaemia (HCL) cells. The antigen is also expressed on cells from cases of B-cell chronic lymphocytic leukaemia (B-CLL), B-cell prolymphocytic leukaemia (B-PLL), splenic marginal zone lymphoma, acute myeloblastic leukaemia (AML), acute myelomonocytic leukaemia (AMML) and acute megakaryoblastic leukaemia (AMKL) (1-3, 5, 6).						
Reagent provided	The Anti-CD11c conjugate, F 0713, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 μ L of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.						
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 0713</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 0713	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927
Antibody Code No.	Fluorochrome	Negative Control Code No.					
F 0713	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927					
Immunogen	Hairy cell leukaemia cells (7).						
Specificity	Anti-CD11c, KB90, was included in the Sixth International Workshop and Conference on Leucocyte Cell Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with the CD11c antigen (5). Anti-CD11c, KB90, is reactive to the C-terminal region of CD11c (8). Anti-CD11c, KB90, reacts with a small population of lymphocytes, strongly with monocytes and weakly with granulocytes in peripheral blood (7). In chronic lymphoproliferative disorders, Anti-CD11c, KB90, has been shown to label HCL (21/21 cases (3) and 14/14 cases (6)), B-CLL (52/106 cases (3) and 8/39 cases (6)), B-PLL (9/9 cases (3) and 2/6 cases (7)) and CD5-negative non-Hodgkin's lymphoma (11/19 cases (3)). In acute leukaemias, Anti-CD11c, KB90, has been shown to label AMML (36/36 cases), AMKL (4/9 cases) and AML (16/30 cases) (6).						
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. 						
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.						
Staining procedure	<ol style="list-style-type: none"> Collect venous blood into a test tube containing an anticoagulant. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. In one test tube, mix 100 μL cell suspension with 10 μL F 0713. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). Incubate in the dark at 4 °C for 30 minutes. Wash twice with PBS containing 2% BSA. Resuspend the cells in each test tube in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. 						

8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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

- van Dongen JJM, Adriaansen HJ. Immunobiology of leukemia. In: Henderson ES, Lister TA, Greaves MF, editors. Leukemia. Philadelphia, London, Toronto, Montreal, Sydney, Tokyo: WB Saunders Company; 1996. p. 83-130.
- Sánchez ML, Almeida J, Vidriales B, López-Berges MC, García-Marcos MA, Moro MJ, et al. Incidence of phenotypic aberrations in a series of 467 patients with B chronic lymphoproliferative disorders: basis for the design of specific four-color stainings to be used for minimal residual disease investigation. Leukemia 2002;16:1460-9.
- Marotta G, Raspadori D, Sestigiani C, Scalia G, Bigazzi C, Lauria F. Expression of the CD11c antigen in B-cell chronic lymphoproliferative disorders. Leuk Lymphoma 2000;37:145-9.
- Hogg N. AS4. CD11c workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 347-348.
- Luk J, Springer TA. AS5.3. CD11c cluster report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 2. p. 1590-1591.
- Erber WN, Mynheer LC, Mason DY. APAAP labelling of blood and bone-marrow samples for phenotyping leukaemia. Lancet 1986;i:761-5.
- MacDonald SM, Pulford K, Falini B, Micklem K, Mason DY. A monoclonal antibody recognizing the p150/95 leucocyte differentiation antigen. Immunology 1986;59:427-31.
- Luk J, Luther E, Diamond MS, Springer TA. AS5.9. Subunit specificity and epitope mapping of Mac-1 and p150,95 mAb using chimeric CD11b X CD11c transfectants. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 2. p. 1599-1601.

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