Monoclonal Mouse Anti-Human CD2, Clone MT910

Code No./ Code/ Code-Nr. **F 0767 FITC-Conjugated**
Code No./ Code/ Code-Nr. **R 0807 RPE-Conjugated**

ENGLISH

Intended use For in vitro diagnostic use.

F 0767 and R 0807 are intended for use in flow cytometry. CD2 is present on virtually all normal T lymphocytes, and it can be considered a pan-T cell antigen. CD2 is a useful marker in the assessment of lymphoid malignancies as it is expressed in the majority of precursor and postthymic lymphomas and leukaemias. In some neoplastic T-cell populations, e.g. in peripheral T-cell lymphomas, CD2 may be aberrantly deleted (1). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Synonyms for antigen CD2R, E-rosette receptor, T11, lymphocyte-function antigen-2 (LFA-2) (2).

Introduction CD2 is a 50 kDa single-chain type I transmembrane glycoprotein, originally defined as the sheep erythrocyte receptor. CD2 primarily binds CD58 on antigen-presenting cells and induces costimulatory signals in T cells (2). CD2 is one of the earliest T-cell lineage-restricted antigens to appear during T-cell differentiation, and it is expressed on most thymocytes, peripheral T cells and NK cells. Some B-cell non-Hodgkin lymphomas show aberrant expression of CD2. This may represent clonal expansion of CD2+ B lymphocytes that normally constitute a small fraction of peripheral B lymphocytes (3). Additionally, aberrant CD2 expression has been found in some acute myeloid leukaemias almost exclusively restricted to the M3 and M4Eo variants (4).

Reagent provided The Anti-CD2 conjugates, F 0767 and R 0807, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).
Isotype: IgG1, kappa. Conjugate concentration mg/L: see label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
F 0767	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927
R 0807	RPE (R-Phycoerythrin)	X 0928

Specificity Anti-CD2, MT910, was included in the Second, Third and Fourth International Workshops and Conferences on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD2 (5).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

Storage Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Staining procedure

1. Collect venous blood into a test tube containing an anticoagulant.
2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
4. Mix 100 µL cell suspension with 10 µL of the fluorochrome-conjugated Anti-CD2.
5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
6. Incubate in the dark at 4 °C for 30 minutes.
7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
1. Leong A, Cooper K, Leong F. Manual of diagnostic antibodies for immunohistology; London: Oxford University Press; 1999. p. 45-6.
2. Kato K. CD Guide CD2. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 748-9.
3. Kaleem Z, White G, Zutter MM. Aberrant expression of T-cell-associated antigens on B-cell non-Hodgkin lymphomas. Am J Clin Pathol 2001;115:396-403.
4. Thalhammer-Scherrer R, Mitterbauer G, Simonitsch I, Jaeger U, Lechner K, Schneider B, et al. The immunophenotype of 325 adult acute leukemias: relationship to morphologic and molecular classification and proposal for a minimal screening program highly predictive for lineage discrimination. Am J Clin Pathol 2002;117:380-9.
5. Bernard A, Brown MH, Yang SY, Wallace DL. T2.0. Summary of the CD2 workshop. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sept 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 106.

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