## Monoclonal Mouse Anti-Human Ki-67 Antigen, Clone Ki-67

**Code No./ Code/ Code-Nr. F 0788**      **FITC-Conjugated**  
**Code No./ Code/ Code-Nr. R 0840**      **RPE-Conjugated**

### ENGLISH

#### Intended use

For in vitro diagnostic use.

F 0788 and R 0840 are intended for use in flow cytometry. Anti-Ki-67 antigen is useful for assessing cellular proliferation in acute leukaemia, lymphoproliferative disorders and breast cancer (1-3). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

#### Introduction

The Ki-67 antigen is a nuclear protein defined by its reactivity with the monoclonal antibody from the Ki-67 clone (4, 5). During interphase, the antigen can be exclusively detected within the nucleus, whereas in mitosis most of the protein is relocated to the surface of the chromosomes (4). The complete gene locus of the Ki-67 has been sequenced. The size of the gene is approximately 30,000 base pairs organized in 15 exons with sizes from 67 to 6845 base pairs and in 14 introns with sizes from 87 to 3569 base pairs. The gene is located on chromosome 10 (6).

The Ki-67 protein is expressed in all proliferating cells during late G<sub>1</sub>, S, M and G<sub>2</sub> phases of the cell cycle while cells in the G<sub>0</sub> (non-cycling cells) phase consistently lack the Ki-67 antigen. Un-stimulated normal human cells do not express the Ki-67 antigen (4, 5).

Numerous studies have shown the use of anti-Ki-67 antigen in measuring the Ki-67 labelling index or growth fraction in various solid tumours using immunohistochemical methods (4).

Flow cytometric analysis has demonstrated to be a useful method for detecting Ki-67 antigen and assessing cellular proliferation in tumor cells as an alternative to S-phase cell cycle determination, continuous 3H-thymidine labelling, acridine orange and bromodeoxyuridine staining (7, 8).

#### Reagent provided

The Anti-Ki-67 Antigen conjugates, F 0788 and R 0840, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10<sup>6</sup> RAJI Cells).

Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
F 0788	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927
R 0840	RPE (R-Phycoerythrin)	X 0928

#### Immunogen

Crude nuclear fraction of L428 cells (4).

#### Specificity

Immunohistochemical and flow cytometric analysis has shown that Anti-Ki-67 Antigen reacts negatively or weakly with un-stimulated cells from normal healthy donors (1-3, 5).

Flow cytometric analysis of breast cancer samples demonstrated that Anti-Ki-67 Antigen labelled proliferating tumor cells (2, 3). One study with 438 breast cancer cases showed variable cellular Ki-67 expression levels (1 – 60% above normal level) (2). Infiltrating lobular carcinomas showed lower Ki-67 antigen expression levels compared to ductal carcinomas. Medulary carcinomas had the highest Ki-67 antigen levels. Another study with 154 breast cancer cases showed overall (all cell cycle phases: G0+G1, S+G2M) and G1 cellular Ki-67 antigen expression levels from 5 – 100% and 2 - 95%, respectively (3). In both studies there was a positive correlation between S-phase or S+G<sub>2</sub>M-phase fraction and Ki-67 staining (2, 3).

Anti-Ki-67 Antigen also reacts with proliferating neoplastic cells in acute lymphoblastic leukaemia (ALL), acute myeloid leukaemia (AML), chronic myeloid leukaemia (CML), chronic myeloid leukaemia in blast crisis (CML-BC), B-cell chronic lymphocytic leukaemia (B-CLL), hairy cell leukaemia (HCL), prolymphocytic leukaemia (PLL) and immunocytoma (IC). Results show high cellular Ki-67 antigen levels in ALL (10 - 65%), IC (10 - 44%), B-CLL (up to 20%) and CML-BC cases (10 – 35%). Prolymphocytic leukaemia and hairy cell leukaemia have lower Ki-67 antigen levels (1).

#### Precautions

- For professional users.
- This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
- As with any product derived from biological sources, proper handling procedures should be used.

#### Storage

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

#### Staining procedure

- Harvest cells and determine total number present. Wash twice in wash buffer. (If required, this procedure may also be used with whole blood samples)
- If required perform staining of cell surface antigens using appropriate directly conjugated monoclonal antibodies at this stage. Following staining wash cells once in PBS and discard supernatant. N.B. Phycoerythrin conjugates are not suitable for the detection of cell surface antigens using this method.
- Add DakoCytomation IntraStain Reagent A (Fixation), code No. K 2311, using 50µl per 1x10<sup>6</sup> cells. Mix gently by using a vortex mixer to ensure that the cells are in suspension.
- Incubate at room temperature for 10 minutes.
- Add 500µl cold (-20°C) absolute methanol. Mix gently by using a vortex mixer and incubate for 10 minutes at room temperature.
- Add 2 mL PBS and mix gently by using a vortex mixer.
- Centrifuge at 300 x g for 5 minutes, then aspirate the supernatant, leaving approximately 50 µL of fluid..

(102453-002)

F 0788/R 0840/EFG/SSA/01.03.03 p. 1/4

- Mix thoroughly by using a vortex mixer to ensure that the cells are in suspension and add DakoCytomation IntraStain Reagent B (Permeabilization), code No. K 2311, using 50µl per 1 x 10<sup>6</sup> cells.
- Add 10 µl of the anti-Ki-67 Antigen conjugate. Mix gently by using a vortex mixer to ensure that the cells are in suspension.
- Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
- Incubate in the dark at 4 °C for 30 minutes.
- Repeat steps 6 and 7.
- Resuspend in 0.25 mL 0.5% paraformaldehyde (fixative) in 0,01 mol/L PBS, pH 7,4.PBS.
- Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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