Monoclonal Mouse Anti-Human CD7, Clone DK24

Code No./ Code/ Code-Nr. **F 0789** FITC-Conjugated
Code No./ Code/ Code-Nr. **R 7208** RPE-Conjugated

ENGLISH											
Intended use	For in vitro diagnostic use. F 0789 and R 7208 are intended for use in flow cytometry. Anti-CD7 is considered essential for the initial evaluation of T cell acute lymphoblastic leukemias (T-ALL) and T cell chronic leukemias together with a panel of other antibodies (1). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.										
Synonym for antigen	gp40 (2).										
Introduction	The CD7 antigen is a membrane-bound 40 kDa glycoprotein expressed on thymocytes, mature T cells, a large majority of natural killer cells, pluripotent hematopoietic stem cells, and progenitor cells of lymphoid and myeloid cells (3, 4, 5). CD7 is the earliest T-cell specific antigen to be expressed in lymphocytes, and the only early marker to persist throughout differentiation (4). CD7 is expressed on immature, common thymocytic and mature T-ALL, T cell prolymphocytic leukemia, adult T-cell leukemia, cutaneous T-cell leukemia/lymphoma and large granular leukaemias (1, 6). CD7 is also expressed on special cases of AML of hematopoietic stem cell origin (7). CD7-deficient T cells have been observed in severe combined immunodeficiency (SCID) and rheumatoid arthritis (2, 6). Although the exact function of CD7 is not known, cross-linking studies with monoclonal antibodies suggest that CD7 can function as a costimulatory molecule, induce cytokine secretion, and modify cellular adhesion (2).										
Reagent provided	The Anti-CD7 conjugates, F 0789 and R 7208, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG2b, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.										
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 0789</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0941</td> </tr> <tr> <td>R 7208</td> <td>RPE (R-Phycoerythrin)</td> <td>X 0951</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 0789	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0941	R 7208	RPE (R-Phycoerythrin)	X 0951	
Antibody Code No.	Fluorochrome	Negative Control Code No.									
F 0789	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0941									
R 7208	RPE (R-Phycoerythrin)	X 0951									
Specificity	Anti-CD7, DK24, labels the majority of T cells present in normal peripheral blood and bone marrow.										
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. 										
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.										
Staining procedure	<ol style="list-style-type: none"> Collect venous blood into a test tube containing an anticoagulant. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD7. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). Incubate in the dark at 4 °C for 30 minutes. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. Analyse on a flow cytometer. <p>It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.</p>										

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