Monoclonal Mouse Anti-Human CD45R0, Clone UCHL1

Code No./ Code/ Code-Nr. F 0800 **FITC-Conjugated**
Code No./ Code/ Code-Nr. R 0843 **RPE-Conjugated**

ENGLISH										
Intended use	For in vitro diagnostic use. F 0800 and R 0843 are intended for use in flow cytometry. Anti-CD45R0, UCHL1, labels both normal and neoplastic cells (1, 2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.									
Introduction	CD45 is a transmembrane glycoprotein expressed on all nucleated cells of haematopoietic origin. CD45, encoded by a single gene mapped to chromosome 1, has various isoforms, based on differential splicing of exons 4, 5 and 6. Five different isoforms of CD45, named ABC, AB, BC, B and 0, have been identified on human leucocytes. CD45RC, CD45RB, CD45RA and CD45R0 antibodies recognize these isoforms. The Mr of the isoforms ranges from 220 000 for the ABC isoform to 180 000 for the 0 isoform (3). All the CD45 isoforms share the same intracellular segment, which has been shown to have tyrosine phosphatase activity. Various leucocytes express characteristic CD45 isoforms, and T cells differentially express CD45 isoforms at various stages of their development and activation (3). Most B-lymphocytes and leukaemic B-cell lines predominantly express the CD45RA isoform and lack the CD45R0 isoform, whereas thymocytes and many leukaemic T-cell lines primarily express CD45R0 and CD45RA (4). Peripheral T cells can roughly be divided into two populations: CD45RA-, CD45R0+ "memory" T cells and CD45+, CD45R0- "naïve" T cells. Memory T cells respond to recall antigen and possess B-cell helper function, while naïve T cells respond poorly to recall antigen and lack B cell helper function (4).									
Reagent provided	The Anti-CD45 conjugates, F 0800 and R 0843, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). Isotype: IgG2a, kappa. Conjugate concentration mg/L: see label on vial.									
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 0800</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0933</td> </tr> <tr> <td>R 0843</td> <td>RPE (R-Phycoerythrin)</td> <td>X 0950</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 0800	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0933	R 0843	RPE (R-Phycoerythrin)	X 0950
Antibody Code No.	Fluorochrome	Negative Control Code No.								
F 0800	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0933								
R 0843	RPE (R-Phycoerythrin)	X 0950								
Immunogen	IL-2 dependent T-cell line, CA1 (5).									
Specificity	Anti-CD45R0, UCHL1, was clustered as anti-CD45R0 at the Fourth (6) and Fifth (4) International Workshops and Conferences on Human Leucocyte Differentiation Antigens. In Western blotting of lysate of peripheral blood mononuclear cells, the antibody was found to react specifically with a 180 kDa band corresponding to CD45R0. In leucocyte common antigen (LCA) transfected cell lines, expressing the ABC, the AB, the BC, the B or the 0 CD45 isoform only, the antibody reacted solely with the transfectant expressing CD45R0 (7). Anti-CD45R0, UCHL1, labels most thymocytes, a subpopulation of resting T cells within both the CD4 and CD8 subsets, and mature activated T cells. This antibody also labels cells of the myeloid lineage, whereas most normal B cells and NK cells are consistently negative (5). Weak cytoplasmic staining is, however, seen in cases of centroblastic and immunoblastic lymphoma (8). Anti-CD45R0, UCHL1, recognizes approximately 50% of human T cells in peripheral blood. The antibody is of interest to immunologists studying T cell activation since it recognizes a functionally distinct subset of human helper-inducer T cells. Unprimed T cells, being largely CD45RA-positive and UCHL1-negative, lose the CD45RA antigen and concomitantly acquire CD45R0 on activation by PHA or alloantigen. The CD45R0-positive T cells are most probably a primed population within which memory cells reside (4, 9).									
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. 									
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.									
Staining procedure	<ol style="list-style-type: none"> Collect venous blood into a test tube containing an anticoagulant. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. Mix 100 µL cell suspension with 10 µL of the fluorochrome-conjugated Anti-CD45R0. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). Incubate in the dark at 4 °C for 30 minutes. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. Analyse on a flow cytometer. <p>It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.</p>									

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
1. Norton AJ, Ramsay AD, Smith SH, Beverley PCL, Isaacson PG. Monoclonal antibody (UCHL1) that recognises normal and neoplastic T cells in routinely fixed tissues. J Clin Pathol 1986;39:399-405.
2. Davey FR, Elghetany MT, Kurec AS. Immunophenotyping of hematologic neoplasms in paraffin-embedded tissue sections. Am J Clin Pathol 1990;93 Suppl 1:S17-S26.
3. Sewell WA, Cooley MA, Hegen M. NL6. CD45 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 499-502.
4. Morimoto C. T18. CD45 cluster report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 386-9.
5. Smith SH, Brown MH, Rowe D, Callard RE, Beverley PCL. Functional subsets of human helper-inducer cells defined by a new monoclonal antibody, UCHL1. Immunology 1986;58:63-70.
6. Schmidt RE. Non-lineage/natural killer section report: new and previously defined clusters. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 517-42.
7. Poppema S, Lai R, Visser L. Monoclonal antibody OPD4 is reactive with CD45R0, but differs from UCHL1 by the absence of monocyte reactivity. Am J Pathol 1991;139:725-9.
8. Hall PA, D'Ardenne AJ, Butler MG, Habeshaw JR, Stansfeld AG. New marker of B lymphocytes, MB2: comparison with other lymphocyte subset markers in conventionally processed tissue sections. J Clin Pathol 1987;40:151-6.
9. Akbar AN, Terry L, Timms A, Beverley PCL, Janossy G. Loss of CD45R and gain of UCHL1 reactivity is a feature of primed T cells. J Immunol 1988;140:2171-8.

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