Monoclonal Mouse Anti-Human CD25, Interleukin-2 Receptor, Clone ACT-1

Code No./ Code/ Code-Nr. **F 0801 FITC-Conjugated**
Code No./ Code/ Code-Nr. **R 0811 RPE-Conjugated**

ENGLISH										
Intended use	For in vitro diagnostic use. F 0801 and R 0811 are intended for use in flow cytometry. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.									
Synonyms for antigen	Interleukin-2 receptor (IL-2R), interleukin-2 receptor α chain, Tac antigen (1, 2, 3).									
Introduction	The interleukin-2 receptors consist of three polypeptides (α , β , γ), which in different combinations bind IL-2 (1). The 55 kDa α chain is called CD25 (3). The α chain alone, the $\beta\gamma$ heterodimer, and the $\alpha\beta\gamma$ heterotrimer form the low-, intermediate- and high affinity IL-2R, respectively. The intermediate- and high affinity IL-2R play a role in intracellular signal transduction (1). CD25 is strongly expressed on the surface of T cells activated by antigens, mitogens (phytohaemagglutinin A (PHA), concanavalin A (Con-A), pokeweed mitogen (PWM)), <i>Staphylococcus aureus</i> , viruses (HTLV-I, HTLV-II, Epstein-Barr virus) (4, 5), and is also expressed on B cells stimulated with anti-IgM antibody, and on monocytes/macrophages stimulated with lipopolysaccharide (6). In normal peripheral blood a variable staining of lymphoid cells (about 3-25%) is observed.									
Reagent provided	The Anti-CD25 conjugates, F 0801 and R 0811, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 μ L of conjugate for up to 10 ⁶ PHA-stimulated lymphocytes from normal peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> see label on vial.									
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 0801</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> <tr> <td>R 0811</td> <td>RPE (R-Phycoerythrin)</td> <td>X 0928</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 0801	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927	R 0811	RPE (R-Phycoerythrin)	X 0928
Antibody Code No.	Fluorochrome	Negative Control Code No.								
F 0801	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927								
R 0811	RPE (R-Phycoerythrin)	X 0928								
Immunogen	Phytohemagglutinin-stimulated peripheral blood lymphocytes.									
Specificity	The specificity of Anti-CD25, ACT-1, is identical to that of CD25-clustered antibodies as indicated by analysis of immunoprecipitates and the immunocytochemical labelling pattern (7, 8).									
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. 									
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.									
Staining procedure	<ol style="list-style-type: none"> Collect venous blood into a test tube containing an anticoagulant. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. Mix 100 μL cell suspension with 10 μL fluorochrome-conjugated Anti-CD25. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome as a negative control (see table). Incubate in the dark at 4 °C for 30 minutes. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. Analyse on a flow cytometer. <p>It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.</p>									

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
Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p.1855-64.

- Uchiyama T, Broder S, Waldmann TA. A monoclonal antibody (anti-Tac) reactive with activated and functionally mature human T cells. I. Production of anti-Tac monoclonal antibody and distribution of Tac (+) cells. J Immunol 1981;126:1393-7.
- Uchiyama T, Nelson DL, Fleisher TA, Waldmann TA. A monoclonal antibody (anti-Tac) reactive with activated and functionally mature T cells. II. Expression of Tac antigen on activated cytotoxic killer T cells, suppressor cells, and on one of two types of helper T cells. J Immunol 1981;126:1398-1403.
- Sugamura K. CD25. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 1127-8.
- Schwartz R, Gerdes J, Ziegler A, Stein H. Immunoprecipitation of the interleukin-2-receptor from Hodgkin's disease derived cell lines by monoclonal antibodies. Hematol Oncol 1987;5:57-64.
- Aziz KE, McCluskey PJ, Wakefield D. Expression of selectins (CD62 E,L,P) and cellular adhesion molecules in primary Sjögren's syndrome: questions to immunoregulation. Clin Immunol Immunopathol 1996;80:55-66.

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- Ishii N, Kondo M, Takeshita T, Sugamura K. CR2.1. mAb specific for the γ chain of the IL-2 receptor. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leukocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p.1867-8.
- Leonard WJ, Depper JM, Robb RJ, Waldmann TA, Greene WC. Characterization of the human receptor for T-cell growth factor. Proc Natl Acad Sci (USA) 1983;80:6957-61.
- Kikutani H, Kishimoto T. CR1. The cytokine receptors: section report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leukocyte typing V. White cell differentiation antigens.