**Monoclonal Mouse  
Anti-Human  
HLA-DP, DQ, DR Antigen/FITC**  
Clone CR3/43  
**Code No./ Code/ Code-Nr. F0817**

ENGLISH							
<b>Intended use</b>	For in vitro diagnostic use. F0817 is intended for use in flow cytometry. Anti-HLA-DP, DQ, DR antigen is considered essential for the initial evaluation of acute leukaemias, chronic T and myeloid leukaemias together with a panel of other antibodies (1-3). F0817 is not intended for use in tissue typing. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.						
<b>Synonym for antigen</b>	HLA-DR (5-7).						
<b>Introduction</b>	The human leucocyte antigen (HLA) system, originally discovered as the result of a transfusion reaction, is now known to play a crucial role in many areas of clinical medicine. The HLA molecules are encoded by a cluster of tightly linked genes located on the short arm of chromosome 6. Based on some of the structural and functional characteristics of the genes, the region has been divided into three: HLA class I, Class II and class III regions (4). The A and B genes of the HLA class II, DP, DQ and DR encode a heterodimer formed by two non-covalently associated $\alpha$ and $\beta$ chains of approximately 34 and 28 kDa respectively (4). The main function of the HLA-DP, DQ and DR molecules is to present antigenic peptides, mostly of exogenous nature, to CD4+ T-cells. HLA molecules are also known to be associated with a variety of autoimmune, non-autoimmune and infectious diseases and to restrict the antibody response to certain antigens and vaccines (4). HLA-DP, DQ and DR molecules are constitutively expressed on antigen-presenting cells (APC) such as B lymphocytes, monocytes and dendritic cells but can also be detected on cytotoxic/suppressor T lymphocytes and activated granulocytes (4, 5). It is uncertain whether HLA-DP, DQ and DR antigens are also expressed on activated platelets. HLA class II expression can also be induced on cells and tissues such as fibroblasts and endothelial cells as a result of activation and/or by certain cytokines such as $\gamma$ -interferon, tumor necrosis factor and interleukin-10 (4). The antigen has been found on the cell surface of leukaemic blasts from cases of B-cell acute lymphoblastic leukaemia (ALL), T-cell pre-ALL, acute myeloid leukaemia (AML) except AML-M3, and chronic B and T cell leukaemia, chronic myeloid leukaemia (CML) in blast crisis and lymphomas of B cell and T cell type (1-3, 6, 7). HLA-DP, DQ, DR antigen is normally not present on non-haematopoietic tumors and multiple myeloma (6).						
<b>Reagent provided</b>	The Anti-HLA-DP, DQ, DR Antigen conjugate, F 0817, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2. Each vial contains 100 tests (10 $\mu$ L of conjugate for up to 10 <sup>5</sup> leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG1, Kappa. <u>Coniugate concentration mg/L:</u> See label on vial.						
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F0817</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X0927</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F0817	FITC (Fluorescein Isothiocyanate Isomer 1)	X0927
Antibody Code No.	Fluorochrome	Negative Control Code No.					
F0817	FITC (Fluorescein Isothiocyanate Isomer 1)	X0927					
<b>Specificity</b>	Anti-Human HLA-DP, DQ, DR Antigen, CR3/43, reacts with the $\beta$ -chain of the $\alpha$ $\beta$ heterodimer of all products of the gene families DP, DQ and DR. The antibody was included in the First International Workshop and Conference on Monoclonal Antibodies to Human MHC Class II Antigens (1983) and its specificity and other characteristics were ascertained by a variety of techniques, including reactivity with isolated antigen, immunoblotting, and labelling of transfected cells (8). In normal peripheral blood the antibody stains B cells and most monocytes but is unreactive with normal T cells and polymorphs. It will, however, stain activated T cells in peripheral blood (5). Anti-HLA-DP, DQ, DR Antigen does not react with erythrocytes and megakaryocytes (6). Immunohistochemical analysis demonstrated that Anti-HLA-DP, DQ, DR, Antigen, CR3/43 labels AML (5/5 cases (6)), B cell ALL (3/3 cases (6)), chronic leukaemias and lymphomas of B and T cell type (3/3 cases (6) and 45/46 cases (7)) and CML in myeloid blast crisis (1/1 case (6)). The antibody does not label multiple myeloma (0/3 cases) but shows weak staining of a minority of cells in metastatic breast carcinomas (2/5 cases) (6).						
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.						
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.						
<b>Staining procedure</b>	<ol style="list-style-type: none"> <li>Collect venous blood into a test tube containing an anticoagulant.</li> <li>Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.</li> <li>Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.</li> <li>Mix 100 <math>\mu</math>L cell suspension with 10 <math>\mu</math>L F0817.</li> <li>Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).</li> <li>Incubate in the dark at 4 °C for 30 minutes.</li> <li>Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.</li> </ol>						
(102455-003)	F0817/EFG/MBH/04.03.05 p. 1/4						
DakoCytomation Denmark A/S · Produktionsvej 42 · DK-2600 Glostrup · Denmark · Tel. +45 44 85 95 00 · Fax +45 44 85 95 95 · CVR No. 33 21 13 17							

8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
- Van Dongen JJM, Adriaansen HJ. Chapter 6. Immunobiology of leukaemia. In: Henderson ES, Lister TA, Greaves MF. Leukemia. Philadelphia, London, Toronto, Montreal, Sydney, Tokyo: WB Saunders Company; 1996. p. 83-130.
- Mittelman M, Karcher DS, Kammerman LA, Lessin LS. High Ia (HLA-DR) and low CD11b (Mo1) expression may predict early conversion to leukemia in myelodysplastic syndromes. Am J Hematol 1993;43:165-71.
- De Zen L, Orfao A, Cazzaniga G, Masiero L, Cocito MG, Spinelli M, et al. Quantitative multiparametric immunophenotyping in acute lymphoblastic leukemia: correlation with specific genotype. I. ETV6/AML1 ALLs identification. Leukemia 2000;14:1225-31.
- Naverrete CV. The HLA system in blood transfusion. Baillière's Clin Haematol 2000;13:511-32.
- Erber WN, Pinching AJ, Mason DY. Immunocytochemical detection of T and B cell populations in routine blood smears. Lancet 1984;1:1042-5.
- Falini B, Martelli F, Tarallo F, Moir DJ, Cordell JL, Gatter KC, et al. Immunohistological analysis of human bone marrow trephine biopsies using monoclonal antibodies. Br J Haematol 1984;56:365-86.
- Smith MEF, Holgate CS, Williamson JMS, Grigor I, Quirke P, Bird CC. Major histocompatibility complex class II antigen expression in B and T cell non-Hodgkin's lymphoma. J Clin Pathol 1987;40:34-41.
- Appendix. Table of workshop monoclonal antibodies and a synopsis of their principal characteristics. In: Steel CM, editor. Disease markers. Special Issue. Human MHC Class II Antigens: Genetics, Structure and Function. Sussex: John Wiley & Sons, Ltd., 1984;2:363-69.

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