## Monoclonal Mouse Anti-Human CD10, Clone SS2/36

**Code No./ Code/ Code-Nr. F 0826 FITC-Conjugated**  
**Code No./ Code/ Code-Nr. R 0848 RPE-Conjugated**

### ENGLISH

**Intended use** For in vitro diagnostic use.  
 F 0826 and R 0848 are intended for use in flow cytometry. CD10 remains useful for the subclassification of B-lineage leukaemias, but it can also be found on other types of leukaemic cells (1). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

**Synonym for antigen** Common acute lymphoblastic leukaemia antigen (CALLA) (1).

**Introduction** CD10 is a 100 kDa protein belonging to a family of type II transmembrane metalloproteases which also includes the leukocyte antigens CD13, CD26, CD73 and aminopeptidase. The gene encoding CD10 is located on chromosome 3. CD10 is a zinc-dependent enzyme and is thought to down-regulate cellular responses to peptide hormones (2). On lymphoid cells, CD10 is expressed on immature T and B precursor cells but is lost as the cells reach maturation. CD10 is, however, re-expressed on proliferating B cells and mature neutrophils (3). In lymphoid malignancies, CD10 is expressed in acute lymphoblastic leukaemia (ALL) arising from precursor B cells, but is also observed in a proportion of T cell ALL. Additionally, it is expressed selectively in mature B-cell leukaemia, including multiple myeloma and in lymphomas (2, 3).

**Reagent provided** The Anti-CD10 conjugates, F 0826 and R 0848, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10<sup>6</sup> leucocytes from normal human peripheral blood).  
**Isotype:** IgG1, kappa. **Conjugate concentration mg/L:** See label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
F 0826	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927
R 0848	RPE (R-Phycoerythrin)	X 0928

**Immunogen** Acute lymphoblastic leukaemia cells.

**Specificity** Anti-CD10, SS2/36, was included in the Third International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD10 (4).

**Precautions**

- For professional users.
- This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
- As with any product derived from biological sources, proper handling procedures should be used.

**Storage** Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Staining procedure**

- Collect venous blood into a test tube containing an anticoagulant.
- Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
- Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
- Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD10.
- Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). The negative control may be used at a volume of 10 µL when adjusted to the same conjugate concentration as the test conjugate. However, the optimal amount of the negative control may vary, e.g. depending on the specimen and preparation method, and should be established by each individual laboratory.
- Incubate in the dark at 4 °C for 30 minutes.
- Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
- Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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