Monoclonal Mouse
Anti-Human CD71, Transferrin Receptor/FITC
 Clone Ber-T9
Code No./ Code/ Code-Nr. F 0829
 Edition/ Ausgabe 01.01.03

ENGLISH

Intended use For in vitro diagnostic use.
 F 0829 is intended for use in flow cytometry. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Synonyms for antigen T9, Transferrin Receptor (1)

Introduction CD71 is a 190 kDa disulfide-bonded homodimeric type II transmembrane glycoprotein. CD71 is the receptor for transferrin. CD71 is involved in iron (Fe³⁺) uptake, and expression is regulated by the metabolic demand for iron, which is essential for the growing cell (1). *In vitro* studies of human leukaemia cell lines, breast carcinoma cell line and fresh human leukaemic and glioblastoma cells have shown that CD71 is mainly expressed at the transition from G0/G1 to S-phase (2).
 CD71 is more abundantly expressed in rapidly dividing cells than quiescent cells (3, 4, 7). In normal tissues, high expression of CD71 is seen in the placental trophoblast and haemoglobin-synthesizing reticulocytes (3) but expression is lost in mature erythrocytes (4). Activated, but not resting lymphocytes express CD71 (5, 6). The proliferation status in both normal and neoplastic cells correlate with CD71 expression (3, 4).

Reagent provided The Anti-CD71 conjugate, F 0829, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10⁶ A431 cells).
Isotype: IgG1, kappa. **Conjugate concentration mg/L:** See label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
F 0829	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927

Specificity Anti-CD71, Ber-T9, was included in the Fifth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD71 (8).

Precautions

- For professional users.
- This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
- As with any product derived from biological sources, proper handling procedures should be used.

Storage Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Staining procedure

- Collect venous blood into a test tube containing an anticoagulant.
- Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
- Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
- Mix 100 µL cell suspension with 10 µL F 0829.
- Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
- Incubate in the dark at 4 °C for 30 minutes.
- Wash twice with PBS containing 2% BSA. Resuspend the cells in each test tube in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
- Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
1. Goding J. CD71. CD Guide. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 823-4.
2. Danova M, Riccardi A, Giordano M, Girino M, Mazzini G, Dezza L, et al. Cell cycle-related proteins: a flow cytometric study in human tumors. *Biol Cell* 1988;64:23-8.
3. Newman R, Schneider C, Sutherland R, Vodineli L, Greaves M. The transferrin receptor. [Review]. *Trends Biochem Sci* 1982;1:397-400.
4. Sutherland R, Delia D, Schneider C, Newman R, Kemshead J, Greaves M. Ubiquitous cell-surface glycoprotein on tumor cells is proliferation-associated receptor for transferrin. *Proc Natl Acad Sci (USA)* 1981;78:4515-9.

5. Cotner T, Williams JM, Christenson L, Shapiro HM, Strom TB, Strominger J. Simultaneous flow cytometric analysis of human T cell activation antigen expression and DNA content. *J Exp Med* 1983;157:461-72.
6. Klein J, ed. Natural history of the major histocompatibility complex. London: John Wiley & Sons, 1986:297-300.
7. Gatter KC, Brown G, Trowbridge IS, Woolston RE, and Mason DY. *J Clin Pathol* 1983;36:539-45
8. Boeker M, Werfel T. AA12.1. Binding of activation antigen panel mAb to mononuclear cells after antigen-specific stimulation in vitro. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 1181-2.

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