Monoclonal Mouse
Anti-Human CD15/FITC
 Clone C3D-1
Code No./ Code/ Code-Nr. F 0830

ENGLISH	
Intended use	For in vitro diagnostic use. F 0830 is intended for use in flow cytometry. In flow cytometry, anti-CD15 is considered essential for the initial evaluation of acute myeloid leukaemias together with a panel of other antibodies (1-3). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.
Synonyms for antigen	Lewis X, X-hapten, lacto- <i>N</i> -neo-fucopentaose III (LNFP III), stage-specific embryonic antigen (SSEA-1), 3-FAL and 3-FL (4,5).
Introduction	CD15 is a carbohydrate antigen carried by both glycoproteins and glycolipids. The structure of the antigen was determined to be Gal β 1-4[Fuc α 1-3]GlcNAc β 1-3Gal β 1-R. CD15 serves as a ligand for selectins and might also be involved in cell adhesion through a direct Lewis X-Lewis X interaction (5). CD15 is expressed mainly on mature granulocytes and monocytes but also on immature bone marrow cells and on leukaemic cells of the myelo-monocytic lineage (5). Additionally, CD15 is expressed in rare cases of acute lymphoblastic leukaemias (ALL) (5,6), and on adenocarcinoma cells and immature epithelial cells (5).
Reagent provided	The Anti-CD15 conjugate, F 0830, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 μ L of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Isotype:</u> IgM, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.
Immunogen	Purified human neutrophils (7).
Specificity	Anti-CD15, C3D-1, was included in the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD15 (5). Anti-CD15, C3D-1, reacts with peripheral blood granulocytes and monocytes (5). In acute myeloid leukaemia (AML), Anti-CD15, C3D-1, has been shown to label 21/35 (61.7%) of cases expressing CD15 in association with CD117 and/or CD34 (3). In acute lymphoblastic leukaemia (ALL) Anti-CD15, C3D-1, labelled 8/128 (6.25%) of B-lineage ALL and 0/26 T-lineage ALL cases (6).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Staining procedure	1. Collect venous blood into a test tube containing an anticoagulant. 2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. 3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. 4. In one test tube, mix 100 μ L cell suspension with 10 μ L F 0830. 5. In another test tube, mix 100 μ L cell suspension with an appropriate volume, e.g. 10 μ L, negative control reagent. The recommended negative control for F 0830 is DakoCytomation code No. X 0934. 6. Incubate in the dark at 4 °C for 30 minutes. 7. Wash twice with PBS containing 2% BSA. Resuspend the cells in each test tube in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. 8. Analyse on a flow cytometer. It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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