## Monoclonal Mouse Anti-Human CD13, Clone WM-47

Code No./ Code/ Code-Nr. **F 0831**      **FITC-Conjugated**  
 Code No./ Code/ Code-Nr. **R 0715**      **RPE-Conjugated**  
 Code No./ Code/ Code-Nr. **C 7073**      **RPE-Cy5-Conjugated**

### ENGLISH

**Intended use** For in vitro diagnostic use.  
 F 0831, R 0715 and C 7073 are intended for use in flow cytometry. Anti-CD13 is considered essential for the initial evaluation of acute myeloid leukaemias together with a panel of other antibodies (1, 2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

**Synonym for antigen** Aminopeptidase N (APN), EC 3.4.11.2, gp150 (3).

**Introduction** CD13 is a peptidase M1, categorized as a zinc metallopeptidase (3-5). It is a 150 kDa Type II integral membrane protein composed of 967 amino acids (3-6). CD13 is identical in structure to aminopeptidase N, a membrane bound zinc-binding metalloprotease, which degrades regulatory peptides produced by diverse cell types. The CD13 gene is located on chromosome 15 (3, 4). CD13 has many different functions. The molecule catalyzes the removal of single amino acids from the amino terminus of small peptides and probably plays a role in their digestion (3-6). CD13 participates in trimming peptides to MHC class II (3). It cleaves the MIP-1 chemokine that alters the target cell specificity of basophils to eosinophils (3). CD13 acts as a receptor for specific strains of RNA viruses (coronaviruses) that cause a relatively large percentage of upper respiratory tract infections (3-6). CD13 serves an important function for the early events in the interaction between human cytomegalovirus and the target cells. CMV incorporates the cellular CD13 protein in its envelope that may induce CD13-specific antibody response in certain immunocompromised patients (3).

CD13 is expressed on the surface of committed granulocyte-monocyte progenitors (CFU-GM) and by cells of the granulocyte and monocyte lineages at all stages of differentiation, as well as by malignant counterparts of these cells, including leukaemic blasts from a high percentage of acute and chronic myeloid leukaemias and a smaller percentage of lymphoid leukaemias (1-7). The antigen is also expressed on a small proportion of large granular lymphocytes but not on other lymphocytes as well as on a spectrum of non-haematopoietic tissues including fibroblasts, endothelial cells, epithelial cells from renal proximal tubules and intestinal brush border, bone marrow stromal cells, osteoclasts, and cells lining the bile duct canaliculi (3, 5-7).

**Reagent provided** The Anti-CD13 conjugates, F 0831, R 0715 and C 7073, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10<sup>6</sup> leucocytes from normal human peripheral blood).

Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

| Antibody Code No. | Fluorochrome                               | Negative Control Code No. |
|-------------------|--|---------------------------|
| F 0831            | FITC (Fluorescein Isothiocyanate Isomer 1) | X 0927                    |
| R 0715            | RPE (R-Phycoerythrin)                      | X 0928                    |
| C 7073            | RPE-Cy5 (R-Phycoerythrin-Cyanine 5)        | X 0955                    |

**Immunogen** Leukaemic cells from a case of chronic myeloid leukaemia (CML).

**Specificity** Anti-CD13, WM-47, was included in the Fourth and Fifth International Workshop and Conference on Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with the CD13 antigen (5, 6).

Anti-CD13, WM-47, reacts with CFU-GM cells, normal granulocytic and monocytic cells at all stages of differentiation (5-7). It does not react with lymphocytes and platelets (7).

Flow cytometric analysis demonstrated that Anti-CD13, WM-47, labels acute myeloid leukaemia (9/9 cases) and chronic myeloid leukaemia in myeloid blast crisis (7/7 cases). It does not react with common acute lymphoblastic leukaemia (0/8 cases) or T-cell acute lymphoblastic leukaemia (0/1 case) (7).

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

**Storage** Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Staining procedure**

1. Collect venous blood into a test tube containing an anticoagulant.
2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
4. Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD13.
5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
6. Incubate in the dark at 4 °C for 30 minutes.
7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

#### Product-specific limitations

Additionally, it has been observed that RPE-Cy5-conjugates may bind to monocytes, resulting in background staining (8).

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
1. van Dongen JJM, Adriaansen HJ. Immunobiology of leukemia. In: Henderson ES, Lister TA, Greaves MF, editors. Leukemia. Philadelphia, London, Toronto, Montreal, Sydney, Tokyo: WB Saunders Company; 1996. p. 83-130.
2. Thalhammer-Scherrer R, Mitterbauer G, Simonitsch I, Jaeger U, Lechner K, Schneider B, et al. The immunophenotype of 325 adult acute leukemias. Relationship to morphologic and molecular classification and proposal for a minimal screening program highly predictive for lineage discrimination. Am J Clin Pathol 2002;117:380-9.
3. Turni L, Shaw S, Watson B, Mason D. CD guide. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 761.
4. Riemann D, Kehlen A, Langner J. CD13 – not just a marker in leukemia typing. Immunology Today 1999;20:83-8.
5. Ashmun RA, Holmes KV, Shapiro LH, Favaloro EJ, Razak K, de Crom RPG, et al. M3. CD13 (aminopeptidase N) cluster workshop report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the Fifth International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 1. p. 771-5.
6. Look AT, Ashmun RA, Shapiro LH, O'Connell PJ, Gerkis V, D'Apice AJ, et al. M2.1. Report on the CD13 (aminopeptidase N) cluster workshop. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 784-7.
7. Favaloro EJ, Bradstock KF, Kabral A, Grimsley P, Zowtyj H, Zola H. Further characterization of human myeloid antigens (gp160,95; gp150; gp67): investigation of epitopic heterogeneity and non-haemopoietic distribution using panels of monoclonal antibodies belonging to CD-11b, CD-13 and CD-33. Br J Haematol 1988;69:163-71.
8. van Vugt MJ, van den Herik-Oudijk IE, van de Winkel JGJ. Binding of PE-CY5 conjugates to the human high-affinity receptor for IgG (CD64). Blood 1996;88:2358-61.

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Agilent Technologies Singapore (International) Pte Ltd.