Monoclonal Mouse Anti-Human CD14, Clone TÜK4

Code No./ Code/ Code-Nr. **F 0844** FITC-Conjugated
Code No./ Code/ Code-Nr. **R 0864** RPE-Conjugated

ENGLISH										
Intended use	For in vitro diagnostic use. F 0844 and R 0864 are intended for use in flow cytometry. Anti-CD14 is important for excluding CD14+ monocytes in relation to T lymphocyte subset analysis (1, 2). In flow cytometry, Anti-CD14 is considered essential for the initial evaluation of acute myeloid leukaemias together with a panel of other antibodies (3). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.									
Synonym for antigen	LPS receptor, LPS-R and the MO2, My4 and LeuM3 antigen (4, 7).									
Introduction	The mature CD14 molecule, which is a 53-55 kDa membrane protein, is composed of 356 amino acids and is anchored to the cell membrane by linkage to glycosyl phosphatidylinositol (GPI). In addition to the cell surface form, two soluble forms have been described. One soluble form is produced by shedding the cell surface form, which results in an approximately 48 kDa molecule (4-7). A second soluble form is released from cells before addition of the GPI anchor, yielding a higher molecular weight molecule (6). CD14 functions as a receptor for endotoxin (lipopolysaccharide (LPS)) and is the major receptor involved in the lethal response to both endotoxin and Gram-negative bacteria (4, 6). Antibodies against CD14 inhibit LPS-induced tumor necrosis factor alpha (TNF α) production (4, 6). CD14 is expressed strongly on the surface of monocytes and weakly on the surface of granulocytes. Most tissue macrophages and Langerhans' cells also express the antigen (5-7). The CD14 antigen is expressed on leukemic cells of the myelo-monocytic lineage especially AML FAB Type M4, M5a and M5b (3).									
Reagent provided	The Anti-CD14 conjugates, F 0844 and R 0864, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L Na ₂ S ₂ O ₃ , pH 7.2. Each vial contains 100 tests (10 μ L of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG2a, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.									
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 0844</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0933</td> </tr> <tr> <td>R 0864</td> <td>RPE (R-Phycoerythrin)</td> <td>X 0950</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 0844	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0933	R 0864	RPE (R-Phycoerythrin)	X 0950
Antibody Code No.	Fluorochrome	Negative Control Code No.								
F 0844	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0933								
R 0864	RPE (R-Phycoerythrin)	X 0950								
Specificity	Anti-CD14, TÜK4, was included in the Fourth, Fifth and Sixth International Workshop and Conference on Leukocyte Cell Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with the CD14 antigen (5-7). Anti-CD14, TÜK4 reacts strongly with peripheral blood monocytes and shows high reactivity to granulocytes. It also reacts with perivascular macrophages. The antibody labels 50% of Langerhans' cells in epidermal cell suspensions and is therefore characterized as weakly reactive (7).									
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.									
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.									
Staining procedure	1. Collect venous blood into a test tube containing an anticoagulant. 2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. 3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. 4. Mix 100 μ L cell suspension with 10 μ L fluorochrome-conjugated Anti-CD14. 5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). 6. Incubate in the dark at 4 °C for 30 minutes. 7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. 8. Analyse on a flow cytometer.									

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
- Centers for disease control. Guidelines for the performance of CD4+ T-cell determinations in persons with human immunodeficiency virus infection. MMWR 1992;41 (No. RR-8);001. CDC Publication date 05/08/1992.
- National Committee for Clinical Laboratory Standards (NCCLS) Document. Clinical applications of flow cytometry: Quality assurance and immunophenotyping of peripheral blood lymphocytes. H42-T 1992;12:1-76.
- van Dongen JJM, Adriaansen HJ. Immunobiology of leukemia. In: Henderson ES, Lister TA, Greaves MF, editors. Leukemia. Philadelphia, London, Toronto, Montreal, Sydney, Tokyo: WB Saunders Company; 1996. p. 83-130.
- Turni L, Shaw S, Watson B, Mason D. CD guide. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, De Haas M, et al, editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 June; Harrogate, UK. Oxford, New York: Oxford University Press Inc.; 2002. p. 761-3.
- Goyert SM, Haziot A, Jiao D, Katz IR, Kruger C, Ross J, Schütt C. M4. CD14 cluster workshop report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, Ritz J, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the Fifth International Workshop and Conference; Held in Boston, USA 3-7 November 1993. Volume One. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 778-82.
- Goyert SM, Cohen L, Gangloff SC, Ashmun R, Haefner-Cavaillon N. MC4. CD14 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997.p. 963-5.
- Goyert SM, Tesio L, Ashman LK, Ball E, Bazil V, Garrido F, Hogg N, Horejsi V, et al. M3.1. Report on the CD14 cluster workshop. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, von dem Borne AEG, editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference on Human Leucocyte Differentiation Antigens; 1989 Februarv 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 789-94.

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