Monoclonal Mouse
Anti-Human CD30/FITC
 Clone Ber-H2
Code No./ Code/ Code-Nr. F 0849
 Edition/ Ausgabe 01.01.03

ENGLISH							
Intended use	For in vitro diagnostic use. F 0849 is intended for use in flow cytometry. The antibody labels anaplastic large-cell lymphoma (ALCL) and Reed-Sternberg cells, and is a useful tool for the identification of ALCL and as a secondary marker for Hodgkin's disease (1). Differential identification is aided by the results from a panel of antibodies. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.						
Synonym for antigen	Ki-1 Antigen (1)						
Introduction	CD30 is a transmembrane cytokine receptor belonging to the tumour necrosis factor (TNF) receptor superfamily. Mature CD30 has a molecular mass of 120 kDa and is derived from a 90 kDa precursor protein (2). The extracellular domain of CD30, is homologous to that of TNF receptor superfamily members, whereas there is no homology in the cytoplasmic domain, suggesting major differences in signalling mechanisms (3). The intracellular part of CD30 possesses kinase activity, indicating that CD30 plays a role in regulating the function, differentiation and/or proliferation of normal lymphoid cells (2). A soluble 85 kDa form of CD30, sCD30, released from the membrane-bound molecule by proteolytic cleavage, can be detected in the sera of patients with CD30-expressing neoplasms (3). CD30 expression is found on Hodgkin and Reed-Sternberg (H-RS) cells, anaplastic large-cell lymphoma cells, and on activated B and T lymphocytes (2, 5, 6).						
Reagent provided	The Anti-CD30 conjugate, F 0849, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ L540 cells). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.						
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 0849</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 0849	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927
Antibody Code No.	Fluorochrome	Negative Control Code No.					
F 0849	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927					
Immunogen	Co cells, a Hodgkin cell line (1).						
Specificity	The antibody from the Ber-H2 hybridoma cell line was included in the Fourth International Workshop on Human Leucocyte Differentiation Antigens (Vienna 1989) and was found to belong to the CD30 cluster (4). The epitope recognized by anti CD30, Ber-H2, is located extracellularly between amino acid residues 112 and 412 (3). The antibody labels: Cell lines derived from Hodgkin's disease, L428, L540, L591, Co, Ho and KM-H2; HTLV-1 transfected T-cell lines, Hut-102 and MT-2; EBV-transformed B-cell lines (non-Burkitt), B95-8 (monkey), BJA-B and Cess; and the myeloid cell line, K 562 (1). Outside the lympho-haematopoietic system, the anti CD30, Ber-H2, labels embryonal carcinomas and embryonal elements of mixed germ cell tumours (7).						
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.						
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.						
Staining procedure	1. Collect venous blood into a test tube containing an anticoagulant. 2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. 3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. 4. Mix 100 µL cell suspension with 10 µL F 0849. 5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). 6. Incubate in the dark at 4 °C for 30 minutes. 7. Wash twice with PBS containing 2% BSA. Resuspend the cells in each test tube in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.						

8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
- Schwarting R, Gerdes J, Dürkop H, Falini B, Pileri S, Stein H. Ber-H2: a new anti-Ki-1 (CD30) monoclonal antibody directed at a formol-resistant epitope. Blood 1989;74:1678-89.
- de Bruin PC, Gruss H-J, van der Valk P, Willemze R, Meijer CJLM. CD30 expression in normal and neoplastic lymphoid tissue: biological aspects and clinical implications [review]. Leukemia 1995;9:1620-7.
- Falini B, Pileri S, Pizzolo G, Dürkop H, Flenghi L, Stirpe F, et al. CD30 (Ki-1) molecule: a new cytokine receptor of the tumor necrosis factor receptor superfamily as a tool for diagnosis and immunotherapy [Review]. Blood 1995;85:1-14.
- Schwarting R, Stein H. A3. Cluster report: CD30. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p 419-22.
- Stein H, Gerdes J, Schwab U, Lemke H, Mason DY, Ziegler A, et al. Identification of Hodgkin and Sternberg-Reed cells as a unique cell type derived from a newly-detected small-cell population. Int J Cancer 1982;30:445--59.
- Stein H, Mason DY, Gerdes J, O'Connor N, Wainscoat J, Pallesen G, et al. The expression of the Hodgkin's disease associated antigen Ki-1 in reactive and neoplastic lymphoid tissue: evidence that Reed-Sternberg cells and histiocytic malignancies are derived from activated lymphoid cells. Blood 1985;66:848-58.
- Pallesen G, Hamilton-Dutoit SJ. Ki-1 (CD30) antigen is regularly expressed by tumor cells of embryonal carcinoma. Am J Pathol 1988;133:446-50.

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