Monoclonal Mouse
Anti-Human Epithelial Antigen/FITC
 Clone Ber-EP4
 Code No./ Code/ Code-Nr. **F 0860**

ENGLISH

Intended use For in vitro diagnostic use.
 F 0860 is intended for use in flow cytometry. Anti-Human Epithelial Antigen can be used for detecting circulating tumor cells in peripheral blood, serous effusions and peritoneal washings together with a panel of other antibodies (1-3). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Synonym for antigen Ber-EP4 (1-7).

Introduction Cytologic diagnosis of reactive or malignant effusions using immunohistochemistry or flow cytometry is often difficult, especially in cases involving mesothelioma versus adenocarcinoma with overlapping histological features (3-7). The sensitivity of conventional cytology for the detection of malignant cells in effusions has been improved by using antibody combinations that include markers that are positive for mesothelioma and markers negative for mesothelioma and positive for carcinoma (4-7). Vimentin, HBME-1 and thrombomodulin are positive markers for mesothelioma (4, 5). The most commonly used epithelial markers of carcinomas include CEA, B72.3 (TAG-72), Leu-M1 (CD15), MOC-31, AH-6, HB-TN1, CA-125, and Ber-EP4 (3-7). In order to achieve optimal sensitivity and specificity of immunohistochemical or flow cytometric analysis it has been recommended to use antibody panels composed of 2, 3 or 4 carcinoma markers that includes Ber-EP4 (3-6).

Reagent provided The Anti-Human Epithelial Antigen conjugate, F 0860, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).
Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
F 0860	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927

Immunogen MCF-7 cell line (8).

Specificity Anti-Human Epithelial Antigen, Ber-EP4, is directed against a partially formalin-resistant epitope of a human epithelial antigen consisting of two glycoproteins of 34 and 39 kDa. The antigen is present on the surface membrane and in the cytoplasm of all epithelial cells except the superficial layers of squamous epithelia, hepatocytes and parietal cells (8).
 Flow cytometric analysis of malignant effusions has shown that Anti-Human Epithelial Antigen, Ber-EP4, labels carcinoma samples (16/17 cases (2) and 36/44 cases (3)). A few mesothelioma samples (2/3 cases) were also stained (3).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

Storage Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Staining procedure

1. Collect venous blood into a test tube containing an anticoagulant.
2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
4. Mix 100 µL cell suspension with 10 µL F 0860.
5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
6. Incubate in the dark at 4 °C for 30 minutes.
7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
1. Kularatne BY, Lorigan P, Browne S, Suvarna SK, Smith MO, Lawry J. Monitoring tumour cells in the peripheral blood of small cell lung cancer patients. *Cytometry* 2002;50:160-7.
2. Risberg B, Davidson B, Dong HP, Nesland JM, Berner A. Flow cytometric immunophenotyping of serous effusions and peritoneal washings: comparison with immunocytochemistry and morphological findings. *J Clin Pathol* 2000;53:513-7.
3. Davidson B, Dong HP, Berner A, Christensen J, Nielsen S, Johansen P, et al. Detection of malignant epithelial cells in effusions using flow cytometric immunophenotyping. *Am J Clin Pathol* 2002;118:85-92.
4. Lozano MD, Panizo A, Toledo GR, Sola JJ, Pardo-Mindán J. Immunocytochemistry in the differential diagnosis of serous effusions. A comparative evaluation of eight monoclonal antibodies in Papanicolaou stained smears. *Cancer* 2001;93:68-72.
5. Garcia-Prats MD, Ballestin C, Sotelo T, Lopez-Encuentra A, Mayordomo JI. A comparative evaluation of immunohistochemical markers for the differential diagnosis of malignant pleural tumours. *Histopathology* 1998;32:462-72.
6. Motherby H, Kube M, Friedrichs N, Nadjari B, Knops K, Donner A, et al. Immunocytochemistry and DNA-image cytometry in diagnostic effusion cytology. I. Prevalence of markers in tumour cell positive and negative smears. *Anal Cell Pathol* 1999;19:7-20.
7. Roberts F, Harper CM, Downie I, Burnett RA. Immunohistochemical analysis still has a limited role in the diagnosis of malignant mesothelioma. *Am J Clin Pathol* 2001;116:253-62.
8. Latza U, Niedobitek G, Schwarting R, Nekarda H, Stein H. Ber-EP4: new monoclonal antibody which distinguishes epithelia from mesothelia. *J Clin Pathol* 1990;43:213-9.

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Agilent Technologies Singapore (International) Pte Ltd.