Monoclonal Mouse Anti-Human CD235a, Glycophorin A, Clone JC159

Code No./ Code/ Code-Nr. F 0870 **FITC-Conjugated**
Code No./ Code/ Code-Nr. R 7078 **RPE-Conjugated**

ENGLISH										
Intended use	For in vitro diagnostic use. F 0870 and R 7078 are intended for use in flow cytometry. Anti-CD235a, JC159, is a useful tool for the flow cytometric discrimination of the later erythroid cells from lymphoid cells (1), and, generally, antibodies to CD235a have been used for the identification of erythroleukaemia (2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.									
Synonyms for antigen	Sialoglycoprotein, GPA (2, 3).									
Introduction	Glycophorin A, CD235a, is a single-pass membrane sialoglycoprotein and carrier of the blood group specificities M and N. The gene is located on chromosome 4q28-q31 and share sequence homology with genes encoding glycophorin B and glycophorin E (3). Glycophorin A is expressed by erythroid cells beginning on morphologically recognizable erythroid precursors, just after the CFU-E stage, to the mature erythrocyte. Once glycophorin A is maximally expressed, the amount in each erythroid cell remains constant and does not appear to change during further maturation (1). The majority of cases of erythroleukaemia express glycophorin A on neoplastic erythroblasts, whereas acute myeloid leukaemia and acute lymphoblastic leukaemia only very rarely express glycophorin A (2, 4).									
Reagent provided	The Anti-CD235a conjugates, F 0870 and R 7078, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ erythrocytes from normal human peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.									
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 0870</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> <tr> <td>R 7078</td> <td>RPE (R-Phycoerythrin)</td> <td>X 0928</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 0870	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927	R 7078	RPE (R-Phycoerythrin)	X 0928
Antibody Code No.	Fluorochrome	Negative Control Code No.								
F 0870	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927								
R 7078	RPE (R-Phycoerythrin)	X 0928								
Immunogen	Membrane extracts of hairy cell leukemia prepared from human spleen (2).									
Specificity	Anti-CD235a, JC159 strongly labels normal erythroid cells at all stages of differentiation from the erythroblast to mature red cell (2).									
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. 									
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.									
Staining procedure	<ol style="list-style-type: none"> Collect venous blood into a test tube containing an EDTA anticoagulant. Isolate erythrocytes by centrifugation. Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD235a. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). Incubate in the dark at 4 °C for 30 minutes. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. Analyse on a flow cytometer. <p>It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis</p>									

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