Monoclonal Mouse Anti-Human CD16, Fc Gamma Receptor III/FITC, Clone DJ130c Code F7011
Monoclonal Mouse Anti-Human CD16, Fc Gamma Receptor III/RPE, Clone DJ130c Code R7012

ENGLISH

Intended use	For in vitro diagnostic use. F7011 and R7012 are intended for use in flow cytometry. The immunological status of patients can be determined by assessing the level of CD16+ natural killer (NK) cells, and a simultaneous detection of other NK-cell markers and enumeration of other lymphocyte subpopulations. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.									
Synonym for antigen	Low affinity Fc receptor for IgG (FcγRIII).									
Summary and explanation	CD16 is a 50-70 kDa glycoprotein which occurs in two isoforms, CD16a and CD16b. CD16a is a transmembrane molecule expressed on about 90% of NK cells and also found on macrophages and subsets of monocytes and T cells. CD16b is glycosyl phosphatidyl inositol-anchored and is expressed on virtually all neutrophils. CD16b is polymorph with two alleles termed NA1 and NA2 (1, 2). CD16 represents the functional receptor for antibody-dependent cellular cytotoxicity (ADCC). Thus, CD16 ^{bright} , CD56 ^{dim} NK cells have high ADCC, whereas CD16 ^{dim} , CD56 ^{bright} cells have high lymphokine-activated killer (LAK) activity (3). Phenotyping and quantitation of NK cells, together with evaluation of other lymphocyte subpopulations can be of great importance in various disorders, and may provide prognostic information (4-6).									
Reagent provided	The Anti-CD16 conjugates, F7011 and R7012, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.									
	<table border="1"> <thead> <tr> <th>Antibody Code</th> <th>Fluorochrome</th> <th>Negative Control Code</th> </tr> </thead> <tbody> <tr> <td>F7011</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X0927</td> </tr> <tr> <td>R7012</td> <td>RPE (R-Phycoerythrin)</td> <td>X0928</td> </tr> </tbody> </table>	Antibody Code	Fluorochrome	Negative Control Code	F7011	FITC (Fluorescein Isothiocyanate Isomer 1)	X0927	R7012	RPE (R-Phycoerythrin)	X0928
Antibody Code	Fluorochrome	Negative Control Code								
F7011	FITC (Fluorescein Isothiocyanate Isomer 1)	X0927								
R7012	RPE (R-Phycoerythrin)	X0928								
Specificity	Anti-CD16, DJ130c, was included in the Fifth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD16 (7). The antibody reacts with CD16a and both variants of CD16b. The epitope has been located to the first membrane-distal domain of CD16 (2).									
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. Unused solution should be disposed of according to local, State and Federal regulations. 									
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.									
Staining procedure	<ol style="list-style-type: none"> Collect venous blood into a test tube containing an anticoagulant. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD16. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). Incubate in the dark at 4 °C for 30 minutes. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. Analyse on a flow cytometer. <p>It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.</p>									

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
- Tetteroo PAT, Van Der Schoot CE, Visser FJ, Bos MJE, Von Dem Borne AEGKR. M5.3. Three different types of Fcγ receptors on human leucocytes defined by Workshop antibodies; FcγR_{LOW} of neutrophils, FcγR_{LOW} of K/NK lymphocytes, and FcγRII. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 702-6.
- Tamm A, Schmidt RE. The binding epitopes of human CD16 (FcγRIII) monoclonal antibodies. Implications for ligand binding. J Immunol 1996;157:1576-81.
- Cooper MA, Fehniger TA, Caligiuri MA. The biology of human natural killer-cell subsets. Trends Immunol 2001;22:633-40.
- Hu PF, Hultin LE, Hultin P, Hausner MA, Hirji K, Jewett A, et al. Natural killer cell immunodeficiency in HIV disease is manifest by profoundly decreased numbers of CD16⁺CD56⁺ cells and expansion of a population of CD16^{dim} CD56⁺ cells with low lytic activity. J Acquir Immune Defic Syndr Hum Retrovirol 1995;10:331-40.
- Bruunsgaard H, Pedersen C, Skinhøj P, Pedersen BK. Clinical progression of HIV infection: role of NK cells. Scand J Immunol 1997;46:91-5.
- Emmer PM, Nelen WL, Steegers EA, Hendriks JC, Veerhoek M, Joosten I. Peripheral natural killer cytotoxicity and CD56^{pos}CD16^{pos} cells increase during early pregnancy in women with a history of recurrent spontaneous abortion. Hum Reprod 2000;15:1163-9
- Schmidt RE. M6. CD16 cluster workshop report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 2. p. 805-6.

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