## Monoclonal Mouse Anti-Human CD23, Clone MHM6

**Code No./ Code/ Code-Nr. F 7062**      **FITC-Conjugated**  
**Code No./ Code/ Code-Nr. R 7108**      **RPE-Conjugated**

### ENGLISH

**Intended use** For in vitro diagnostic use.  
 F 7062 and R 7108 are intended for use in flow cytometry. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

**Synonym for antigen** FcεRII, BLAST-2, low affinity IgE receptor (3).

**Introduction** CD23 is a 45 kDa type II membrane glycoprotein belonging to the C-type lectin family, and identical to the low affinity IgE receptor (FcεRII) found on B cells (2, 4). CD23 is primarily expressed on B cells and monocytes including a strong expression on EBV-transformed B lymphoblasts (1, 3). Interleukin-4 and Interleukin-13 stimulation increase CD23 expression on B cells (2). On B cells, CD23 is selectively expressed on sIgM/sIgD double-bearing cells, and is lost upon differentiation into Ig-secreting cells (3).  
 CD23 is also present on a large variety of other cells such as T cells, eosinophils, platelets, Langerhans cells and a subset of thymic epithelial cells (3). However, CD23 is absent in a large number of normal non-lymphoid human tissues, such as, central nervous system, endocrine glands, gastrointestinal, respiratory, urinary and female/male reproductive system, skin, bone, cartilage, connective tissue, muscles, peripheral nerves, vessels, mesothelium, synovium, placenta and umbilical cord (5).

**Reagent provided** The Anti-CD23 conjugates, F 7062 and R 7108, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10<sup>6</sup> leucocytes from normal human peripheral blood).

Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
F 7062	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927
R 7108	RPE (R-Phycoerythrin)	X 0928

**Immunogen** EBV-transformed B-lymphoblastoid cells (5).

**Specificity** Anti-CD23, MHM6, was included in the Second, Third, Fourth and Fifth International Workshop and Conference on Human Leucocyte Differentiation Antigens (Boston 1984, Oxford 1986, Vienna 1989, Boston 1993), and studies by a number of laboratories confirmed its reactivity with the CD23 antigen (3). The CD23 antigen is expressed on neoplastic cells from cases of B-cell chronic lymphocytic leukaemia and some cases of centroblastic/centrocytic lymphoma, but not in other types of lymphoid neoplasms (5).

**Precautions** 1. For professional users.  
 2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.  
 3. As with any product derived from biological sources, proper handling procedures should be used.

**Storage** Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Staining procedure**

1. Collect venous blood into a test tube containing an anticoagulant.
2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
4. Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD23.
5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
6. Incubate in the dark at 4 °C for 30 minutes.
7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
8. Analyze on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
- Sarfati M. BC7. CD23 Workshop Panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p.144-7.
- Kikutani H, Suemura M, Owaki H, Yamasaki K, Barsumian EL, Nakamura H, et al. B3.3. CD23 is a low-affinity Fce receptor (FceR): regulation of FceR expression. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 419-22.
- Sarfati M, Ishihara H, Delespesse G. B7. CD23 Workshop Panel report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 530-3.
- Aubry JP, Bonnefoy JY, De Vries JE, Banchereau J. B3.2. The CD23 antigen is the human lymphocyte receptor for IgE. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 417-9.
- Pallesen G. B2.5. The distribution of CD23 in normal human tissues and in malignant lymphomas. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 383-6.

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