Monoclonal Mouse Anti-Human CD41, Platelet Glycoprotein IIb, Clone 5B12

Code No./ Code/ Code-Nr. F 7088 FITC-Conjugated
Code No./ Code/ Code-Nr. R 7058 RPE-Conjugated

ENGLISH

Intended use	For in vitro diagnostic use. F 7088 and R 7058 are intended for use in flow cytometry. CD41 is a selective marker of platelets and platelet precursors, and antibodies to CD41 may be of value for the immunophenotyping of megakaryoblastic leukaemias, and for the identification of Glanzmann's thrombasthenia (1). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.									
Synonyms for antigen	Glycoprotein IIb (GPIIb), α IIb integrin (1).									
Introduction	CD41 is the α subunit of the CD41/CD61 complex (GPIIb-IIIa), a calcium-dependent, non-covalently associated heterodimer. CD41 is a 135 kDa molecule, consisting of 2 peptides: a 120 kDa α chain and a 23 kDa β chain. The α chain and the β chain are linked by a single disulphide bond. The α chain contains four calcium-binding sites and is entirely extracellular, while the β -chain has extracellular, transmembrane and cytoplasmic domains (1). The CD41/CD61 complex appears early in megakaryocyte maturation (2). The activated CD41/CD61 complex is a receptor for von Willebrand factor, soluble fibrinogen and fibronectin (1) and plays a central role in platelet activation and aggregation (1, 3). CD41 is expressed to a variable degree in megakaryoblastic/cytic leukaemias (4). It is absent from, or defective on the platelets of patients suffering from Glanzmann's thrombasthenia (5).									
Reagent provided	The Anti-CD41 conjugates, F 7088 and R 7058, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN_3 , pH 7.2. Each vial contains 100 tests (10 μL of conjugate for up to 10^6 platelets from normal peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> see label on vial.									
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 7088</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> <tr> <td>R 7058</td> <td>RPE (R-Phycoerythrin)</td> <td>X 0928</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 7088	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927	R 7058	RPE (R-Phycoerythrin)	X 0928
Antibody Code No.	Fluorochrome	Negative Control Code No.								
F 7088	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927								
R 7058	RPE (R-Phycoerythrin)	X 0928								
Immunogen	Normal human platelets.									
Specificity	Anti-CD41, 5B12, was included in the Fifth International Workshop and Conference on Human Leucocyte Differentiation Antigens (Boston 1993) and studies by a number of laboratories confirmed its reactivity with CD41. Erroneously the clone was listed as SB12 (6).									
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN_3), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. 									
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.									
Staining procedure	<ol style="list-style-type: none"> Collect venous blood into a test tube containing EDTA as an anticoagulant. Within 5 minutes centrifuge at 200 x g for 5 minutes at room temperature. Collect 100 μL of the upper platelet-rich plasma (this is sufficient for 20 tests) and mix into 1 mL of 1% paraformaldehyd in 0.01 mol/L PBS, pH 7.4. Incubate at 4 °C for ½-1 hour. Centrifuge at 1200 x g for 5 minutes at room temperature. Gently aspirate the supernatant and discard it leaving approximately 50 μL of fluid. Add 2 mL 0.01 mol/L PBS, pH 7.4, and resuspend the platelets by using a vortex mixer. Centrifuge at 1200 x g for 5 minutes at room temperature. Gently aspirate the supernatant and discard it leaving approximately 50 μL of fluid. Add 1 mL 2% foetal calf serum in 0.01 mol/L PBS and resuspend the platelets by using a vortex mixer. Mix 50 μL platelet suspension with 10 μL fluorochrome-conjugated Anti-CD41. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). 									

10. Incubate in the dark at 4 °C for 20 minutes.
11. Wash twice with PBS containing 2% BSA. Resuspend the platelets in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
12. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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Agilent Technologies Singapore (International) Pte Ltd.