## Monoclonal Mouse Anti-Human CD38, Clone AT13/5

Code No./ Code/ Code-Nr. **F 7101 FITC-Conjugated**  
 Code No./ Code/ Code-Nr. **R 7144 RPE-Conjugated**

### ENGLISH

<b>Intended use</b>	For in vitro diagnostic use. F 7101 and R 7144 are intended for use in flow cytometry. CD38 is a useful marker in the immuno-phenotyping of acute leukaemias. Additionally, antibodies to CD38 are valuable for the identification of plasma cells, as poorly differentiated plasma cells may mimick other blastic lymphoid cells (1). Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.									
<b>Synonym for antigen</b>	T10 (1).									
<b>Introduction</b>	<p>CD38 is a 45 kDa type II integral membrane glycoprotein consisting of an N-terminal 20 amino acid cytoplasmic domain, a 23 amino acid transmembrane domain, and a 257 amino acid extracellular domain. CD38 is associated with different molecules in different cell types and these unique combinations cause unique cell-type specific signalling. One of the known ligands that activates CD38-mediated signalling is CD31 (2).</p> <p>CD38 is expressed by a number of different cell-types and the pattern of expression often changes with differentiation and maturation. Thus, CD38 expression is found during early stages of B and T-cell differentiation, is lacking during the later intermediate stages, and is regained near the end of maturation. CD38 is expressed at variable levels on the majority of haematopoietic cells, prevalently during early differentiation and activation, and it is expressed strongly on plasma cells (3-5). CD38-positive cell types outside the haematopoietic system include subsets from the gut, brain, prostate, pancreas, bone, eye and muscle (2).</p> <p>CD38 is expressed on approximately 60% of CD34-positive peripheral blood mononuclear cells. The least mature CD34-positive cells are characterized by a lack of CD38 (6).</p>									
<b>Reagent provided</b>	<p>The Anti-CD38 conjugates, F 7101 and R 7144, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10<sup>6</sup> leucocytes from normal peripheral blood).</p> <p><u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> see label on vial.</p>									
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 7101</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> <tr> <td>R 7144</td> <td>RPE (R-Phycoerythrin)</td> <td>X 0928</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 7101	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927	R 7144	RPE (R-Phycoerythrin)	X 0928
Antibody Code No.	Fluorochrome	Negative Control Code No.								
F 7101	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927								
R 7144	RPE (R-Phycoerythrin)	X 0928								
<b>Immunogen</b>	Namalwa cells (B lymphocyte cell line from a Burkitt's lymphoma) (7).									
<b>Specificity</b>	Anti-CD38, AT13/5, was included in the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens (Kobe 1996), and studies by a number of laboratories confirmed its reactivity with CD38 (5).									
<b>Precautions</b>	<ol style="list-style-type: none"> <li>For professional users.</li> <li>This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>As with any product derived from biological sources, proper handling procedures should be used.</li> </ol>									
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.									
<b>Staining procedure</b>	<ol style="list-style-type: none"> <li>Collect venous blood into a test tube containing an anticoagulant.</li> <li>Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.</li> <li>Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.</li> <li>Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD38.</li> <li>Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).</li> <li>Incubate in the dark at 4 °C for 30 minutes.</li> <li>Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.</li> <li>Analyse on a flow cytometer.</li> </ol> <p>It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.</p>									

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
1. Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 1999. p. 87-8.
2. Deaglio S, Mehta K, Malavasi F. Human CD38: a (r)evolutionary story of enzymes and receptors (review). Leuk Res 2001;25:1-12.
3. Jackson DG, Bell JI. Isolation of a cDNA encoding the human CD38 (T10) molecule, a cell surface glycoprotein with an unusual discontinuous pattern of expression during lymphocyte differentiation. J Immunol 1990;144:2811-5.
4. Terstappen LW, Hollander Z, Meiners H, Loken MR. Quantitative comparison of myeloid antigens on five lineages of mature peripheral blood cells. J Leukoc Biol 1990;48:138-48.
5. Morra M, Deaglio S, Mallone R, di Rosa G, Tibaldi E, Zini G, et al. BC10. CD38 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 151-4.
6. Herbein G, Sovalat H, Wunder E, Baerenzung M, Bachorz J, Lewandowski H, et al. Isolation and identification of two CD34+ cell subpopulations from normal human peripheral blood. Stem Cells 1994;12:187-97.
7. Ellis JH, Barber KA, Tutt A, Hale C, Lewis AP, Glennie MJ, et al. Engineered anti-CD38 monoclonal antibodies for immunotherapy of multiple myeloma. J Immunol 1995;155:925-37.

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Agilent Technologies Singapore (International) Pte Ltd.