Monoclonal Mouse
Anti-Human CD43/FITC
 Clone DF-T1
Code No./ Code/ Code-Nr. F 7102
 Edition/ Ausgabe 01.03.03

ENGLISH							
Intended use	For in vitro diagnostic use. F 7102 is intended for use in flow cytometry. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.						
Synonyms for antigen	Leucosialin, gpL 115, sialophorin, leucocyte sialoglycoprotein or LSGP (3, 5).						
Introduction	CD43 is a human heavily O-glycosylated transmembrane protein, which is expressed by T cells and myeloid lineage cells (1, 3). Reactivity with plasma cells has been observed (2). The molecule is known to be defective in the Wiskott-Aldrich syndrome and may participate in T-cell activation (3). Immunoblotting studies have revealed molecular mass heterogeneity displayed by the antigen extracted from white blood cells of different lineage. Thus, two major bands of 140 and 130 kDa are seen for the promyelocytic cell line HL60, 110 kDa for the myeloblast cell line KG1 and a major 115 kDa band for peripheral blood lymphocytes (3).						
Reagent provided	The Anti-CD43 conjugate, F 7102, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.						
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 7102</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 7102	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927
Antibody Code No.	Fluorochrome	Negative Control Code No.					
F 7102	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927					
Immunogen	Myeloblast cell line KG1 (3).						
Specificity	CD43 was first clustered at the Third International Workshop and Conference on Human Leucocyte Differentiation Antigens (Oxford 1986) (4). Molecular mass determinations, cross blocking studies, tissue staining specificities and CD43 COS transfectant studies confirmed that Anti-CD43, DF-T1 belongs to the group of antibodies reacting with CD43 at the Fourth International Workshop and Conference on Human Leucocyte Differentiation Antigens (Vienna 1989) where Anti-CD43, DF-T1, was included (5). Anti-CD43, DF-T1, recognizes an antigen expressed by normal and neoplastic T cells. A small number of B-cell leukaemias and lymphomas usually of the small cell type, e.g. B-cell chronic lymphocytic leukaemia and/or centrocytic lymphoma, may also express the antigen (3).						
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.						
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.						
Staining procedure	<ol style="list-style-type: none"> Collect venous blood into a test tube containing an anticoagulant. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. Mix 100 µL cell suspension with 10 µL F 7102. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). Incubate in the dark at 4 °C for 30 minutes. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. Analyse on a flow cytometer. <p>It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.</p>						

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
1. Borche L, Lozano F, Vilella R, Vives J. CD43 monoclonal antibodies recognize the large sialoglycoprotein of human leukocytes. Eur J Immunol 1987;17:1523-6.
2. Beschoner R, Horny H-P, Petruch UR, Kaiserling E. Frequent expression of haemopoietic and non-haemopoietic antigens by reactive plasma cells: an immunohistochemical study using formalin-fixed, paraffin-embedded tissue. Histo Histopathol 1999;14:805-12.
3. Stross WP, Warnke RA, Flavell DJ, Flavell SU, Simmons D, Gatter KC, et al. Molecule detected in formalin fixed tissue by antibodies MT1, DF-T1, and L60 (Leu-22) corresponds to CD43 antigen. J Clin Pathol 1989;42:953-61.
4. Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1 family, and leucocyte common antigens: new and previously defined clusters. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leukocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 788-803.
5. Stoll M, Dalchau R, Schmidt RE. N5: Cluster report: CD43. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leukocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 604-8.

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