**Monoclonal Mouse**  
**Anti-Human B Cell/FITC**  
 Clone FMC7  
**Code No./ Code/ Code-Nr. F 7110**  
 Edition/ Ausgabe 01.03.03

ENGLISH							
<b>Intended use</b>	For in vitro diagnostic use. F 7110 is intended for use in flow cytometry. Anti-B Cell, FMC7 is considered essential for the initial evaluation of B-cell chronic lymphoproliferative disorders together with a panel of other antibodies (1-2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.						
<b>Synonym for antigen</b>	FMC7 (1-8).						
<b>Introduction</b>	The 105 kDa FMC7 antigen is present on a subpopulation of functionally mature B cells (3-5). The antigen is an integral membrane protein, not affected by cell cycle status, which is present in similar quantities in normal (resting) B cells, in leukaemic cells, and in rapidly cycling cell lines. It is able to interact with the cytoskeleton (6). The expression of the antigen does not correlate with activation markers such as CD23 (5). Studies have shown that FMC7 antigen may be a conformational epitope of the CD20 B cell antigen (7). The FMC7 antigen is present on B cells in normal peripheral blood and isolated tonsil samples (4, 5). The antigen is expressed on chronic lymphoproliferative disorders including B-cell prolymphocytic leukaemia (B-PLL), hairy cell leukaemia (HCL), mantle cell lymphoma (MCL), follicle center lymphoma (FCL), marginal zone lymphoma (MZL), mucosa-associated lymphoid tissue lymphoma (MALT), diffuse large cell lymphoma (DLCL), Burkitt's lymphoma (BL) and Waldenström macroglobulinaemia (WM) (1-3, 8). Low and negative expression of the antigen is found in B-cell chronic lymphocytic/small lymphocytic lymphoma (B-CLL/SLL), null-cell acute lymphoblastic leukaemia (Null-ALL), B-cell acute lymphoblastic leukaemia (B-ALL) and multiple myeloma (MM) (1-4, 8). The antigen is not expressed on acute myeloid leukaemia (AML), chronic myeloid leukaemia (CML), T- ALL and T-PLL (4).						
<b>Reagent provided</b>	The Anti-B Cell, FMC7 conjugate, F 7110, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 <sup>6</sup> leucocytes from normal human peripheral blood). <u>Isotype:</u> IgM, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.						
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 7110</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0934</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 7110	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0934
Antibody Code No.	Fluorochrome	Negative Control Code No.					
F 7110	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0934					
<b>Immunogen</b>	Human B lymphoblastoid cell line HRIK (5).						
<b>Specificity</b>	Anti-B Cell, FMC7, reacts with 3-6% of normal peripheral blood mononuclear cells. These represent a subset of B cells that also expresses surface membrane immunoglobulin (Smlg). The antibody does not stain B cells without Smlg, T cells and monocytes (4). In chronic lymphoproliferative disorders, Anti-B Cell, FMC7, has been shown to label B-PLL (2/2 cases (1)), HCL (14/14 cases (1)), MCL (15/15 cases (1) and 12/13 cases (2)), FCL (24/24 cases (1) and 59/62 cases (2)), MZL (20/20 cases (1) and 13/14 cases (2)), MALT (14/17 cases (2)), DLCL (46/64 cases (2)), BL (2/2 cases (1) and 1/1 case (2)) and WM (5/6 cases (1) and 10/11 cases (8)). Anti-B Cell, FMC7, shows weak or no reactivity to B-CLL/SSL (27/27 cases (2)), 105/121 cases (1); and 17/17 cases (4)), Null-ALL (19/19 cases (4)), B-ALL (1/1 case (2)), T-ALL (5/5 cases (4)), AML/CML (6/6 cases (4)), T-PLL (1/1 case (4)) and MM (5/5 cases (8)).						
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.						
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.						
<b>Staining procedure</b>	1. Collect venous blood into a test tube containing an anticoagulant. 2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. 3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. 4. Mix 100 µL cell suspension with 10 µL F 7110.						
(103507-001)	F 7110/EFG/SSA/01.03.03 p. 1/4						

- Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
- Incubate in the dark at 4 °C for 30 minutes.
- Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
- Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
1. Ahmad E, Garcia D, Davis BH. Clinical utility of CD23 and FMC7 antigen coexistent expression in B-cell lymphoproliferative disorder subclassification. Cytometry 2002;50:1-7.
2. Garcia DP, Rooney MT, Ahmad E, Davis BH. Diagnostic usefulness of CD23 and FMC-7 antigen expression patterns in B-cell lymphoma classification. Am J Clin Pathol 2001;115:258-65.
3. Zola H, Neoh SH, Potter A, Melo JV, De Oliveria MSP, Catovsky D. Markers of differentiated B cell leukaemia: CD22 antibodies and FMC7 react with different molecules. Dis Markers 1987;5:227-35.
4. Brooks DA, Beckman IGR, Bradley J, McNamara PJ, Thomas ME, Zola H. Human lymphocyte markers defined by antibodies derived from somatic cell hybrids. IV. A monoclonal antibody reacting specifically with a subpopulation of human B lymphocytes. J Immunol 1981;126:1373-7.
5. Ferro LM, Zola H. Modulation of expression of the antigen identified by FMC7 upon human B-lymphocyte activation: evidence for differences between activation *in vivo* and *in vitro*. Immunology 1990;69:373-8.
6. Zola H, Moore HA, Hohmann A, Hunter IK, Nikoloutsopoulos A, Bradley J. The antigen of mature human B cells detected by the monoclonal antibody FMC7: studies on the nature of the antigen and modulation of its expression. J Immunol 1984;133:321-6.
7. Serke S, Schwane I, Yordanova M, Szczepek A, Huhn D. Monoclonal antibody FMC7 detects a conformational epitope on the CD20 molecule: evidence from phenotyping after rituxan therapy and transfectant cell analyses. Cytometry 2001;46:98-104.
8. San Miguel JF, Caballero MD, Gonzalez M, Zola H, Lopez Borrascas A. Immunological phenotype of neoplasms involving the B cell in the last step of differentiation. Br J Haematol 1986;62:75-83.

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