## Monoclonal Mouse Anti-Human CD66abce, Clone Kat4c

**Code No./ Code/ Code-Nr. F 7112**      **FITC-Conjugated**  
**Code No./ Code/ Code-Nr. C 7113**      **RPE-Cy5-Conjugated**

ENGLISH											
<b>Intended use</b>	For in vitro diagnostic use. F 7112 and C 7113 are intended for use in flow cytometry. Anti-CD66abce is considered essential for the initial evaluation of adult cases of acute myeloid leukaemia (AML) and chronic myelomonocytic leukaemia (CML), and for detecting minimal residual disease in acute lymphoblastic leukaemia (ALL) together with a panel of other antibodies (1-3). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.										
<b>Synonym for antigen</b>	NCA-160, biliary glycoprotein (BGP), CD67, CGM6, NCA-95, NCA, NCA-50/90, and CEA (4).										
<b>Introduction</b>	Antibodies against CD66 recognize the carcinoembryonic antigen (CEA) family of proteins and may be further characterized by reactivity with one or more family members indicated by a low letter after CD66: e.g. CD66a (BGP), CD66b (CGM6, CD67), CD66c (NCA), CD66d (CGM1), CD66e (CEA), CD66f (pregnancy-specific glycoprotein (PSG)) (4, 5). The human CEA gene family consists of 22 closely related genes belonging to the CEA superfamily that is located on chromosome 19. The gene family can be divided into two main subgroups, with 9 belonging to the CEA subgroup and 13 to the PSG subgroup (6). The molecular weight of the different CD66 types varies between 35 kDa and 200 kDa: CD66a (140-180 kDa), CD66b (95-100 kDa), CD66c (90 kDa), CD66d (35 kDa), CD66e (180-200 kDa), and CD66f (found in serum 54-72 kDa) (4). In peripheral blood samples CD66a, CD66b, CD66c and CD66d are only expressed on granulocytes. CD66a, CD66c and CD66e are expressed on epithelial cells as well (4). Anti-CD66a, CD66abce, CD66ace, CD66acd, CD66acde, and CD66c show binding to lymphoid cell lines of B and T cell origin, childhood ALLs, CML, and AML (3, 7). Anti-CD66abce has in conjunction with other antibodies been found useful for initial evaluation of adult cases of AML and chronic myeloproliferative disorders, and for detecting minimal residual disease in ALL (1-3). The function of the different CD66 molecules is unclear but they are capable of homophilic and heterophilic adhesion, and of activating neutrophils (4, 6).										
<b>Reagent provided</b>	The Anti-CD66abce conjugates, F 7112 and C 7113, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 <sup>6</sup> leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.										
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 7112</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> <tr> <td>C 7113</td> <td>RPE-Cy5 (R-Phycoerythrin-Cyanine 5)</td> <td>X 0955</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 7112	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927	C 7113	RPE-Cy5 (R-Phycoerythrin-Cyanine 5)	X 0955	
Antibody Code No.	Fluorochrome	Negative Control Code No.									
F 7112	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927									
C 7113	RPE-Cy5 (R-Phycoerythrin-Cyanine 5)	X 0955									
<b>Specificity</b>	Anti-CD66abce, Kat4c, reacts with a, b, c and e transfectants from the 6 subgroups of the carcinoembryonic antigen (CEA) group within the Ig superfamily, and is therefore classified as CD66abce (5). Anti-CD66abce, Kat4c labels mature myelopoietic cells and maturing myelopoietic cells in normal bone marrow (1, 7). Flow cytometric analysis of bone marrow samples using Anti-CD66abce, Kat4c, anti-CD13 and anti-CD14 allows for division of myelopoietic cells into three subsets: CD13+CD14-CD66abce- cells (immature myelopoietic cells), CD13+CD14-CD66abce+ cells (maturing myelopoietic cells), and CD13-CD14-CD66abce+ cells (mature myelopoietic cells) (1). This can be useful in the evaluation of patients with malignant myeloid disorders (1). Flow cytometric analysis of bone marrow samples demonstrates that Anti-CD66abce, Kat4c labels AML (2/29 cases) and ALL (8/12 cases) (3). Anti-CD66abce, Kat4c labels 19/108 cases of AML and 9/20 cases of chronic myeloproliferative disorder show CD66abce and CD14 co-expression (2).										
<b>Precautions</b>	<ol style="list-style-type: none"> <li>For professional users.</li> <li>This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>As with any product derived from biological sources, proper handling procedures should be used.</li> </ol>										
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.										
<b>Staining procedure</b>	<ol style="list-style-type: none"> <li>Collect venous blood into a test tube containing an anticoagulant.</li> <li>Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.</li> <li>Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.</li> <li>Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD66abce.</li> <li>Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).</li> <li>Incubate in the dark at 4 °C for 30 minutes.</li> <li>Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 0.01 mol/L PBS, pH 7.4.</li> <li>Analyse on a flow cytometer.</li> </ol>										
(102395-001)	F 7112/C 7113/EFG/SSA/01.03.03 p. 1/4										

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

**Product-specific limitations**      Additionally, it has been observed that RPE-Cy5-conjugates may bind to monocytes, resulting in background staining (8).

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
- Bonde J, Meyer K, Broe MK, Hokland M, Turley H, Hokland P. Improved flow cytometric identification of myelopoiesis by the simultaneous labeling with CD13, CD14 and CD66 monoclonal antibodies. Br J Haematol 1996;92:269-79.
- Hansen I, Meyer K, Hokland P. Flow cytometric identification of myeloid disorders by asynchronous expression of the CD14 and CD66 antigens. Eur J Haematol 1998;61:339-46.
- Carrasco M, Munoz L, Bellido M, Bernat S, Rubiol E, Ubeda J, et al. CD66 expression in acute leukaemia. Ann Hematol 2000;79:299-303.
- Turni L, Shaw S, Watson B, Mason D. CD guide. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, De Haas M, et al, editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002 .p. 816-21.
- Skubitz KM, Micklem K, van der Schoot E. M14. CD66 and CD67 cluster workshop report. In: Schlossman SF, Bousmell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume One .p. 889-99.
- Thompson JA, Grunert F, Zimmermann W. Carcinoembryonic antigen gene family: molecular biology and clinical perspectives. J Clin Lab Anal 1991;5:344-66.
- Hanenberg H, Quentin I, Baumann M, Grunert F, Burdach S. M14.8. Expression of CEA gene family members on haemopoietic cells. In: Schlossman SF, Bousmell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume One. p. 919-21.
- van Vugt MJ, van den Herik-Oudijk IE, van de Winkel JGJ. Binding of PE-CY5 conjugates to the human high-affinity receptor for IgG (CD64). Blood 1996;88:2358-61.

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