Monoclonal Mouse
Anti-Human CD24/FITC
 Clone SN3
Code No./ Code/ Code-Nr. F 7134
 Edition/ Ausgabe 01.01.03

ENGLISH							
Intended use	For in vitro diagnostic use. F 7134 is intended for use in flow cytometry. Anti-CD24 is considered essential for the initial evaluation of non-T-cell acute lymphoblastic leukaemias (ALL) and B-cell lymphoproliferative disorders together with a panel of other antibodies (1-5). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.						
Synonym for antigen	BA-1 and heat stable antigen (HSA) (6).						
Introduction	The CD24 antigen, which has a molecular weight of 35-45 kDa, is a glycosyl phosphatidylinositol (GPI)-anchored surface molecule. It is a member of a family that also includes HSA and CD52. Cross-linking of CD24 can induce increased intracellular Ca ⁺⁺ in mature B lymphocytes (6). CD24 is expressed at multiple stages of B cell development, beginning with the bone marrow CD34+ pro-B-cell compartment and continuing through mature surface IgM+/IgD+ B cells. Plasma cells express very low to negligible levels of CD24 (7). CD24 is not expressed on the cell surface of normal T cells or monocytes, but is expressed on mature granulocytes (4, 7). The antigen is also present on a variety of epithelial types (6). CD24 is expressed on the vast majority of B-ALL, B-cell chronic lymphocytic leukaemia (B-CLL) and B-cell non-Hodgkin's lymphoma (1-5). The antigen is also present on acute myeloid leukaemias of FAB subtype M4 and M5 and on CD13+/CD33+ myeloid cells in the bone marrow of patients with chronic myeloid leukaemia (CML) (3, 7). CD24 is present at high levels in small cell lung carcinoma and in ovarian cancer (6, 8).						
Reagent provided	The Anti-CD24 conjugate, F 7134, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L Na ₃ N, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.						
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 7134</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 7134	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927
Antibody Code No.	Fluorochrome	Negative Control Code No.					
F 7134	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927					
Immunogen	Isolated membrane fraction of the pre-B leukaemic cell line NALM-1 (5).						
Specificity	Anti-CD24, SN3, was included in the Fifth International Workshop and Conference on Leucocyte Cell Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with the CD24 antigen (7). Anti-CD24, SN3, reacts with B cells and granulocytes (4, 7). It shows no reactivity to T cells, monocytes, erythrocytes or platelets. In acute leukaemias, Anti-CD24, SN3, has been shown to label 17/23 of non-T/non-B-cell ALL cases including pre-B-ALL cases, 1/1 case of B-ALL, 1/5 of acute myelocytic leukaemia and 1/3 cases of acute monocytic leukaemia (4). In chronic lymphoproliferative disorders, Anti-CD24, SN3, has been shown to label 6/6 cases of B-CLL, 1/1 case of B-cell lymphoma (4), 16/20 cases of B-cell lymphoma (5) and 2/2 cases of Hodgkin's disease. The antibody also labelled 5/5 cases of CML (4).						
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.						
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.						
Staining procedure	1. Collect venous blood into a test tube containing an anticoagulant. 2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. 3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. 4. Mix 100 µL cell suspension with 10 µL F 7134. 5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). 6. Incubate in the dark at 4 °C for 30 minutes.						

- Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
- Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
- van Dongen JJM, Adriaansen HJ. Immunobiology of leukemia. In: Henderson ES, Lister TA, Greaves MF, editors. Leukemia. Philadelphia, London, Toronto, Montreal, Sydney, Tokyo: WB Saunders Company; 1996. p. 83-130.
- Ratei R, Sperling C, Karawajew L, Schott G, Schrappe M, Harbott J, Riehm H, et al. Immunophenotype and clinical characteristics of CD45-negative and CD45-positive childhood acute lymphoblastic leukemia. Ann Hematol 1998;77:107-14.
- Raife TJ, Lager DJ, Kemp JD, Dick FR. Expression of CD24 (BA-1) predicts monocytic lineage in acute myeloid leukemia. Am J Clin Pathol 1994;101:296-9.
- Fukukawa T, Matsuzaki H, Haruta Y, Hara H, Seon BK. New monoclonal antibodies SN3, SN3a, and SN3b directed to sialic acid of glycoprotein on human non-T leukemia cells. Exp Hematol 1986;14:850-5.
- Barcos M, Pollard C, Fukukawa T, Seon BK. Follicular mantle zone cell subpopulations detected by monoclonal antibody SN3. Hematol Oncol 1986;4:251-9.
- Turni L, Shaw S, Watson B, Mason D. CD guide. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 769-70.
- Le Bien TW. B8. CD24 workshop panel report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the Fifth International Workshop and Conference; 1993 Niv 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. Volume 1 p. 539-41.
- Kristiansen G, Denkert C, Schlüns K, Dahl E, Pilarsky C, Hauptmann S. CD24 is expressed in ovarian cancer and is a new independent prognostic marker of patient survival. Am J Pathol 2002;161:1215-21.

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