**Monoclonal Mouse
Anti-Human CD68/FITC**

Clone KP1

Code No./ Code/ Code-Nr. F 7135
ENGLISH

Intended use For in vitro diagnostic use.
F 7135 is intended for use in flow cytometry. Anti-CD68, Clone KP1, is a monocytes/macrophage associated antigen and can be used, together with a panel of other antibodies, to discriminate cells of lymphoid and monocytes/macrophage origin (1, 2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Synonyms for antigen Macrosialin, gp110 (8).

Introduction CD68 is a 110 kDa highly glycosylated lysosomal membrane protein. The CD68 protein belongs to a family of lysosomal glycoproteins/plasma membrane shuttling proteins that play a role in endocytosis and/or lysosomal trafficking (3, 7). Unlike many other CD leucocyte antigens, the CD68 molecule is antigenically very heterogeneous, and different antibodies to CD68 show different cellular reactivities (4).

CD68 is expressed strongly in cytoplasmic granules, and weakly on the surface of macrophages, monocytes, neutrophils, basophils and NK-cells (3). CD68 is also found in neoplastic cells of myeloid background. Additionally, CD68 is expressed by approximately 40% of peripheral blood B cells and is weakly expressed in 50% of B cell type acute lymphoblastic leukaemia cells (3). CD68 can also be found in the cytoplasm of non-haematopoietic tissues, especially the liver, glomeruli and renal tubules (3).

Reagent provided The Anti-CD68 conjugate, F 7135, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).

Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
F 7135	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927

Immunogen Lysosomal fraction of human lung macrophages (6).

Specificity Anti-CD68, KP1, was clustered as anti-CD68 at the Fourth International Workshop and Conference on Human Leucocyte Differentiation Antigens held in Vienna in 1989 (7).

Anti-CD68, KP1, labels the cytoplasm in type M1-M5 AML, giving a strong granular staining (1, 4).

Mast cells, both normal and neoplastic, and 1/1 true histiocytic neoplasia, were labelled by the antibody (2, 5).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

Storage Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Staining procedure

1. Transfer 50 µL (up to 10⁶ cells) of the cell suspension to be analysed (whole blood, bone marrow or mononuclear cells) to each test tube.
2. Add 100 µL Dako IntraStain, Reagent A (Fixation), code no. K 2311. Mix gently by using a vortex mixer to ensure that cells are in suspension.
3. Incubate at room temperature for 15 minutes.
4. Add 2 mL PBS and mix gently by using a vortex mixer.
5. Centrifuge at 300 x g for 5 minutes, then aspirate the supernatant, leaving approximately 50 µL of fluid.
6. Mix thoroughly by using a vortex mixer to ensure that cells are in suspension and add 100 µL Dako IntraStain, Reagent B (Permeabilisation), code no. K 2311. Add 10 µL F 7135. Mix gently by using a vortex mixer to ensure that the cells are in suspension.
7. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
8. Incubate in the dark at room temperature for 15 minutes.
9. Repeat steps 4 and 5.

10. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
11. Analyse on a flow cytometer

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
1. Horny H-P, Wehrmann M, Steinke B, Kaiserling E. Assessment of the value of immunohistochemistry in the subtyping of acute leukemia on routinely processed bone marrow biopsy specimens with particular reference to macrophage-associated antibodies. *Hum Pathol* 1994;25:810-4.
2. Warnke RA, Pulford KAF, Pallesen G, Ralfkiaer E, Brown DC, Gatter KC, et al. Diagnosis of myelomonocytic and macrophage neoplasms in routinely processed tissue biopsies with monoclonal antibody KP1. *Am J Pathol* 1989;135:1089-95.
3. Goyert SM. MC12. CD68 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc; 1997. p. 1015-6.
4. Falini B, Flenghi L, Pileri S, Gambacorta M, Bigerna B, Durkop H, et al. PG-M1: a new monoclonal antibody directed against a fixative-resistant epitope on the macrophage-restricted form of the CD68 molecule. *Am J Pathol* 1993;142:1359-72.
5. Horny H-P, Schaumburg-Lever G, Bolz S, Geerts ML, Kaiserling E. Use of monoclonal antibody KP1 for identifying normal and neoplastic human mast cells. *J Clin Pathol* 1990;43:719-22.
6. Pulford KAF, Rigney EM, Micklem KJ, Jones M, Stross WP, Gatter KC, et al. KP1: a new monoclonal antibody that detects a monocyte/macrophage associated antigen in routinely processed tissue sections. *J Clin Pathol* 1989;42:414-21.
7. Micklem K, Cordell J, Rigney E, Simmons D, Pulford K, Stross P, et al. M13.1. A macrophage-associated antigen defined by five mAB. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 843-6.
8. Goyert SM. CD Guide. CD68. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 1166.

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Agilent Technologies Singapore (International) Pte Ltd.