**Monoclonal Mouse Anti-Human CD79 $\beta$ /FITC, Clone SN8**  
**Monoclonal Mouse Anti-Human CD79 $\beta$ /RPE, Clone SN8**

**Code No. F 7137**  
**Code No. R 7272**

**ENGLISH**

**Intended use** For in vitro diagnostic use.  
 F 7137 and R 7272 are intended for use in flow cytometry. F 7137 and R 7272 are intended for use in identification of cells expressing CD79 $\beta$  and react with an epitope on the extracellular portion of the  $\beta$  chain of CD79 (1). In flow cytometry, anti-CD79 $\beta$ , together with a panel of other antibodies, is useful for the immunophenotyping of leukaemias and lymphomas (2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

**Synonym for antigen** Ig- $\beta$ , B29 (3).

**Introduction** Together with the CD79 $\alpha$  polypeptide, the CD79 $\beta$  polypeptide forms a disulphide-linked, transmembrane heterodimer, which is non-covalently associated with membrane-bound immunoglobulins on B cells (4). These 3 components constitute the B-cell antigen receptor complex. Antigen-binding by the immunoglobulin in this complex leads to CD79 $\alpha$ /CD79 $\beta$ -mediated signal transduction resulting in activation of multiple biochemical pathways. CD79 $\alpha$  and CD79 $\beta$  have a molecular mass of 32-33 kDa and 37-39 kDa, respectively, and their cytoplasmic determinants are named CD79 $\alpha$ cy and CD79 $\beta$ cy (3). CD79 $\beta$  shows high identity (96%) between man and mouse in the cytoplasmic and transmembrane domains, whereas only 59% homology is found in the extracellular domain (2).  
 CD79 $\beta$  is restricted to B-cell lineage, and is found in pro-B cells before Ig heavy-chain gene rearrangement but after CD79 $\alpha$  expression. In contrast to CD79 $\alpha$ , CD79 $\beta$  expression persists to the latest stage of B-cell development (1, 2), and furthermore, CD79 $\beta$  expression increases upon B-cell activation (2).  
 CD79 $\beta$  is expressed in B-cell neoplasms, with an incidence of approximately 80%, with the exception of CLL/SLL (chronic lymphocytic leukaemia/small lymphocytic lymphoma) where CD79 $\beta$  is negative in the majority of the cases (2).

**Reagent provided** The Anti-CD79 $\beta$  conjugates, F 7137 and R 7272, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 100 tests (10  $\mu$ L of conjugate for up to 10<sup>6</sup> leucocytes from normal human peripheral blood).  
Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
F 7137	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927
R 7272	RPE (R-Phycoerythrin)	X 0928

**Immunogen** Membrane antigen preparation from human B prolymphocytic leukaemia cells (5).

**Specificity** Anti-CD79 $\beta$ , SN8, was included in the Fifth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with the CD79 $\beta$  (3, 4). SN8 reacts with an epitope on the extracellular portion of the  $\beta$  chain of the human CD79 antigen (1).

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

**Storage** Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Staining procedure**

1. Collect venous blood into a test tube containing an anticoagulant.
2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
4. Mix 100  $\mu$ L cell suspension with 10  $\mu$ L Anti-CD79 $\beta$ .
5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
6. Incubate in the dark at 4 °C for 30 minutes.
7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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
1. Yoshida M, Vasile S, Coligan JE, Okazaki M. B23.3. mAb SN8 defines a unique extracellular epitope of the B29/mb1 heterodimer of the human B-cell antigen receptor complex. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leukocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 674-5.
2. Chu PG, Arber DA. CD79: A review. Appl Immunohistochem Mol Morphol 2001;9:97-106.
3. Engel P, Wagner N, Tedder TF. B23. CD79 workshop report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 667-70.
4. Persin C, Jensen WA, Cambier JC. B23.1. Definition of the submolecular specificity of mAb against mb1 (Igα) and B29 (Igβ/γ). In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leukocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 670-3.
5. Okazaki M, Luo Y, Han T, Yoshida M, Seon BK. Three new monoclonal antibodies that define a unique antigen associated with prolymphocytic leukemia/non-Hodgkin's lymphoma and are effectively internalized after binding to the cell surface antigen. Blood 1993;81:84-94.

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Agilent Technologies Singapore (International) Pte Ltd.