## Monoclonal Mouse Anti-Human CD103, Mucosa Lymphocyte Antigen, Clone Ber-ACT8

Code No./ Code/ Code-Nr. F 7138      FITC-Conjugated  
 Code No./ Code/ Code-Nr. R 7188      RPE-Conjugated

### ENGLISH

**Intended use** For in vitro diagnostic use.  
 F 7138 and R 7188 are intended for use in flow cytometry. CD103 is expressed on mucosal T cells, on activated CD8+ T cells and on hairy cell leukaemia cells. In several malignant conditions, such as enteropathy-associated T-cell lymphomas and hairy cell leukaemia, the cells express CD103. The antibody is well-suited for the immunophenotyping of leukaemias and lymphomas (1, 2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

**Synonyms for antigen** ITGAE, Human mucosal lymphocyte antigen 1 (HML-1), Integrin  $\alpha$ E chain (3).

**Introduction** CD103 is the  $\alpha^E$  integrin subunit of the heterodimeric  $\alpha^E\beta 7$  integrin belonging to the small  $\beta 7$  integrin subfamily. Antibodies against CD103 immunoprecipitate an antigenic protein composed of two chains of 105 kDa ( $\beta 7$ ) and 175 kDa ( $\alpha^E$ ) under non-reduced conditions. Upon reduction, the 175 kDa chain dissociates into two fragments of 150 and 25 kDa (4).

**Reagent provided** The Anti-CD103 conjugates, F 7138 and R 7188, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L  $\text{NaN}_3$ , pH 7.2. Each vial contains 100 tests (10  $\mu\text{L}$  of conjugate for up to  $10^6$  PHA-stimulated lymphocytes isolated from normal human peripheral blood).

Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
F 7138	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927
R 7188	RPE (R-Phycoerythrin)	X 0928

**Immunogen** HTLV-1+ MAPS-16 cell line (1).

**Specificity** Anti-CD103, Ber-ACT8, was included in the Fifth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with the CD103 antigen (4). The antibody reacts with an epitope localized on the 150 kDa polypeptide chain, and detects the same molecule target as monoclonal antibodies HML-1, LF61 and B-Ly7 (1). In the gut, the antibody labels more than 90% of intestinal intraepithelial lymphocytes and about 40% of lamina propria T lymphocytes (3). In other normal human tissues the reactivity of the antibody is restricted to a few (0.5-5%) resting CD8+ T cells in peripheral blood, lymph nodes and tonsil it stains some scattered small lymphoid cells present in the T zones and the follicle mantels. In the red pulp of the spleen it stains 8-10% of cells (1, 3).

**Precautions**

1. For professional users.
2. This product contains sodium azide ( $\text{NaN}_3$ ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

**Storage** Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Staining procedure**

1. Collect venous blood into a test tube containing an anticoagulant.
2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
4. Mix 100  $\mu\text{L}$  cell suspension with 10  $\mu\text{L}$  fluorochrome-conjugated Anti-CD103.
5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
6. Incubate in the dark at 4 °C for 30 minutes.
7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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Agilent Technologies Singapore (International) Pte Ltd.