**Monoclonal Mouse Anti-Human CD1a/FITC, Clone NA1/34**  
**Monoclonal Mouse Anti-Human CD1a/RPE, Clone NA1/34**

**F7141**  
**R7189**

ENGLISH										
<b>Intended use</b>	For in vitro diagnostic use. F7141 and R7189 are intended for use in flow cytometry. Antibodies to CD1a may be a useful tool for identification of Langerhans' cell histiocytosis and for the classification of thymomas and malignancies of T-cell precursors (1). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.									
<b>Synonym for antigen</b>	R4 (2).									
<b>Introduction</b>	The human CD1 locus contains five genes, CD1A-E. CD1 molecules consist of glycosylated polypeptid $\alpha$ chains that span the cell membrane and are non-covalent associated with $\beta$ 2-microglobulin (3). Anti-CD1a recognizes a heavy chain molecule of 49 kDa. This antigen is expressed strongly on cortical thymocytes and on a variety of antigen-presenting cells including Langerhans' cells and interdigitating dendritic cells. CD1a is not expressed on early thymocytes, and is absent on mature peripheral blood T cells but intracytoplasmic expression is detected on activated T lymphocytes (1-3).									
<b>Reagent provided</b>	The Anti-CD1a conjugates, F7141 and R7189, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2. Each vial contains 100 tests (10 $\mu$ L of conjugate for up to 10 <sup>6</sup> JM cells). <u>Isotype</u> : IgG2a, kappa. <u>Conjugate concentration mg/L</u> : See label on vial.									
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Control Reagent Code No.</th> </tr> </thead> <tbody> <tr> <td>F7141</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X0933</td> </tr> <tr> <td>R7189</td> <td>RPE (R-Phycoerythrin)</td> <td>X0950</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Control Reagent Code No.	F7141	FITC (Fluorescein Isothiocyanate Isomer 1)	X0933	R7189	RPE (R-Phycoerythrin)	X0950
Antibody Code No.	Fluorochrome	Control Reagent Code No.								
F7141	FITC (Fluorescein Isothiocyanate Isomer 1)	X0933								
R7189	RPE (R-Phycoerythrin)	X0950								
<b>Immunogen</b>	Human thymocytes (4).									
<b>Specificity</b>	Anti-CD1a, NA1/34, was included in the First, Second, Third and Fourth International Workshop and Conferences on Human Leucocyte Differentiation Antigens, and studies by different laboratories confirmed its reactivity with the CD1a antigen (2).									
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.									
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.									
<b>Staining procedure</b>	<ol style="list-style-type: none"> <li>Transfer 100 <math>\mu</math>L of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.</li> <li>Add 10 <math>\mu</math>L of fluorochrome-conjugated anti-CD1a and mix gently by using a vortex mixer.</li> <li>Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes.</li> <li>Add 100 <math>\mu</math>L of Uti-Lyse™ Reagent A (DakoCytomation code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.</li> <li>Add 1 mL of Uti-Lyse™ Reagent B (DakoCytomation code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.</li> <li>Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 <math>\mu</math>L of fluid in the tube.</li> <li>Add 2 mL of PBS (DakoCytomation code No. S3024) to the tube and resuspend the cells by using a vortex mixer.</li> <li>Repeat step 6.</li> <li>Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.</li> <li>Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining.</li> </ol> <p>Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.</p> <p><b>Procedural notes</b></p> <p>Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.</p> <p>Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.</p> <p>It is optional to include a control reagent test tube. The control reagent should match the isotype and fluorochrome of the conjugated antibody. Recommended control reagents are shown in the table above.</p> <p>Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, such as DakoCytomation EasyLyse™, code No. S2364, the PBS in step 9 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.</p>									

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis.

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


1. Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 2003. p. 59-60.
2. Fainboim L, Salamone MC. CD guide. CD1a-e. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 747-8.
3. Fainboim L, Salamone MC. TC2. CD1 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 33-7.
4. McMichael AJ, Pilch JR, Galfré G, Mason DY, Fabre JW, Milstein C. A human thymocyte antigen defined by a hybrid myeloma monoclonal antibody. Eur J Immunol 1979;9:205-10.

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