Monoclonal Mouse
Anti-Human Plasma Cell/FITC
 Clone VS38c
Code No./ Code/ Code-Nr. F 7149
 Edition/ Edition/ Ausgabe 01.11.02

ENGLISH							
Intended use	For in vitro diagnostic use. F 7149 is intended for use in flow cytometry. The antibody reacts strongly with human plasma cells (1). The VS38c antibody has been shown to identify the human p63 protein (2, 3). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.						
Synonym for antigen	p63 (2).						
Introduction	Anti-Human Plasma Cell recognizes the intracytoplasmic antigen p63 of 64 kDa present in normal and neoplastic cells (4). p63 is an intracellular type II transmembrane protein which is localized in the rough endoplasmic reticulum. The function of p63 is unknown, but its abundance in secretory cells and localization to the endoplasmic reticulum suggests a conserved role in protein processing or secretion (2).						
Reagent provided	The Anti-Human Plasma Cell conjugate, F7149, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 ⁶ HS-Sultan cells). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.						
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 7149</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 7149	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927
Antibody Code No.	Fluorochrome	Negative Control Code No.					
F 7149	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927					
Immunogen	MCF-7 AdR, a multi-drug resistant breast carcinoma cell line (1).						
Specificity	Anti-Human Plasma Cell, VS38c, was submitted to and characterized at the sixth International Workshop and Conference on Human Leucocyte Differentiation Antigen (2). Intracellular markers are not clustered. Anti-Human Plasma Cell, VS38c, labels normal and neoplastic plasma cells in various tissues, but does not label peripheral blood lymphocytes, although some weak labelling of monocytes and neutrophils could occasionally be detected (1). The VS38c antibody also labels melanocytic cells, particularly malignant melanoma cells, a number of epithelial cells from various organs (5) and osteoblasts from bone marrow biopsies (4). Myeloma, plasmacytoma, and lymphoplasmacytoid lymphomas are strongly labelled and about a third of large B-cell lymphomas are weakly labelled. T-cell lymphomas are negative. Occasional weak labelling of endothelial cells may occur (2).						
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. 						
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.						
Staining procedure	<ol style="list-style-type: none"> Transfer 50 µL (up to 10⁶ cells) of the cell suspension to be analysed (whole blood, bone marrow or mononuclear cells) to a test tube. Add 100 µL DakoCytomation IntraStain Reagent A (Fixation), code No. K 2311. Mix gently by using a vortex mixer to ensure that the cells are in suspension. Incubate at room temperature for 15 minutes. Add 2 mL PBS and mix gently by using a vortex mixer. Centrifuge at 300 x g for 5 minutes, then aspirate the supernatant, leaving approximately 50 µL of fluid. Mix thoroughly by using a vortex mixer to ensure that the cells are in suspension and add 100 µL DakoCytomation IntraStain Reagent B (Permeabilization), code No. K 2311. Add 10 µL of F 7149. Mix gently by using a vortex mixer to ensure that the cells are in suspension. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). Incubate in the dark at room temperature for 15 minutes. Repeat steps 4 and 5. 						

- Resuspend the pellet in an appropriate fluid for flow cytometric analysis, e.g. 0,3 mL 1% paraformaldehyde (fixative) in 0,01 mol/L PBS, pH 7,4.
- Analyse on a flow cytometer

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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