Monoclonal Mouse Anti-Human HLA-DR Antigen/FITC, Clone AB3
Monoclonal Mouse Anti-Human HLA-DR Antigen/RPE, Clone AB3
Code No./ Code/ Code-Nr. F 7266
Code No./ Code/ Code-Nr. R 7267

ENGLISH										
Intended use	For in vitro diagnostic use. F 7266 and R 7267 are intended for use in flow cytometry. Antibodies to HLA-DR are considered essential for the initial evaluation of acute leukaemia, chronic T- and B-cell leukaemia and myeloid leukaemia when used together with a panel of other antibodies (1-3). F7266 and R7267 are not intended for tissue typing. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.									
Introduction	The main function of human leucocyte antigen (HLA) molecules is to present antigenic peptides to the T-cell receptor, thereby regulating the induction of the immune response (4). The HLA molecules are encoded by a cluster of tightly linked genes located on the short arm of chromosome 6. Three classes of HLA molecules (I, II and III) have been denoted. Human class II genes are located in the HLA-D region, consisting of three families called <i>DQ</i> , <i>DP</i> and <i>DR</i> (4). The products of class II genes form a heterodimeric transmembrane protein, consisting of a heavy (~34 kDa) α -chain and a light (~28 kDa) β -chain (4). The DR α -chain is expressed from one non-polymorphic gene, whereas the DR β -chain originates from nine highly polymorphic genes (4). HLA-DR antigen is constitutively expressed on antigen-presenting cells, such as B lymphocytes, monocytes and dendritic cells but can also be detected on activated T lymphocytes and activated granulocytes (4, 5). Occasionally, natural killer cells express HLA-DR antigen (1). The antigen has been found expressed in cases of different types of acute lymphoblastic leukaemias, acute myeloid leukaemias except AML-M3, chronic lymphoblastic leukaemias, chronic myeloid leukaemias and B- and T-cell non-Hodgkin's leukaemias (1-3, 6, 7). However, the antigen is normally not present on non-haematopoietic tumours and multiple myelomas (6).									
Reagent provided	The Anti-HLA-DR Antigen conjugates, F 7266 and R 7267, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 100 tests (10 μ L of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood). <u>Isotype:</u> IgG2a, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.									
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 7266</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0933</td> </tr> <tr> <td>R 7267</td> <td>RPE (R-Phycoerythrin)</td> <td>X 0950</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 7266	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0933	R 7267	RPE (R-Phycoerythrin)	X 0950
Antibody Code No.	Fluorochrome	Negative Control Code No.								
F 7266	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0933								
R 7267	RPE (R-Phycoerythrin)	X 0950								
Immunogen	B cells from a follicular lymphoma.									
Specificity	Anti-Human HLA-DR Antigen, AB3, reacts with the DRw52 determinant of the HLA-DR molecule. The specificity has been tested with transfectants expressing HLA-DRw52. Anti-Human HLA-DR Antigen, AB3, does not react with transfectants expressing HLA-DP antigen and HLA-DQ antigen. The antibody has been submitted to the 8th International Workshop on Human Leucocyte Differentiation Antigens 2000-2004.									
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. 									
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.									
Staining procedure	<ol style="list-style-type: none"> Collect venous blood into a test tube containing an anticoagulant. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4. Mix 100 μL cell suspension with 10 μL fluorochrome-conjugated Anti-HLA-DR Antigen. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table). Incubate in the dark at 4 °C for 30 minutes. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4. Analyse on a flow cytometer. 									
(106537-002)	F 7266/R 7267/EFG/TJA/27.11.03 p. 1/4									

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

The escapee phenomenon has been described for FITC-conjugated HLA-DR antigen antibodies when using ammonium chloride (NH₄Cl)-based lysing reagents. The phenomenon can be reduced by inclusion of paraformaldehyde in the lysing reagent (8).

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
- Van Dongen JJM, Adriaansen HJ. Immunobiology of leukaemia. In: Henderson ES, Lister TA, Greaves MF, editors. Leukemia. Philadelphia, London, Toronto, Montreal, Sydney, Tokyo: WB Saunders Company; 1996. p. 83-130.
- Mittelman M, Karcher DS, Kammerman LA, Lessin LS. High Ia (HLA-DR) and low CD11b (Mo1) expression may predict early conversion to leukemia in myelodysplastic syndromes. Am J Hematol 1993;43:165-71.
- De Zen L, Orfao A, Cazzaniga G, Masiero L, Cocito MG, Spinelli M, et al. Quantitative multiparametric immunophenotyping in acute lymphoblastic leukemia: correlation with specific genotype. I. ETV6/AML1 ALLs identification. Leukemia 2000;14:1225-31.
- Naverrete CV. The HLA system in blood transfusion. Baillière's Clin Haematol 2000;13:511-32.
- Erber WN, Pinching AJ, Mason DY. Immunocytochemical detection of T and B cell populations in routine blood smears. Lancet 1984;i:1042-5.
- Falini B, Martelli F, Tarallo F, Moir DJ, Cordell JL, Gatter KC, et al. Immunohistological analysis of human bone marrow trephine biopsies using monoclonal antibodies. Br J Haematol 1984;56:365-86.
- Smith MEF, Holgate CS, Williamson JMS, Grigor I, Quirke P, Bird CC. Major histocompatibility complex class II antigen expression in B and T cell non-Hodgkin's lymphoma. J Clin Pathol 1987;40:34-41.
- Prince HE, York J, Kuttner DK. Reduction of escapee formation in flow cytometric analysis of lymphocyte subsets. J Immunol Methods 1994;177:165-73.

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