MultiMix™

Dual-Colour Reagent

Anti-Human Kappa Light Chains/FITC

Anti-Human Lambda Light Chains/RPE

Code No./ Code/ Code-Nr. FR481

ENGLISH

Intended use

For in vitro diagnostic use.

FR481 is intended for use in flow cytometry for simultaneous detection and enumeration of kappa light chains and lambda light chains. In flow cytometry, antibodies to kappa light chains are useful for the demonstration of cell surface kappa light chains, and antibodies to lambda light chains are useful for the demonstration of cell surface lambda light chains, and, thus, for the identification of monoclonality (clonal excess) in B-cell lymphoproliferative disorders together with a panel of other antibodies (1-4). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Reagent provided

FR481 comprises the following two, carefully matched, fluorescent antibodies:

Polyclonal Rabbit Anti-Human Kappa Light Chains conjugated with fluorescein isothiocyanate isomer 1 (FITC).

Polyclonal Rabbit Anti-Human Lambda Light Chains, conjugated with R-phycoerythrin (RPE).

The two conjugates in FR481, Anti-Kappa Light Chains and Anti-Lambda Light Chains conjugate have been produced from F(ab')₂ fragments of affinity-isolated polyclonal rabbit antibodies. FR481 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 50 tests (10 µL of conjugate for up to 10⁶ leucocytes from normal human peripheral blood).

Antibody Code No.	Fluorochrome	Control Reagent Code No.
FR481	FITC and RPE	X0935

Immunogens

Rabbit Anti-Human Kappa Light Chains: Polyclonal immunoglobulin light chains of kappa type isolated from a pool of human sera.

Rabbit Anti-Human Lambda Light Chains: Polyclonal immunoglobulin light chains of lambda type isolated from a pool of human sera.

Specificity

Anti-Kappa Light Chains reacts with free kappa chains as well as kappa chains in intact immunoglobulin molecules. The specificity has been ascertained as described below. To obtain maximum sensitivity, the crossed immunoelectrophoresis specificity test was performed prior to affinity purification and pepsin degradation.

Crossed immunoelectrophoresis: Only kappa-related precipitates appear when the antibody is run against human plasma. Staining: Coomassie Brilliant Blue.

Flow cytometry: When Anti-Kappa Light Chains/FITC is applied as described in the staining procedure in combination with Anti-CD19/RPE, HD37, on lysed human whole blood, a specific staining of a part of the CD19-positive B lymphocytes is seen corresponding to the expected range of kappa light chain expression.

Anti-Lambda Light Chains reacts with free lambda chains as well as lambda chains in intact immunoglobulin molecules. The specificity has been ascertained as described below. To obtain maximum sensitivity, the crossed immunoelectrophoresis specificity test was performed prior to affinity purification and pepsin degradation.

Crossed immunoelectrophoresis: Only lambda-related precipitates appear when the antibody is run against human plasma. Staining: Coomassie Brilliant Blue.

Flow cytometry: When Anti-Lambda Light Chains/RPE is applied as described in the staining procedure in combination with Anti-CD19/FITC, HD37, on lysed human whole blood, a specific staining of a part of the CD19-positive B lymphocytes is seen corresponding to the expected range of lambda light chain expression.

Flow cytometry: When FR481 is applied as described in the staining procedure, on lysed human whole blood, a specific staining of a part of the lymphocytes is seen, within the expected range of kappa and lambda light chains expression.

Flow cytometric analysis of single suspensions from formalin-fixed, paraffin-embedded tissue samples demonstrated that Anti-Lambda Light Chains labels reactive hyperplastic lymph nodes (10/10 cases), and so does Anti-Kappa Light Chains (10/10 cases). In non-Hodgkin's lymphomas, the antibody labeled 5/10 cases. The remaining cases were positive for kappa light chains (4/10 cases) or showed no expression of light chains (1/10 cases) (4).

Flow cytometric analysis of peripheral blood lymphocytes demonstrated that Anti-Kappa Light Chains and Anti-Lambda Light Chains, respectively, label a proportion of B-cell chronic lymphocytic leukaemias. Thus, in one study of 121 cases (2), 37 were positive for kappa and 17 were positive for lambda. In another study of 165 cases (3), 97 were positive for kappa and 68 were positive for lambda.

Precautions

- For professional users.
- This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
- As with any product derived from biological sources, proper handling procedures should be used.

Storage

Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Staining procedure

Before staining samples of peripheral blood, the mononuclear cells must be isolated by centrifugation on a separation medium or the blood sample washed to remove the soluble serum proteins. As human monocytes bind serum immunoglobulins via their surface Fc receptors, these cells should be removed by depletion or identified.

- Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.
- Add 2 mL of PBS (DakoCytomation code No. S3024). Mix gently by using a vortex mixer.
- Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
- Repeat steps 2 and 3 two more times.
- Add 10 µL Rabbit Immunoglobulin Fraction (DakoCytomation code No. X0903) for blocking. Mix gently by using a vortex mixer and incubate the tube at 37 °C for 30 minutes.
- Add 10 µL of FR481 and mix gently by using a vortex mixer.
- Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes.
- Add 100 µL of Uti-Lyse™ Reagent A (DakoCytomation code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
- Add 1 mL of Uti-Lyse™ Reagent B (DakoCytomation code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
- Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.
- Add 2 mL of PBS to the tube and resuspend the cells by using a vortex mixer.

FR481/EFG/MTH/11.03.05 p. 1/4

- Repeat step 10.
- Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.
- Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Procedural notes

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 6: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the conjugated antibody reagent. The recommended Dual-Colour Control Reagent for FR481 is DakoCytomation code No. X0935.

Steps 8 and 9: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, such as DakoCytomation EasyLyse™, code No. S2364, the PBS in step 13 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.

Step 14: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis..

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