MultiMix™
Dual-Colour Reagent
Anti-Human CD45/FITC
Anti-Human CD14/RPE
Code No./ Code/ Code-Nr. FR700

ENGLISH							
Intended use	For in vitro diagnostic use. FR700 is intended for use in flow cytometry. CD45 is one of the most abundant leucocyte cell surface glycoproteins and is expressed exclusively on cells of the haematopoietic system and their progenitors (1). Anti-CD14 is important for excluding CD14 positive monocytes in relation to T lymphocyte subset analysis (2, 3). The combination of CD45 and CD14 allows simultaneous subdivision of leucocytes into lymphocytes, monocytes and granulocytes. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.						
Reagent provided	FR700 comprises the following two, carefully matched, fluorescent antibodies: Monoclonal Mouse Anti-Human CD45, Clone T29/33, conjugated with fluorescein isothiocyanate isomer 1 (FITC). Monoclonal Mouse Anti-Human CD14, Clone TÜK4, conjugated with R-phycoerythrin (RPE). The two conjugates in FR700 have been produced from purified monoclonal mouse antibodies IgG1, kappa (CD45) and IgG2a, kappa (CD14) isotypes. FR700 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN ₃ , pH 7.2. Each vial contains 50 tests (10 µL of conjugate for up to 10 ⁶ leucocytes from normal human peripheral blood).						
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Control Reagent Code No.</th> </tr> </thead> <tbody> <tr> <td>FR700</td> <td>FITC and RPE</td> <td>X0949</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Control Reagent Code No.	FR700	FITC and RPE	X0949
Antibody Code No.	Fluorochrome	Control Reagent Code No.					
FR700	FITC and RPE	X0949					
Specificity	<p>Anti-CD45, T29/33, was included in the Third and Fourth International Workshops and Conferences on Human Leucocyte Differentiation Antigens. At the Third Workshop, the antibody was clustered as a CD45 antibody reacting with all the known isotypes of the CD45 family, also called the leucocyte common antigen family (4). At the Fourth Workshop, the expression of CD45 in a range of haematopoietic cell lines, and the lack of CD45 in non-haematopoietic cell lines, were demonstrated. Notably, CD45 was expressed at a very low level in myeloma cells and in the U-266 plasmacytoid cell line (5).</p> <p>Anti-CD14, TÜK4, was included in the Fourth, Fifth and Sixth International Workshop and Conference on Leukocyte Cell Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with the CD14 antigen (6-8).</p> <p>Anti-CD14, TÜK4 reacts strongly with peripheral blood monocytes and shows high reactivity to granulocytes. It also reacts with perivascular macrophages. The antibody labels 50% of Langerhans' cells in epidermal cell suspensions and is therefore characterized as weakly reactive (8).</p>						
Precautions	<ol style="list-style-type: none"> For professional users. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. As with any product derived from biological sources, proper handling procedures should be used. 						
Storage	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.						
Staining procedure	<ol style="list-style-type: none"> Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube. Add 10 µL of FR700 and mix gently by using a vortex mixer. Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes. Add 100 µL of Uti-Lyse™ Reagent A (DakoCytomation code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark. Add 1 mL of Uti-Lyse™ Reagent B (DakoCytomation code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube. Add 2 mL of PBS (DakoCytomation code No. S3024) to the tube and resuspend the cells by using a vortex mixer. Repeat step 6. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS. 						

10. Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Procedural notes

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the isotypes and fluorochromes of the conjugated antibody reagent. The recommended Dual-Colour Control Reagent for FR700 is DakoCytomation code No. X0949.

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, such as DakoCytomation EasyLyse™, code No. S2364, the PBS in step 9 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis.

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
1. Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 2003. p. 121-4.
2. Centers for disease control. Guidelines for the performance of CD4+ T-cell determinations in persons with human immunodeficiency virus infection. MMWR 1992;41 (No. RR-8);001. CDC Publication date 05/08/1992.
3. National Committee for Clinical Laboratory Standards (NCCLS) Document. Clinical applications of flow cytometry: Quality assurance and immunophenotyping of peripheral blood lymphocytes. H42-T 1992;12:1-76.
4. Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1, and leucocyte common antigens: new and previously defined clusters. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sep 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 788-803.
5. Schwitzer R. N7. Cluster report: CD45/CD45R. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 628-34.
6. Goyert SM, Haziot A, Jiao D, Katz IR, Kruger C, Ross J, Schütt C. M4. CD14 cluster workshop report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, Ritz J, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the Fifth International Workshop and Conference; Held in Boston, USA 3-7 November 1993. Volume One. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 778-82.
7. Goyert SM, Cohen L, Gangloff SC, Ashmun R, Haeffner-Cavaillon N. MC4. CD14 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 963-5.
8. Goyert SM, Tesio L, Ashman LK, Ball E, Bazil V, Garrido F, Hogg N, Horejsi V, et al. M3.1. Report on the CD14 cluster workshop. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, von dem Borne AEG, editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference on Human Leucocyte Differentiation Antigens; 1989 Februarv 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 789-94.

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