**MultiMix™**  
**Dual-Colour Reagent**  
**Anti-Human CD5/FITC**  
**Anti-Human CD20/RPE**

**Code No./ Code-Nr. FR729**

**ENGLISH**

**Intended use** For in vitro diagnostic use.  
 FR729 is intended for use in flow cytometry. FR729 allows simultaneous detection and enumeration of CD5+ cells and CD20+ cells. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

**Introduction** CD5 is present on all mature T cells and most thymocytes. It is also expressed on a small subpopulation of normal B cells. The majority of T-cell malignancies (76%) express CD5, and almost 85% of T-cell acute lymphoblastic leukaemias are CD5 positive (1).  
 Further, CD5 is expressed in some B-cell derived lymphoproliferative disorders, notably chronic lymphocytic leukaemia (CLL) (>90%) (1, 2), and centrocytic leukaemia (3).

CD20 can be considered a pan-B cell antigen, as it is expressed on the surface of all mature B lymphocytes but not in secreting plasma cells (4). CD20 is expressed early during pre-B cell development, presumably just before the expression of cytoplasmic  $\mu$ -chains, and it persists until plasma cell differentiation (5). A small subset of T cells expresses low levels of CD20 (6). CD20 is expressed in acute lymphoblastic leukaemia, B-cell chronic lymphocytic leukaemia, hairy cell leukaemia and Burkitt's lymphoma (7).

**Reagent provided** FR729 comprises the following two, carefully matched, fluorescent antibodies:  
 Monoclonal Mouse Anti-Human CD5, Clone DK23, conjugated with fluorescein isothiocyanate isomer 1 (FITC).  
 Monoclonal Mouse Anti-Human CD20, Clone B-Ly1, conjugated with R-phycoerythrin (RPE).  
 The two conjugates in FR729 have been produced from purified monoclonal mouse antibodies of isotype IgG1, kappa. FR729 is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 50 tests (10  $\mu$ L of conjugate for up to 10<sup>6</sup> leucocytes from normal human peripheral blood).

Antibody Code No.	Fluorochrome	Control Reagent Code No.
FR729	FITC and RPE	X0932

**Specificity** The specificity of Anti-CD5, DK23, is equivalent to that of the CD5-clustered antibody Leu-1 (3).  
 Anti-CD20, B-Ly1, was included in the Fourth International Workshop and Conference on Human Leucocyte Differentiation Antigens, and studies by a number of laboratories confirmed its reactivity with CD20 (5). Cross-blocking studies indicate that a majority of antibodies in the CD20 cluster studied at the above workshop defines a single epitope (8).

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

**Storage** Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Staining procedure**

1. Transfer 100  $\mu$ L of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.
2. Add 10  $\mu$ L of FR729 and mix gently by using a vortex mixer.
3. Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes.
4. Add 100  $\mu$ L of Uti-Lyse™ Reagent A (Dako code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
5. Add 1 mL of Uti-Lyse™ Reagent B (Dako code No. S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.
6. Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50  $\mu$ L of fluid in the tube.
7. Add 2 mL of PBS (Dako code No. S3024) to the tube and resuspend the cells by using a vortex mixer.
8. Repeat step 6.
9. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.

10. Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining.

Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

#### Procedural notes

Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.

Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.

It is optional to include a control reagent test tube. The control reagent should match the isotypes and fluorochromes of the conjugated antibody reagent. The recommended Dual-Colour Control Reagent for FR729 is Dako code No. X0932.

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, such as Dako EasyLyse™, code No. S2364, the PBS in step 9 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens. The correct use of colour compensation is particularly important in multicolour analysis.

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
---

## References/ Literatur

1. Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 1999. p. 51-2.
2. Lozano F, Simarro M, Calvo J. CD Guide. CD5. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 1112-3.
3. Stein H, Lennert K, Feller AC, Mason DY. Immunohistological analysis of human lymphoma: correlation of histological and immunological categories. *Adv Cancer Res* 1984;42:67-147.
4. Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology; London: Oxford University Press; 1999. p. 69-70.
5. Dörken B, Möller P, Pezzutto A, Schwartz-Albiez R, Moldenhauer G. B3. B-cell antigens: CD20. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 46-8.
6. Hultin LE, Hausner MA, Hultin PM, Giorgi JV. CD20 (pan-B cell) antigen is expressed at a low level on a subpopulation of human T lymphocytes. *Cytometry* 1993;14:196-204.
7. Pezzutto A, Behm F, Callard RE, Clark EA, Genetet N, Goodall AH, et al. B22.3. Flow cytometry analysis of the B-cell blind panel: joint report. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 165-74.
8. Tedder TF, Penta A. B3.1. Structure of the CD20 antigen and gene of human and mouse B-cells: use of transfected cell lines to examine the Workshop panel of antibodies. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 48-50.

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Agilent Technologies Singapore (International) Pte Ltd.