**FluoroSpheres**

Code K0110

5th edition/ 5ème édition/ 5. Ausgabe

Calibration beads for daily monitoring of the flow cytometer.

The kit contains reagents for 40 calibrations.

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Page/ Page/ Seite

### **ENGLISH**

Intended Use .....	4
Summary and Explanation.....	4
Principle of the Assay .....	4
Reagents.....	5
A. Materials provided .....	5
B. Materials required but not provided.....	5
Precautions .....	5
Storage.....	5
Assay Procedure and Data Acquisition.....	5
Data Analysis, Calibration Curve .....	6
A. Definitions.....	7
B. Manual calibration. For instruments expressing MFI in linear values .....	7
C. Manual calibration. For instruments expressing MFI in channel numbers.....	8
Explanation of symbols.....	22

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### **Intended Use**

For In vitro diagnostic use.

FluoroSpheres are intended for monitoring of flow cytometer performance.

### **Summary and Explanation**

The kit provides Blank Beads and Calibration Beads with a size of 3.2 µm. The Calibration Beads are a mixture of 5 bead populations having different fluorescence intensities, and one non-fluorescent bead population. Entrapped in the fluorescent beads is a combination of unique fluorochromes that enables excitation by light of any wavelength from 365 to 650 nm. Since the bead fluorochromes are different from the fluorochromes normally used in flow cytometry, the spectral properties are different. Therefore the beads cannot be used for setting up fluorescence compensation. With assigned molecules of equivalent fluorochrome (MEF) values for the fluorescent bead populations, the use of FluoroSpheres enables the transformation of arbitrary units of mean fluorescence intensity (MFI) into absolute units for quality assurance and reporting of flow cytometry data.

### **Principle of the Assay**

1. Firstly, Blank Beads are analysed to establish optimal instrument settings for all channels of interest.
2. Secondly, Calibration Beads are analysed. The data are used for the construction of the calibration curve, where MFI is plotted against MEF.

The calibration curve and associated statistical parameters are used to evaluate the linearity of the instrument. Furthermore, the detection threshold and dynamic range of the instrument can also be determined from the calibration curve.

Calibration with FluoroSpheres can be performed in two ways:

- Constant instrument setting.  
When using this method all instrument settings remain unchanged over time. The MFI for the different peaks are monitored and the results are displayed as date of calibration versus MFI.
- Constant instrument response.  
By this method the PMT voltage is adjusted so the MFI of one of the high fluorescent peaks appear in exactly the same channel every time. The results are then displayed as date of calibration versus PMT voltage.

## Reagents

### A. Materials provided

**Vial 1**            **BLANK BEADS**  
1.7 mL  
In 0.02% sodium azide, 0.01% NP40

**Vial 2**            **CALIBRATION BEADS**  
1.7 mL  
Contains Blank Beads and 5 bead populations bearing different quantities of fluorochromes expressed in MEF values for FITC, RPE, RPE-TxRed, RPE-Cy5 and APC, respectively. For lot-specific MEF values, please see the enclosed Analytical Value Sheet. In 0.02% sodium azide, 0.01% NP40.

### B. Materials required but not provided

General laboratory equipment for flow cytometry procedures.

Semi-logarithmic and double logarithmic 5-cycle graph paper.

Calculator.

Software programs for flow cytometry data analysis.

Spreadsheet programs like Microsoft Excel. These programs are optional.

## Precautions

1. For professional users.
2. This product contains sodium azide ( $\text{NaN}_3$ ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.

## Storage

Store kit at 2-8 °C. Do not use after expiration date stamped on vials. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the kit is suspected, contact Dako Technical Services.

## Assay Procedure and Data Acquisition

1. Mix the beads vigorously with either a vortex mixer or by shaking.
2. Add 0.5 mL of a suitable dilution buffer, for example phosphate-buffered saline, pH 7.2, to each of two test tubes. Add one drop from Vial 1, Blank Beads, to the first tube and one drop from Vial 2, Calibration Beads, to the second tube. Mix with a vortex mixer to ensure that the beads are in suspension.
3. Prior to data acquisition, logarithmic amplification for the fluorescence detectors of interest is selected, and the fluorescence compensation for all parameters is set to zero.
4. Apply the tube with the Blank Beads on the flow cytometer.
5. Adjust PMT voltages to allow the peak of the Blank Beads to appear in the first decade in the histogram.

6. Using forward scatter versus side scatter, set a live gate on bead singlets (see Figure 1).
7. Acquire a minimum of 5 000 gated events of the Blank Beads.
8. Apply the tube with Calibration Beads on the flow cytometer and acquire data for a minimum of 5 000 gated events.
9. Six different fluorescence intensity peaks should be visible in the fluorescence channel of interest (See Figure 2).

## Data Analysis, Calibration Curve

Analysis of results is confined to gated bead singlets clear of debris, as defined on a forward scatter versus side scatter dot plot (see Figure 1).

Individual markers or regions are used to define each of the 6 bead populations of the Calibration Beads (see Figure 2). The regions are used to determine the MFI of each bead population.


Figure 1. Forward scatter versus side scatter of Blank Beads. The gate has been set to collect bead singlets.


Figure 2. Histogram of FluoroSpheres Bead populations.

## A. Definitions

### Log offset

Log offset represents the fluorescence value in channel 0. Channel 0 is usually represented in manufacturers' software by arbitrary values, e.g. 1 for most Becton Dickinson software, and 0.1 for most Coulter software.

### AvgRes%

The average residual percent (AvgRes%) is a non-weighted measure of the instrument linearity of response. It is calculated by determining the average absolute percentage the regression line varies from the actual points.

### $r^2$

Coefficient of determination ( $r^2$ ) is a weighted measure of linearity. It is most useful when analysing log data.

### Detection threshold

Even though Blank Beads are without entrapped fluorochromes, they can give rise to a signal. This signal represents the background noise, and, thus, the detection level of the instrument.

### MEF

Molecules of equivalent fluorochrome is the amount of fluorochrome per bead.

### Log decade

Log decade is the calibrated full-scale range, or dynamic range, of the fluorescence parameter that is being calibrated. Flow cytometers are, generally, provided with amplifiers having a dynamic range of 3.5 or 4 log decades, depending on manufacturer. Decreased performance of the flow cytometer will affect the log decade.

## B. Manual calibration. For instruments expressing MFI in linear values

Construction of the calibration curve:

1. Determine the MFI for each of the 6 populations of the FluoroSpheres Calibration Beads.
2. Plot log MFI against the corresponding log MEF for the different peaks. Calculate the "best-fit" straight line according to the following formula:

$$\log(\text{MEF}) = a \times \log(\text{MFI}) + b$$

where a is the slope and b is the y-axis intercept, also called “log offset” or “zero channel value”.

3. Calculate the associated statistical parameters AvgRes% and  $r^2$  of the calibration curve.
4. The dynamic range of the amplifier is expressed in log decades. Log decades can be calculated from the regression line by subtraction of  $\log(\text{MEF}_{\min})$  from  $\log(\text{MEF}_{\max})$ .  $\log(\text{MEF}_{\min})$  and  $\log(\text{MEF}_{\max})$  are derived from the lowest and highest MFI values, respectively. In the present case MFI for  $\text{MEF}_{\min}$  and  $\text{MEF}_{\max}$  are  $10^0$  and  $10^4$ , respectively.


Figure 3. Calibration curve (regression line) based on FluoroSpheres.

Example of calculation:

$$\begin{aligned} \log \text{ decades} &= \log(\text{MEF}_{\max}) - \log(\text{MEF}_{\min}) \\ \log \text{ decades} &= \log(188\,735) - \log(30.8) \\ \log \text{ decades} &= 5.28 - 1.49 \\ \log \text{ decades} &= 3.79 \end{aligned}$$

Channels per decade:

$$\begin{aligned} \text{Since the number of channels are 1024, the number of channels per decade are:} \\ 1024/3.79 = 270 \end{aligned}$$

### C. Manual calibration. For instruments expressing MFI in channel numbers

Construction of the calibration curve:

1. Determine the MFI for each of the 6 populations of the FluoroSpheres Calibration Beads.
2. Plot MFI against  $\log(\text{MEF})$  and calculate the “best-fit” straight line according to the following formula:  

$$\log(\text{MEF}) = a \times \text{MFI} + b$$
 where a is the slope and b is the y-axis intercept.
3. Calculate the associated statistical parameters AvgRes% and  $r^2$  of the calibration curve.
4. The dynamic range of the amplifier is expressed in log decades. Log decades can be calculated from the regression line by subtraction of  $\log(\text{MEF}_{\min})$  from  $\log(\text{MEF}_{\max})$ . Log

( $MEF_{min}$ ) and  $\log(MEF_{max})$  are derived from the lowest and highest MFI values, respectively. In the present case MFI for  $MEF_{min}$  and  $MEF_{max}$  are 0 and 1023, respectively.


Figure 4. Calibration curve (regression line) based on FluoroSpheres.

Example of calculation:

$$\begin{aligned} \log \text{ decades} &= \log(MEF_{max}) - \log(MEF_{min}) \\ \log \text{ decades} &= \log(188735) - \log(30.8) \\ \log \text{ decades} &= 5.28 - 1.49 \\ \log \text{ decades} &= 3.79 \end{aligned}$$

Channels per decade:

Since the number of channels are 1024, the number of channels per decade are:  
 $1024/3.79 = 270$

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Dako Denmark A/S