IntraStain

Code No./ Code/ Code-Nr. K 2311

3rd edition / 3ème édition/ 3. Ausgabe

Fixation and permeabilization kit for flow cytometry.

The kit contains reagents for 100 tests.

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Page/ Page/ Seite

ENGLISH

Intended Use	3
Reagents	3
A. Materials provided	3
B. Materials required but not provided.....	3
Precautions	4
Storage.....	4
Sample Preparation and Staining Procedure	4
Reaction Patterns When Using the IntraStain Procedure	5
Limitations	5
References.....	12
Explanation of symbols.....	12

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Intended Use

For In vitro diagnostic use.

IntraStain is intended for fixation and permeabilization of single-cell suspensions for use in flow cytometry. The IntraStain procedure allows immunological detection of intracellular antigens, while the cellular structure, morphology scatter, and cell surface immunoreactivity remains intact. Cells treated with IntraStain can, therefore, be identified in flow cytometry by their light scatter properties and surface marker expression, while simultaneously being analysed for intracellular antigens (1, 2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Reagents

A. Materials provided

Vial 1

FIXATIVE

IntraStain Reagent A, Fixative

10 mL ready-to-use reagent. For 100 tests.

Vial 2

PERMEABILIZATION

IntraStain Reagent B, Permeabilization

10 mL ready-to-use reagent. For 100 tests.

B. Materials required but not provided

1. Fluorochrome-Conjugated Antibodies

Monoclonal or polyclonal fluorochrome-conjugated antibodies directed against intracellular or surface antigens.

2. Negative Controls

For monoclonal antibodies, DakoCytomation offers a range of isotype-matched FITC, RPE, RPE-Cy5 and APC-conjugated negative controls.

For polyclonal antibodies, DakoCytomation provides FITC and RPE-conjugated negative controls for the F(ab')₂ fragment of rabbit antibodies.

The optimal concentration of the negative control may vary depending on the specimen and preparation method, and should be determined by each individual laboratory.

3. Phosphate-Buffered Saline (PBS) Stock Solution (10 x concentrated)

80.0 g NaCl

0.2 g KH₂PO₄

14.4 g Na₂HPO₄, 2H₂O

2.0 g KCl

Add distilled water to 1 litre. Check pH. At a dilution of 1:10 the pH should be 7.4 ± 0.1.

4. PBS Working Solution (0.01 mol/L PBS)

Add 100 mL of PBS Stock Solution to 900 mL of distilled water.

5. 1% Paraformaldehyde in PBS

Add 1.0 g paraformaldehyde to 90 mL of 70 °C distilled water. Leave at 70 °C for 15-30 minutes in a fume hood. Add 10 mL of PBS Stock Solution and allow to cool. Store at 4 °C. Stable for no longer than three weeks.

General laboratory equipment for flow cytometry procedures.

Precautions

1. For professional users.
2. IntraStain Reagent A, Fixative, contains 5-10% formaldehyde, and is labelled:
Harmful.
R20/21/22 Harmful by inhalation, in contact with skin and if swallowed.
R36/37/38 Irritating to eyes, respiratory system and skin.
R40 Limited evidence of a carcinogenic effect.
R43 May cause sensitisation by skin contact.
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
S35 This material and its container must be disposed of in a safe way.
S36/37 Wear suitable protective clothing and gloves.
As a main rule, persons under 18 years of age are not allowed to work with this product. Users must be carefully instructed in the proper working procedure, the dangerous properties of the product and the necessary safety instructions. Please refer to the Material Safety Data Sheet (MSDS) for additional information.

Storage

Store kit at room temperature. Do not use after expiration date stamped on vials. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the product is suspected, contact our Technical Services.

Sample Preparation and Staining Procedure

1. Transfer 50 μL (up to 10^6 cells) of the cell suspension to be analysed (whole blood, bone marrow or mononuclear cells) to each of two test tubes, labelled 1 and 2 (test and control).
2. *This step is omitted if only intracellular staining is performed.* To test tube 1, add an appropriate volume of fluorochrome-conjugated antibody specific for **the cell surface antigen** to be stained. To test tube 2, add an appropriate volume of a suitably matched fluorochrome-conjugated negative control. Incubate for 15 minutes in the dark at room temperature.
3. Add 100 μL IntraStain Reagent A, Fixative, to each test tube. Mix gently with a vortex mixer to ensure that the cells are in suspension.
4. Incubate at room temperature for 15 minutes.
5. Add 2 mL PBS Working Solution to each test tube and mix gently.
6. Centrifuge at 300 x g for 5 minutes, then aspirate the supernatant, leaving approximately 50 μL of fluid. Mix thoroughly to ensure that the cells are in suspension.
7. Add 100 μL IntraStain Reagent B, Permeabilization, to each test tube. To test tube 1, add an appropriate volume of fluorochrome-conjugated antibody specific for **the intracellular antigen** to be stained. To test tube 2, add an appropriate volume of a suitably matched fluorochrome-conjugated negative control. Mix gently to ensure that the cells are in suspension.
8. Incubate in the dark at room temperature for 15 minutes.
9. Repeat steps 5 and 6.
10. Resuspend the pellet in an appropriate fluid for flow cytometric analysis.
11. Analyse on a flow cytometer. If the samples are not analysed immediately, they may be stored in the dark at 2-8 $^{\circ}\text{C}$ for up to 8 hours.

NOTE: If the cell-surface antigen is formalin resistant, it is possible to incubate the cells with the cell-surface-reactive and the intracellular-reactive antibodies simultaneously. Thus, omit step 2 only, and add in step 7 both antibodies to test tube 1, and both negative controls to test tube 2.

Reaction Patterns When Using the IntraStain Procedure


Figure 1. Forward scatter versus side scatter representation of a normal blood sample.


Figure 2. Surface (CD19) versus intracellular (CD79 α cy) antigen. Cells from a normal blood sample labelled with Anti-CD19/FITC, Clone HD37, code No. F 0768, and Anti-CD79 α cy/RPE, Clone HM57, code No. R 7159.


Figure 3. Intracellular antigen (MPO) versus side scatter. Cells from a normal blood sample labelled with Anti-MPO/FITC, Clone MPO-7, code No. F 0714.


Figure 4. A mixture of 70% CEM and 30% RAJI cells (which both contain intracellular CD3) labelled with Anti-CD3/RPE, clone UCHT1, code No. R 0810.

Limitations

Most fluorochrome-conjugated antibodies directed against intracellular or surface antigens are suitable for use together with IntraStain. Some antigen determinations, however, are sensitive to fixation with formaldehyde. In the test procedure, particular attention should be paid to the formalin resistance of the intracellular antigen.

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
References / Références / Literatur

1. Kappelmayer J, Gratama JW, Karászi É, Menéndez P, Ciudad J, Rivas R, et al. Flow cytometric detection of intracellular myeloperoxidase, CD3 and CD79a. Interaction between monoclonal antibody clones, fluorochromes and sample preparation protocols. J Immunol Methods 2000;242:53-65.
2. Verdier M, Jayat C, Ratinaud M-H, Troutaud D. Optimization of cell permeabilization for multiparametric flow cytometric analysis with lectin staining. Cytometry 2000;41:55-61.

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DakoCytomation Denmark A/S