**Monoclonal Mouse  
Anti-Human  
CD45, Leucocyte Common Antigen/PerCP,  
Clone 2D1**

**Code PR701**

ENGLISH							
<b>Intended use</b>	For in vitro diagnostic use. PR701 is intended for use in flow cytometry. The antibody is intended for use in the identification of cells expressing CD45. CD45 is one of the most abundant leucocyte cell surface glycoproteins and is expressed exclusively on cells of the haematopoietic system and their progenitors (1). In flow cytometry, anti-CD45, together with a panel of other antibodies, is considered essential for the initial evaluation of chronic lymphoproliferative disorders and acute leukaemias (2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.						
<b>Summary and explanation</b>	CD45 is a single chain type I transmembrane protein typically expressed at high levels on nucleated cells of haematopoietic origin (3). Thus CD45 is present on T and B lymphocytes, granulocytes, monocytes and macrophages, with the exception of maturing erythrocytes and megakaryocytes (1). There are five different isoforms of CD45, based on differential splicing of exons 4, 5 and 6, named ABC, AB, BC, B and 0. The Mr of the isoforms ranges from 220 000 for the ABC isoform to 180 000 for the 0 isoform. All the CD45 isoforms share the same intracellular segment which has been shown to have tyrosine phosphatase activity and which has a functional role in lymphocyte activation and differentiation (3). Antibodies that recognize all five isoforms are known as anti-CD45 (3).						
<b>Reagent provided</b>	The Anti-CD45 conjugate, PR701, has been produced from a purified monoclonal mouse antibody. The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 <sup>6</sup> leucocytes from normal human peripheral blood). <u>Isotype</u> : IgG1, kappa. <u>Conjugate concentration mg/L</u> : See label on vial.						
	<table border="1"> <thead> <tr> <th>Antibody Code</th> <th>Fluorochrome</th> <th>Control Reagent Code</th> </tr> </thead> <tbody> <tr> <td>PR701</td> <td>Peridinin Chlorophyll Protein (PerCP)</td> <td>X7909</td> </tr> </tbody> </table>	Antibody Code	Fluorochrome	Control Reagent Code	PR701	Peridinin Chlorophyll Protein (PerCP)	X7909
Antibody Code	Fluorochrome	Control Reagent Code					
PR701	Peridinin Chlorophyll Protein (PerCP)	X7909					
<b>Immunogen</b>	Ficoll-Triosil-separated human peripheral blood mononuclear cells (4, 5).						
<b>Specificity</b>	Anti-CD45, 2D1, was included at the Third International Workshop and Conference on Human Leucocyte Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with CD45 (6, 7).						
<b>Precautions</b>	<ol style="list-style-type: none"> <li>For professional users.</li> <li>This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>As with any product derived from biological sources, proper handling procedures should be used.</li> <li>Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.</li> <li>Unused solution should be disposed of according to local, State and Federal regulations.</li> </ol>						
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.						
<b>Staining procedure</b>	<ol style="list-style-type: none"> <li>Transfer 100 µL of anticoagulated (EDTA) blood to a 12 mm x 75 mm polystyrene test tube.</li> <li>Add 10 µL of PR701 and mix gently by using a vortex mixer.</li> <li>Incubate the tube in the dark at 2-8 °C for 30 minutes or at room temperature (20-25 °C) for 15-30 minutes.</li> <li>Add 100 µL of Uti-Lyse™ Reagent A (Dako Code S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.</li> <li>Add 1 mL of Uti-Lyse™ Reagent B (Dako Code S3325) to the tube and mix gently by using a vortex mixer. Incubate the tube for 10 minutes at room temperature in the dark.</li> <li>Centrifuge the tube at 300 x g for 5 minutes. Gently aspirate the supernatant and discard it leaving approximately 50 µL of fluid in the tube.</li> <li>Add 2 mL of PBS to the tube and resuspend the cells by using a vortex mixer.</li> <li>Repeat step 6.</li> <li>Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL PBS.</li> <li>Analyse the sample on a flow cytometer or store it in the dark at 2-8 °C until analysis. Samples should be analysed within 24 hours after staining.</li> </ol> <p>Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.</p>						
<b>Procedural notes</b>	<p>Step 1: It is optional to include a suitable positive and negative control sample with each run for reagent and preparation control.</p> <p>Step 2: The volume of conjugate recommended is a guideline only. Optimal staining conditions may vary depending on specimen and preparation method and should be determined by each individual laboratory.</p> <p>It is optional to include a control reagent test tube. The control reagent should match the isotype and fluorochrome of the conjugated antibody. The recommended control reagent is shown in the table above.</p>						

(119154-001)

PR701/EFG/LHP/2008.09.17 p. 1/4

Steps 4 and 5: If another cell-lysing reagent is used, please follow the recommendations for that reagent. Note that if the alternative lysing reagent does not contain fixative, such as Dako EasyLyse™, Code S2364, the PBS in step 9 should contain 1% paraformaldehyde unless the sample is analysed within the time frame recommended for the lysing reagent.

Step 10: Abnormal numbers of cells expressing the target antigen or cells showing aberrant expression levels of the antigen can be expected in some disease states. This may result in an altered pattern of staining. It is important to understand the normal expression pattern of an antigen and its relationship to the expression of other relevant antigens in order to perform appropriate analysis.

Multicolour reagents are preferable to single-colour reagents for the comprehensive analysis of flow cytometry specimens.

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

- Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 1999. p. 95-8.
- Braylan RC, Orfao A, Borowitz MJ, Davis BH. Optimal number of reagents required to evaluate hematology neoplasias: results of an international consensus meeting. Cytometry 2001;46:23-7.
- Sewell WA, Cooley MA, Hegen M. NL6. CD45 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 499-502.
- Beverley PCL. Production and use of monoclonal antibodies in transplantation immunology. In: Touraine JL, Trager J, Betuel H, editors. Transplantation and clinical immunology XI. Excerpta Medica Amsterdam. 1980. p. 87-94.
- Bradstock KF, Janossy G, Pizzolo G, Hoffbrand AV, McMichael A, Pilch JR, et al. Subpopulations of normal and leukemic human thymocytes: an analysis with the use of monoclonal antibodies. J Natl Cancer Inst 1980;65:33-42.
- Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1 family, and leucocyte common antigens: new and previously defined clusters. In: McMichael AJ, Beverley PCL, Cobbold S, Crumpton MJ, Gilks W, Gotch FM, et al., editors. Leucocyte typing III. White cell differentiation antigens. Proceedings of the 3rd International Workshop and Conference; 1986 Sept 21-26; Oxford, England. Oxford, New York, Tokyo: Oxford University Press; 1987. p. 788-803.
- Sewell WA, Cooley MA, Katz KS. CD Guide. CD45. In Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 793-4.

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