## Monoclonal Mouse Anti-Human CD33, Clone WM-54

Code No./ Code/ Code-Nr. **F 0832**      **FITC-Conjugated**  
 Code No./ Code/ Code-Nr. **R 0745**      **RPE-Conjugated**

ENGLISH										
<b>Intended use</b>	For in vitro diagnostic use. F 0832 and R 0745 are intended for use in flow cytometry. Anti-CD33 antibody is considered essential for the initial evaluation of acute myeloid leukaemia (AML) together with a panel of other antibodies (1, 2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.									
<b>Synonym for antigen</b>	gp67, siglec-3 (3, 5).									
<b>Introduction</b>	CD33 is a member of the siglec family (sialic acid binding immunoglobulin-like lectins) referred to as Siglec-3 (5). It is the smallest member of the siglec family having a molecular weight of 67 kDa (4). CD33 mediates protein-carbohydrate cellular interactions in addition to protein-protein recognition (5). The CD33 gene is located on chromosome 19 (4). CD33 is expressed by subsets of monocytes and neutrophils. Committed hemopoietic progenitors, including multipotential stem cells (CFU-GEMM) and precursors of granulocytes and macrophages (CFU-GM) and early erythroid progenitors (BFU-E) express CD33 (3, 4). The antigen has been found on the cell surface of leukaemic blasts from the vast majority of cases of AML, including megakaryoblastic, but not erythroblastic variants. Leukaemic lymphoblasts do not show any expression of CD33 (1-4, 6).									
<b>Reagent provided</b>	The Anti-CD33 conjugates, F 0832 and R 0745, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN <sub>3</sub> , pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 10 <sup>6</sup> U937 cells). <u>Isotype:</u> IgG1, kappa. <u>Conjugate concentration mg/L:</u> See label on vial.									
	<table border="1"> <thead> <tr> <th>Antibody Code No.</th> <th>Fluorochrome</th> <th>Negative Control Code No.</th> </tr> </thead> <tbody> <tr> <td>F 0832</td> <td>FITC (Fluorescein Isothiocyanate Isomer 1)</td> <td>X 0927</td> </tr> <tr> <td>R 0745</td> <td>RPE (R-Phycoerythrin)</td> <td>X 0928</td> </tr> </tbody> </table>	Antibody Code No.	Fluorochrome	Negative Control Code No.	F 0832	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927	R 0745	RPE (R-Phycoerythrin)	X 0928
Antibody Code No.	Fluorochrome	Negative Control Code No.								
F 0832	FITC (Fluorescein Isothiocyanate Isomer 1)	X 0927								
R 0745	RPE (R-Phycoerythrin)	X 0928								
<b>Immunogen</b>	Leukaemic cells from a patient with AML (6).									
<b>Specificity</b>	Anti-CD33, WM-54, was included in the Fourth International Workshop and Conference on Leukocyte Cell Differentiation Antigens and studies by a number of laboratories confirmed its reactivity with the CD33 antigen (3). Anti-CD33, WM-54, reacts with cells of the monocyte-macrophage lineage and granulocyte progenitors (3). The antibody also shows a faint to weak reactivity with mature granulocytes (6). The antibody does not react with erythrocytes, platelets or lymphoid cells (6). Bone marrow progenitors CFU-GM, CFU-Mix and BFU-E demonstrate positive staining with Anti-CD33, WM-54 (6). In AML, Anti-CD33, WM-54, has been shown to label 14/21 cases (6). In chronic myeloid leukaemia with myeloid blast crisis Anti-CD33, WM-54, labelled 3/3 cases (6). In common acute lymphoblastic leukaemia Anti-CD33, WM-54, has been shown to label 1/8 cases (6).									
<b>Precautions</b>	<ol style="list-style-type: none"> <li>For professional users.</li> <li>This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>As with any product derived from biological sources, proper handling procedures should be used.</li> </ol>									
<b>Storage</b>	Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.									
<b>Staining procedure</b>	<ol style="list-style-type: none"> <li>Collect venous blood into a test tube containing an anticoagulant.</li> <li>Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.</li> <li>Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.</li> <li>Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD33.</li> <li>Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).</li> <li>Incubate in the dark at 4 °C for 30 minutes.</li> <li>Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.</li> <li>Analyse on a flow cytometer.</li> </ol>									

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

**Product-specific limitations** The fluorescence intensity of the antibody is reduced by the addition of 1% paraformaldehyde.

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
1. van Dongen JJM, Adriaansen HJ. Immunobiology of leukemia. In: Henderson ES, Lister TA, Greaves MF, editors. Leukemia. Philadelphia, London, Toronto, Montreal, Sydney, Tokyo: WB Saunders Company; 1996. p. 83-130.
2. Baer MR, Stewart CC, Dodge RK, Leget G, Sulé N, Mrózek K, et al. High frequency of immunophenotype changes in acute myeloid leukemia at relapse: implications for residual disease detection (Cancer and Leukemia Group B Study 8361). Blood 2001;97:3574-80.
3. Peiper SC, Leboeuf RD, Hughes CB, Prasthofer EF, Borowitz MJ, Dewutter-Dambuyant C, et al. M6.1. Report on the CD33 cluster workshop: biochemical and genetic characterization of gp67. In: Knapp W, Dörken B, Gilks WR, Rieber EP, Schmidt RE, Stein H, et al., editors. Leucocyte typing IV. White cell differentiation antigens. Proceedings of the 4th International Workshop and Conference on Human Leucocyte Differentiation Antigens; 1989 Feb 21-25; Vienna, Austria. Oxford, New York, Tokyo: Oxford University Press; 1989. p. 814-6.
4. Peiper SC, Guo HH. MC6. CD33 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyet SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 972-4.
5. Crocker PR, Cornish A, Floyd H, Jones C, Nicoll G, Zhang J. AS6.1. Siglecs-an emerging family of sialic acid binding receptors of the immunoglobulin superfamily. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 31-3.
6. Favaloro EJ, Bradstock KF, Kabral A, Grimsley P, Berndt MC. Characterization of monoclonal antibodies to the human myeloid-differentiation antigen "gp67" (CD-33). Disease Markers 1987;5:215-25.

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Agilent Technologies Singapore (International) Pte Ltd.