ENGLISH

Intended use For in vitro diagnostic use.
 F0801 and R0811 are intended for use in flow cytometry. The antibody labels CD25-expressing leucocyte subtypes (1, 2). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Synonyms for antigen Interleukin-2 receptor (IL-2R), interleukin-2 receptor α chain, Tac antigen (3, 4, 5).

Summary and explanation The interleukin-2 receptors consist of three polypeptides (α , β , γ), which in different combinations bind IL-2 (1). The 55 kDa α chain is called CD25 (5). The α chain alone, the $\beta\gamma$ heterodimer, and the $\alpha\beta\gamma$ heterotrimer form the low-, intermediate- and high affinity IL-2R, respectively. The intermediate- and high affinity IL-2R play a role in intracellular signal transduction (3).

CD25 is strongly expressed on the surface of T cells activated by antigens, mitogens (phytohemagglutinin (PHA), concanavalin A (Con-A), pokeweed mitogen (PWM)), *Staphylococcus aureus*, viruses (HTLV-I, HTLV-II, Epstein-Barr virus) (6, 7), and is also expressed on B cells stimulated with anti-IgM and on monocytes/macrophages stimulated with lipopolysaccharide (2). In normal peripheral blood a variable staining of T lymphocytes is observed (8).

Reagent provided The Anti-CD25 conjugates, F0801 and R0811, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains conjugate for 100 tests (10 μ L of conjugate for up to 10⁶ cells).

Isotype: IgG1, kappa. Conjugate concentration: see label on vial.

Table 1

Antibody Code	Fluorochrome	Negative Control Code
F0801	FITC (Fluorescein Isothiocyanate Isomer 1)	X0927
R0811	RPE (R-Phycoerythrin)	X0928

Immunogen Phytohemagglutinin-stimulated peripheral blood lymphocytes.

Specificity The specificity of Anti-CD25, ACT-1, is identical to that of CD25-clustered antibodies as indicated by analysis of immunoprecipitates and the immunohistochemical labeling pattern (1, 9).

- Precautions**
1. For in vitro diagnostic use.
 2. For professional users.
 3. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
 4. As with any product derived from biological sources, proper handling procedures should be used.
 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
 6. Unused solution should be disposed of according to local, State and Federal regulations.

Storage Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

- Staining procedure**
1. Collect venous blood into a test tube containing an anticoagulant.
 2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
 3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
 4. Mix 100 μ L cell suspension with 10 μ L fluorochrome-conjugated Anti-CD25.
 5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome as a negative control (see table 1).
 6. Incubate in the dark at 4 °C for 30 minutes.
 7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.
 8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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

References/ Bibliographie/ Literaturnachweise

1. Schwarting R, Gerdes J, Ziegler A, Stein H. Immunoprecipitation of the interleukin-2-receptor from Hodgkin's disease derived cell lines by monoclonal antibodies. *Hematol Oncol* 1987;5:57-64.
2. Sugamura K. CD25. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. *Leukocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan.* New York, London: Garland Publishing Inc.; 1997. p. 1127-8.
3. Ishii N, Kondo M, Takeshita T, Sugamura K. CR2.1. mAb specific for the γ chain of the IL-2 receptor. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. *Leukocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA.* Oxford, New York, Tokyo: Oxford University Press; 1995. p.1867-8.
4. Leonard WJ, Depper JM, Robb RJ, Waldmann TA, Greene WC. Characterization of the human receptor for T-cell growth factor. *Proc Natl Acad Sci (USA)* 1983;80:6957-61.
5. Kikutani H, Kishimoto T. CR1. The cytokine receptors: section report. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. *Leukocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA.* Oxford, New York, Tokyo: Oxford University Press; 1995. p.1855-64.
6. Uchiyama T, Broder S, Waldmann TA. A monoclonal antibody (anti-Tac) reactive with activated and functionally mature human T cells. I. Production of anti-Tac monoclonal antibody and distribution of Tac (+) cells. *J Immunol* 1981;126:1393-7.
7. Uchiyama T, Nelson DL, Fleisher TA, Waldmann TA. A monoclonal antibody (anti-Tac) reactive with activated and functionally mature T cells. II. Expression of Tac antigen on activated cytotoxic killer T cells, suppressor cells, and on one of two types of helper T cells. *J Immunol* 1981;126:1398-1403.
8. Baecher-Allan C, Brown JA, Freeman GJ, Hafler DA. CD4+ CD25high regulatory cells in human peripheral Blood. *J Immunol* 2001; 167:1245-1253.
9. Aziz KE, McCluskey PJ, Wakefield D. Expression of selectins (CD62 E,L,P) and cellular adhesion molecules in primary Sjögren's syndrome: questions to immunoregulation. *Clin Immunol Immunopathol* 1996;80:55-66.

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Agilent Technologies Singapore (International) Pte Ltd.